


Owner's Manual



➤ A903 BATH LIFT

Read this manual before
using the Bath Lift and
save for future reference.

Table of Contents

0. IEC SYMBOLS	2
1. CAUTIONS	2
2. WARNINGS	3
3. PRECAUTIONS	4
4. INTENDED USE	4
5. DESCRIPTION	5
6. SPECIFICATIONS	6
7. HANDLING	8
7.1 Unpacking	8
8. OPERATION	9
8.1 Assemble & separate the backrest	9
8.2 Install the bath lift	9
8.3 Connect the battery Pack to the bath lift	10
8.4 Correct seating	11
8.5 Operating the handheld controller	11
8.6 Removing the bath lift from the bathtub	12
9. RECHARGING THE BATTERY PACK	13
10. MAINTENANCE	14
10.1 General caring	14
10.2 Fitting and removing the soft mats	14
10.3 Cleaning and washing instructions	17
11. TROUBLESHOOTING	17
Instruction Regarding EMC	18
Warranty	22

0. IEC SYMBOLS



Caution, Attention or Consult accompanying documents.



Consult Instructions for Use



Alternating Current



Direct Current




Type B Equipment



May contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.



Weight Capacity. DO NOT overload the Bath Lift.



Battery Recharge Indicator



1. CAUTIONS

- Bathing safely with the Merits Bath Lift is only ensured when the operating instructions are read carefully and the given information is strictly observed.
- Never use the bath lift as an aid to get in or out of anything except for the bathtub, nor as a support, a ladder, lifting platform or for other similar purposes.



2. WARNINGS

- Persons with severe physical handicaps or mentally disturbed persons should never be left unattended in the bath.
- The seat may slip if the Suction Pads do not adhere properly. Clean the bathtub and the Suction Pads before use. AND DO NOT move the Bath Lift once you have installed it in the bathtub. Replace the Suction Pads whenever they are damaged or worn out.
- DO NOT USE the Bath Lift without putting on the Soft Mat.
- Always sit in the middle of the Bath Lift and keep your weight centering in the seat while moving in or out of the bathtub and bathing.
- Only raise or lower the Actuator when the backrest is upright and securely locked.
- There is a risk of injury when the lifting scissors move. Do not reach between the lifting scissors during assembling it.
- While raising or lowering the Bath Lift, do not place your hands between the edge of the bathtub and the Bath Lift or under the seat of it, as this can cause injuries.
- The Bath Lift must be inaccessible to children and animals.
- Use only the Battery Charger that comes with your Bath Lift.
- The Battery Charger must only be used outside the bathroom in a dry environment.
- DO NOT connect the Battery Charger to the power sources until the Battery Pack has been correctly connected with it.
- Do not place the Battery Pack near hot, sparking or burning objects.
- Always disconnect the Battery Pack before transporting the Bath Lift.
- Batteries containing in the Battery Packs are rechargeable and replaceable, but can only be replaced by qualified personnel. DO NOT try to open or replace them by yourself. The used Batteries must be disposed of in the proper manner.
- Repairs must only be undertaken by authorized specialists or the manufacturer.
- The Actuator is not designed to be used under water level. In order not to cause any danger or damage. DO NOT operate the Bath Lift while the Actuator is submerged in the water.
- Do NOT use a damaged or badly worn Bath Lift.
- DO NOT overload the Bath Lift. Otherwise, serious injury may occur.


3. PRECAUTIONS

- The batteries in the Battery Pack have to be charged for at least 8 hours before its first use. And should be recharged after each use.
- The Bath Lift must be cleaned directly after use (preferably under a hot shower) in order to avoid unforeseen defects and to maintain it in a perfect, hygienic condition.
- If your bathtub has grasp bars, make sure to fit the Guide Straps on them. The Bath Lift can be damaged if Guide Straps are not fitted to the bathtub grasp bars properly.
- There is a risk of damaging the Actuator. Never clean the backrest with a high pressure cleaner. Or submerge the entire Actuator into the water.
- Do not use any detergents containing acetic acid for cleaning the Bath Lift.

4. INTENDED USE

The Merit bath lift is intended for lowering and lifting a person in a bathtub. Never use the bath lift as an aid to get in or out of anything except for the bathtub, nor as a support, a ladder, lifting platform or for other similar purposes. Any other uses which this bath lift does not intend for may cause serious injury or damage to personal property.


5. DESCRIPTION





The Merit Bath Lift is basically consists of three parts, the first part is the Base Unit, the second part is the Backrest and the third is the Battery Pack. The Base Unit with sturdy collapsible frame enables the chair to be raised and lowered from 6.5 to 42 cm. An actuator is located on the back of the Backrest. As this actuator stretching and retracting, the Bath Lift can be lifted up, lowered down or adjusted the slope of the Backrest. The movement of this actuator is controlled by the Handheld Controller attached with the Backrest. The Battery Pack which has Ni-MH batteries inside offers the power needed to drive the actuator. Due to its integrated power control, the Merits Bath Lift only lowers if there is still sufficient power available to raise it. The device requires no maintenance tasks other than the routine care explained in this manual. Only a maintenance check is recommended every 2 years to ensure the safety of usage.

6. SPECIFICATIONS

Dimensions:		Inches	mm
Total height including backrest (B)		40.7	1,034
Backrest height (F)		21.5	545
Backrest width (E)		15.6	397
Seat Height (A)	Highest	16.5	422
	Lowest	2.6	65
Seat width (C)	Excluding side flaps	17.4	442
	Including side flap opened	27.6	700
Sitting area length (D)		18.2	462
Cover mat thickness maximum		0.12	3



Technical Specs:	
Weight Capacity	136 kg (300lbs)
Total Weight (with Handheld Controller)	14 kg (30.8 lbs)
Weight of Seat	8.5 kg (18.7 lbs)
Weight of Backrest	5.0 kg (11.0 lbs)
Weight of Handheld Controller	0.5 kg (1.1 lbs)
Shipping Weight	17 kg (37.4 lbs)
Rated Current Consumption	4.8 A (12V )
Battery	Ni-MH / 14.4V / 1300mAh
Fuse (at Battery Pack & Main Control Board)	10A/250V
Operating Force for Hand Controller (Max)	5N
AC/DC Adapter (Use only the original adapter)	Model: MR101 Input: 100 - 240V~ 50-60Hz (0.2A) Output.: 17.6V  170mA
Environmental Protection	Handheld Controller: IP67 Actuator: IP66
Sound Level (Average)	< 45dB(A)
Actuator Duty Cycle	MAXIMUM 10% = 1 Bath every 15 minutes
Operating Ambient Temperature	10°C ~ 40°C
Hydrostatical Pressure	860~1060hpa (12.5~15.4psi)
Recommended Storage and Shipping Temperature	Dry (Up to 95% non-condensing), well ventilated area -20°C ~ 60°C
Regulatory Listing	CSA 601/ UL 2601/ IEC/EN 60601-1/ EN/ISO 10535
Material Used in Main Components	ABS: Seat, Backrest & Headrest PVC & Polyester Mesh: Anti-slip Mat


- (1). Conditions may vary based on altitude above sea level, barometric pressure and temperature. Specification subjects to change without notice.
- (2). This product complies with the following standards: C22.2 No 601.1-M90, UL Std. No 2601-1, IEC/EN 60601-1 and EN/ISO 10535.
- (3). For electromagnetic compatibility the standards are based as follow: CSA 601.1-2 and IEC/EN 60601-1-2.
- (4). The level of protection against electric shock is **Class II** for the battery charger / internally powered equipment for the Bath Lift.
- (5). Type of protection of applied parts against electric shock is **B**.
- (6). Mode of operation belongs to intermittent mode.
- (7). Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

7. HANDLING

The Bath Lift should always be kept in the upright position to prevent cabinet damage while being transported. The shipping container has been designed to assure maximum protection of the Bath Lift. If the Bath Lift is to be reshipped by common carrier, additional cartons are available from Merits.

7.1 UNPACKING

1. Check for any obvious damage to the carton or its contents. If damage is evident, notify the carrier, or your local dealer.
2. Remove all loose packing from the carton.
3. Carefully remove all the components from the carton. There are:
 - a. The Bath Lift Unit with a Handheld Controller and removable anti-slip soft mat.
 - b. The Battery Pack. The batteries contained inside are rechargeable and replaceable, but can only be replaced by qualified personnel.
 - c. The Battery Charger.
 - d. The Headrest with removable anti-slip soft mat. And two (2) Hand Knobs for fixing it.
 - e. Two (2) Guide Strap.



NOTE: Unless the bath lift is to be used immediately, retain container and packing materials for storage until use of bath lift is required.

8. OPERATION

8.1 ASSEMBLE & SEPARATE THE BACKREST

For ease of handling and storage, the backrest may be separated from the Bath Lift. Lift the backrest straightly by using the actuator as a handle, and lift it away from the base unit, see Fig. 1. In order to prevent from mechanism opening and causing injury, you should always hold the base unit by the front and rear handles while moving it.

The backrest may be re-assembled by aligning the openings on the backrest with the lugs of the hinge on the base unit and rotate the backrest to slide the actuator pin into the locker to fix the position.




Fig1.

8.2 INSTALL THE BATH LIFT

1. Always install the Bath Lift to the bathtub in 2 sections: Backrest & Base unit. Follow the procedures described in Sect. 8.1 to separate the Backrest.
2. Fold the Side Flippers onto the Seat.
3. Fit your Bath Lift into the bathtub as shown in Fig. 2 (Install the longitudinal direction of the Bath Lift parallel with the longitudinal direction of bathtub), and fix its position by pressing the four powerful suction pads. Before every use, make sure:
 - a.) The bathtub is clean, and free from debris. Check to confirm effective adhesion of the suction pads. As the suction pads will not adhere to non-slip or textured surfaces. Be sure to avoid these kinds of surfaces.
 - b.) Allow sufficient space behind the Bath Lift for the backrest to recline. Ensure no obstacles will impede the Bath Lift. Such as taps, overflows, and soap dishes.
 - c.) Fitting guide straps - where bathtub grasp bars are present. The Bath Lift can be damaged if Guide Straps are not fitted to the bathtub grasp bars correctly. (Ref. to Fig. 3.)




Fig. 2. Fit the Bath Lift into the bathtub.




Fig. 3. Fit the guide straps.

4. Follow the procedures described in Sect. 8.1 to reassemble the Backrest.

5. Make sure that the actuator is in its shortest position. Rotate the actuator to set the Backrest to its final reclined position as shown in Fig. 4, and lock the Roller Shaft into the Actuator Locker.




Fig. 4. Set the Backrest to its final reclined position.




Fig. 5.

8.3 CONNECT THE BATTERY PACK TO THE BATH LIFT



WARNINGS

- Do not connect the Battery Pack to the Bath Lift until it is assembled and is set to its reclined position in the bath.
- The Battery Pack **MUST** be fully charged before use.

Align and slide the Battery Pack onto the Battery Seat on the backrest as shown in Fig. 5. to connect the Battery Pack. Lift the clip on the lower portion of the Battery Pack and pull it away from the Battery Seat to disconnect.

8.4 CORRECT SEATING



WARNINGS

- Before using the Bath Lift consider means of requesting help if needed.
- Seek advice from a qualified person such as an Occupational Therapist, if you are in any way concerned over the correct procedure for using the Bath Lift.
- Avoid using the seat back as a means to steady yourself.

1. Facing outwards, sit on to the Bath Lift side panel & then move on to the seat area.
2. Slowly lift your legs separately over into the bathtub, ease yourself into the central area of the seat.
3. Face forwards and have your legs and feet extended. Always keep your hands away from the seat, folding side panels & mechanism, when operating the Bath Lift.

8.5 OPERATING THE HANDHELD CONTROLLER

1. Make sure that the fully charged Battery Pack is installed firmly to the Battery Seat on the Backrest.
2. Press the DOWN button to lower and the UP button to raise the Bath Lift. If you wish to recline the backrest further, when the seat has reached its lower position continue to press the down button, this will recline the backrest further. When the desired angle is reached release the button.
3. To raise the backrest, press the UP button. Continue to press the UP button to raise the seat to its original position - level with the bathtub edges.




Fig. 6. Handheld Controller

【Note】

The Handheld Controller can be fitted on to a smooth surface by using the suction pads underneath. To release just pull the edge rim or side strap of the suction pads.

The Handheld Controller has a built in safety feature.

If the batteries do not have sufficient power to lower and then raise the Bath Lift - the unit WILL NOT DESCEND. A red warning light called Battery Recharge Indicator will be blinking on the Handheld Controller to indicate low battery power. See section 9 for recharging the Battery Pack.

**EMERGENCY STOP**

- If the push buttons on the Handheld Controller does not respond while the Bath Lift needs to be stopped, press the Emergency Stop Button on the Handheld Controller. A red led will be illuminated. And the Bath Lift will stop immediately. To operate the Bath Lift, press this button again to release it.
- **THIS PROCEDURE MUST BE USED ONLY IN EMERGENCIES.**
- **NEVER TRY TO OPEN OR DISASSEMBLE THE HANDHELD CONTROLLER. As this may break the seal of the Handheld Controller and damage the unit.**

8.6 REMOVING THE BATH LIFT FROM THE BATHTUB

1. Let the water leak out of the bathtub.
2. Always remove the Bath Lift from the bathtub in 2 sections: Backrest & Base unit. To remove the backrest ensures the seat is at its fully reclined position before folding and removing the backrest. Gently fold the backrest over until it touches the seat. Raise it again straightly and you can lift it away. (Follow the procedures described in Sect. 8.1 to separate the Backrest.)
3. Fold the Side Flippers onto the Seat.
4. First, lift the Bath Lift at the Front Handle and pull this handle to release the front Suction Pads.
When they are free, gently tilt the Bath Lift by the Front Handle.
DO NOT RELEASE STRAP OF THIS HANDLE.
5. Lift the rear side of the Bath Lift by pulling the Rear Handle. Then, lift the Bath Lift away from the bathtub with 2 hands.

9. RECHARGING THE BATTERY PACK



WARNINGS

- When recharging always use a socket away from the bathroom and recharge in a dry environment.
- The Battery Pack containing batteries inside is rechargeable and replaceable, but can only be replaced by qualified personnel. **DO NOT TRY TO REPLACE IT BY YOURSELF.**
- Disposal of the Battery Pack and batteries: Please conform to local regulations. Your supplier will advise on this matter.

Follow the following procedures to recharge your Battery Pack:

1. Remove the Battery Pack from the Bath Lift, as shown in Fig. 5.
2. Plug the charger into the Charging Socket on the Battery Pack. - see Fig. 7.
3. Plug the charger into a suitable power supply.




Fig. 7. Plug in the Charger

Note:

- a. For a new Battery Pack, the first charge shall be no less than 8 hours.
- b. Recharge the Battery Pack as soon as possible after each use. **DO NOT** overcharge it for more than 24 hours.
- c. Disconnect the Battery Pack, from the Bath Lift while not using it. And store it in a cool dry place.
- d. Use only the battery charger which comes with your Bath Lift.
- e. To increase the life of the battery, recharge it at least once a month if you do not use the device for extended period of time.

10. MAINTENANCE

10.1 GENERAL CARING

1. The Bath Lift should be cleaned thoroughly under the shower after each bath.
2. Wash the product with commercial detergents and disinfectants using a cloth or brush.
3. Lightly grease the track and joints by using Vaseline.
4. Follow the following instructions to remove the anti-slip mat and wash it.

10.2 FITTING AND REMOVING THE SOFT MATS

【Fitting/Removing the Seat Mat】

1. Place the Bath Lift in an accessible area. Separate the Backrest and the Seat Assembly.
2. Put the Seat Mat on the seat and let the openings of the seat mat align with the hinge slot on the Seat, as shown in Fig. 8.
3. Insert the Side Flipper into the side pocket of the Seat Mat one at a time.
4. Cap the rear pocket of the Seat Mat on to the Seat Panel.
5. Unfold both the Side Flipper to spread the Seat Mat, as shown in Fig. 9. Let the front end of the Seat Mat hanging freely.
6. To remove the Seat Mat from the seat, reverse the processes above.




Fig. 8. Fitting the Seat Mat onto the seat panel

Let the Front end hanging Freely




Fig. 9. Unfold to spread the Seat Mat

【Fitting/Removing the Backrest Mat】

1. Assemble the Backrest onto the Seat Assembly.
2. Slide the Backrest Mat on to the Backrest Panel. (With the Velcro Straps opened)
4. Center the Backrest Mat in the middle of the Backrest, as shown in Fig. 10.
5. Tighten the Backrest Mat by clipping the Velcro Straps, as shown in Fig. 11.
6. To remove the Backrest Mat from the Backrest, reverse the processes above.




Fig. 10. Slide the Backrest Mat on to the Backrest




Fig. 11. Tighten the Backrest Mat with the Velcro Straps




Fig. 12. Install the Headrest on to the Backrest




Fig. 13. Slide the Headrest Mat on to the Headrest

【Fitting/Removing the Headrest & the Headrest Mat】

1. Press the Headrest Frame into the recesses of the Headrest Panel.
2. Insert the Headrest into the holes on the top of the Backrest and tighten it with 2 hand knobs, as shown in Fig. 12.
3. Slide the Headrest Mat on to the Headrest.
4. To remove the Headrest Mat from the Headrest, reverse the processes above

10.3 CLEANING AND WASHING INSTRUCTIONS

1. Keep your Bath Lift clean and free of residue. Normal washing powder which contains no enzymes or dishwashing liquid are safe to use on this product. Avoid use of abrasive material and solvent based cleaners, such as white spirit, benzene, and thinner.
2. Ensure all areas are thoroughly rinsed before re-using the Bath Lift.
3. The covers are machine washable at 60°C. (140°F)



WARNINGS

- The Bath Lift should not be used without Soft Mat.
- DO NOT IRON the Soft Mat.
- DO NOT TUMBLE DRY the Soft Mat.

11. TROUBLESHOOTING

SYMPTOM	PROBABLE CAUSE	SOLUTION
No motor sound / The Bath Lift does not move	Main Power Loss: 1. The Battery Pack is not installed properly. 2. The battery is flat. 3. The contacts are dirty. 4. The Emergency Stop Button has been pressed. 5. The cable is damaged.	1. Re-install the Battery Pack. 2. Recharge the battery. 3. Clean the contacts. 4. Press the Emergency Stop Button again to Release it. 5. Call dealer to replace the cable.
RED LIGHT illuminated. The Bath Lift can be raised but not lowered	1. Battery Low. Low-voltage protection is activated 2. Internal repairs required.	1. Recharge the battery for at least 8 hours. 2. Call dealer for service.
The Bath Lift stops while being raised.	1. The load is excessive. 2. Internal repairs required.	1. Check to see if the Max. load capacity is exceeded 2. Call dealer for service.
The battery cannot be recharged	1. The battery is faulty. 2. Internal repairs required.	1. Replace the battery. 2. Call dealer for service.
The Bath Lift slips	The suction pads are old or damaged.	Replace the damaged suction pads.


Instructions regarding EMC

1. The Bath Lift is in compliance with EN 60601-1-2.
2. The Bath Lift should not be used adjacent to or stacked with other equipment.
3. The Bath Lift is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies building used for domestic purpose.
4. Mobile RF communications equipment (transmitters) can affect the Bath Lift. However, the user can help prevent electromagnetic interferences by maintaining distance as recommended in the TECHNICAL DESCRIPTION.

Technical Description Regarding Electromagnetic Compatibility

Guidance and manufacturer's declaration - electromagnetic emissions		
The Bath Lift with the model no. A903 is intended for use in the electromagnetic environment specified below. The customer or the user of the above mentioned product should assure that it is used in such environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11 / EN 55011	Group 1	The product use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11 / EN 55011	Class B	The above product is suitable for use in all establishments, including domestic establishments
Harmonics emissions IEC 61000-3-2 / EN 61000-3-2	Class A	and those directly connected to the public low-voltage power supply network that supplies building used for
Voltage fluctuation / flicker emissions IEC 61000-3-2 / EN 61000-3-3	Complies	domestic purpose.

Guidance and manufacturer's declaration - electromagnetic immunity			
The Bath Lift with the model no. A903 is intended for use in the electromagnetic environment specified below. The customer or the user of the above mentioned product should assure that it is used in such environment.			
Immunity test	IEC/EN 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6KV Contact +/- 8KV Air	+/- 6KV Contact +/- 8KV Air	The above product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies building used for domestic purpose.
Electrical fast transient/burst IEC 61000-4-4	+/- 2KV for power supply lines +/- 1KV for input/output lines	+/- 2KV for power supply lines	
Surge IEC 61000-4-5	+/- 1KV for L-N +/- 2KV for line(s) to earth	+/- 1KV for L-N	
Voltage dip and interruptions IEC 61000-4-11	<5 % U _T (>95 % dip in U _T) for 0.5 cycle 40 % U _T (60 % dip in U _T) for 5 cycles <5 % U _T (>95 % dip in U _T) for 5 s	<5 % U _T (>95 % dip in U _T) for 0.5 cycle 40 % U _T (60 % dip in U _T) for 5 cycles <5 % U _T (>95 % dip in U _T) for 5 s	
Power frequency magnetic field IEC 61000-4-8	3 A/m	N/A Not contain component sensitive to magnetic filed	
NOTE: U _T is the AC mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity			
The Bath Lift with the model no. A903 is intended for use in the electromagnetic environment specified below. The customer or the user of the above mentioned product should assure that it is used in such environment.			
Immunity test	IEC/EN 60601-1-2 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 KHz to 80 MHz	3 Vrms,	<p>Portable and mobile RF communications equipment should be used no closer to any part of the above models, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m, 10 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz <p>where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, (NOTE 3) should be less than the compliance level in each frequency range. (NOTE 4)</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1: At 80MHz and 800MHz, the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>NOTE 3: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the above models are used exceeds the applicable RF compliance level above, the above models should be observed to verify normal operation.</p> <p>If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the above models.</p>			
NOTE 4: Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distance between portable and mobile RF communications equipment and the Bath Lift with the model no. A903

The Bath Lift with the model no. A903 is intended for use in an electromagnetic environment in which the radiated RF disturbances are controlled. Mobile RF communications equipment (transmitters) can affect the Bath Lift. However, the customer or the user of the above mentioned product can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Bath Lift as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 KHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitter rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Warranty

Merits Health Products Co., Ltd. warrants the Merits Bath Lift to be free from defective workmanship and materials for a period of 2 years from the date of purchase. The rechargeable battery is warranted for 180 days from date of purchase. Check with your dealer for more details. This warranty is limited to the original purchaser of the Bath Lift. Any defective part or assembly will be repaired or replaced, at the sole discretion and determination of Merits Health Products Co., Ltd., if the unit has been properly operated during the warranty period. Normal maintenance items and disposable components, such as covers, and suction cups are not covered by this warranty. Shipping charges, if any, shall be paid by the purchaser.

Note: There is no other express warranty. Implied warranties, including those of merchantability and fitness for a particular purpose, are limited to the duration of the express limited warranty and to the extent permitted by law and all implied warranties are excluded. This warranty does not cover providing a loaner Bath Lift, compensating for costs incurred for Bath Lift rental, or labor costs incurred in repairing or replacing defective part(s).





Merits Health Products Co., Ltd.

No. 9 Road 36, Taichung Industrial Park, Taichung, Taiwan, R.O.C.
TEL: (886)-4-2359-4985 FAX: (886)-4-2359-4992