

*Jay Protector™*

*User*

*Instruction*

*Manual &*


*Warranty*



## THE JAY PROTECTOR™

The Jay Protector cushion is designed for up to four hours of skin protection when sitting away from your wheelchair. (For all-day sitting with skin protection, we recommend the Jay® J2, Jay Active™, Jay® Care, or Jay® Duo cushions for wheelchairs.)

### Cut Away View of Protector Cushion



### When using your Protector, follow these simple guidelines:

#### I. Distinguish “front” from “rear” of Cushion

The “rear” (A) of your Protector cushion is higher than the front. This should help you sit straighter, with better balance. Plus, the lower front (B) won’t dig into your legs, or interfere with seating heights. When putting the cushion into its cover, always place the lower front in first, so the taller “rear” is toward the zippered end (C).

#### 2. Center yourself carefully on the cushion

When sitting on the cushion, take a few seconds to ensure that you are sitting well-centered with your ischials (seat bones) centered over the Jay Flow™ fluid portion (D) of the cushion. The cushion sides should extend evenly on each side of your buttocks, and the cushion rear should be flush with the back of your buttocks.

#### 3. How to use the sling attachment

Match the Velcro® strips on the cushion with the Velcro strips on the sling and position in place with the taller end of the cushion to the rear. Transfer onto the cushion and make sure you are well-centered. Attach the waistband (E) first and secure tightly, as this is the most important attachment. Now secure each leg strap (F), but not so tightly that it cuts off blood flow or urine drainage if you use a leg bag. (Preferably, the sling should be worn over pants or shorts to prevent strap abrasion against the legs.) Now, try moving around on the floor and check that the cushion does not slide out of place. (If Protector slides out of place, check for proper adjustment of waistband.)

#### 4. Gauge your skin tolerance

You may be able to use the Protector for up to four hours of skin protection. However, extra-sensitive skin or extra-rough use could decrease your skin tolerances. It’s always best to check your skin periodically to gauge your particular tolerances.

#### 5. Easy cleaning

The cushion itself can be scrubbed clean with soap and water. Both the cushion cover and the sling attachment can be machine washed (no bleach) and dried at average temperature settings. Before washing, attach the cushion cover to the sling to prevent the sling’s Velcro strips (G) from abrading attachment straps and other articles in your wash.

## 24 MONTH LIMITED WARRANTY

Your Protector cushion is guaranteed to be free from defects in materials and workmanship for 24 months from date of purchase, given normal use. Should a defect in materials or workmanship occur within 24 months from the original purchase, Sunrise Medical will, at its option, repair or replace it without charge. This warranty does not apply to punctures, tears or burns, nor the cushion's removable cover or sling attachment.

Claims and repairs should be processed through the nearest authorized supplier. Except for express warranties made herein, all other warranties, including implied warranties of merchantability and warranties of fitness for a particular purpose are excluded. There are not warranties which extend beyond the description on the face hereof. Remedies for breach of express warranties herein are limited to repair or replacement of the goods. In no event shall damages for breach of any warranty include any consequential damages or exceed the cost of non-conforming goods sold.



---

**Sunrise Home Healthcare Group  
Mobility Products Division**

7477 East Dry Creek Parkway • Longmont, Colorado • 80503 USA

(303) 218-4500 or (800) 648-8282

Fax: (303) 442-3855 or Toll Free Fax: (800) 300-7502

In Canada (800) 263-3390

Manufactured under one or more of the following U.S. patents and/or U.S. and foreign patents pending.

U.S. Patent #: 4,588,229; 4,842,330

U.K. Patent #: 2,125,695

Japan Patent #: 1,663,369

Germany Patent #: P3337009.5