

AL215 Axis I Value Inside Lift



FEATURES

- Light duty, economical yet effective solution
- Manual rotation with lockout
- Extremely adjustable base and boom arm strap
- Weather resistant hand control
- Curved, offset post allows mobility device to be delivered to passenger or driver's side

SPECIFICATIONS

Lifting Capacity:	250 lbs
Operation:	Power lifting and manual rotation
Adj. Boom Height:	30" - 40"
Adj. Boom Length:	22" - 34"
Adj. Boom Pitch:	0 - 30 degrees
Adj. Rotation Limit:	Yes
Warranty:	3-year transferable

DESCRIPTION

Harmar's AL215 is an efficient and economical solution for transporting smaller, lighter power chairs and scooters. Features include a curved, offset post and an extremely adjustable base and boom arm strap. The AL215 offers powered lifting with manual rotation.