




## MOBILE SERIES

*Mobility, Hi-Low operation and lift accessibility for assisted living.*



# Assured Comfort® Hi-Low Adjustable Beds

### FEATURES

- Includes quiet, electric Hi-Low operation along with head and knee articulation.
- Best solution when bed mobility is necessary.
- Hi-Low feature allows a 10" height adjustment of the foundation deck. Bed can be ordered with or without footboard.
- Leather upholstered, padded assist side rails are available at head and/or foot of bed.
- Cordless and corded remotes are included.
- Choice of mattresses to your specification. Or, use your mattress if it works with our adjustable frame.
- Available in Twin, Full and Queen.
- For additional information about customization, contact one of our representatives

### SALES

**Assured Comfort® Beds, Inc.**

Toll Free: 866.852.2337

3629 Reed Creek Drive / Bassett, VA 24055

[AssuredComfortBed.com](http://AssuredComfortBed.com)

*(Product features and dimensions subject to change without notice.)*

### SPECIFICATIONS

#### HEADBOARD / FOOTBOARD

- Options for wood, metal and upholstered headboards and footboards, connected with matching side panels.

#### ARTICULATION

- Head section moves from zero to about 70 degrees.
- Knee section moves from zero to about 40 degrees.

#### MATTRESS SIZES

- Twin: 38" x 80" / Full: 53" x 80" / Queen: 60" x 80"  
*(Custom bed configurations and sizes are available)*

#### MATTRESSES

- 10" thick, "Comfort Transfer" premium mattress. Other mattress options are available.

#### WEIGHT CAPACITY

- Foundation support varies depending on bed size and mattress. Contact customer service for details.


*Hundreds of combinations available!*



Model shown in down position  
(Shown without mattress)



## Assured Comfort® Hi-Low Adjustable Beds / MOBILE SERIES PROFILE DIMENSIONS



### HEADBOARD DIMENSIONS

All footboards are about 30" in height.  
Exact measurements can be provided once a specific combination is chosen.



Mission Style




Raised Panel

### WOOD HEADBOARDS

| Standard   | Width                             | Height | Depth |
|------------|-----------------------------------|--------|-------|
| Twin       | 41"                               | 41"    | 3"    |
| Full       | 58"                               | 48"    | 3"    |
| Queen      | 64"                               | 48"    | 3"    |
| Split King | (Custom configurations available) |        |       |



Chinese Chippendale



New Modern

### METAL HEADBOARDS

| Standard   | Width                             | Height | Depth |
|------------|-----------------------------------|--------|-------|
| Twin       | 40"                               | 41"    | 2"    |
| Full       | 55"                               | 48"    | 2"    |
| Queen      | 62"                               | 48"    | 2"    |
| Split King | (Custom configurations available) |        |       |



Camelback



Flat Panel

### UPHOLSTERED HEADBOARDS

| Standard   | Width                             | Height | Depth |
|------------|-----------------------------------|--------|-------|
| Twin       | 41"                               | 41"    | 3"    |
| Full       | 56"                               | 47"    | 3"    |
| Queen      | 63"                               | 47"    | 3"    |
| Split King | (Custom configurations available) |        |       |

### CONSIDERATIONS

#### WHAT IS THE PREFERRED TRANSFER HEIGHT?

- With Hi-Low operation, the MOBILE SERIES offers a transfer height range from 24" to 34" when using a standard 10" mattress.

#### DOES THE BED NEED TO ACCOMODATE LIFT?

- Yes. The MOBILE SERIES will accomodate a lift.

#### WHAT IS THE FLOOR SURFACE?

- If the bed will be on a hard floor surface such as wood or tile and needs casters to be moved often, we recommend the MOBILE SERIES.

#### DO YOU NEED SIDE ASSIST RAILS?

- For help getting in and out of bed, consider side assist rails. Height of rail varies with the thickness of the mattress.



### WARRANTY

#### Full One-Year Warranty

Assured Comfort® Adjustable Beds are warranted against defects in workmanship or materials for a period of one year from the date of purchase. Assured Comfort Beds, Inc. will repair or replace any defective part at no cost to the purchaser.

#### Limited Five Year Warranty

During the second through fifth years from the date of purchase, Assured Comfort Beds, Inc. will replace any part found to be defective. Purchaser shall pay all service and shipping costs related to the replacement of the defective parts.

#### Terms and Conditions

This warranty applies to normal use and does not cover any damage caused by excessive wear, abuse, misuse, mishandling or modification of the product.

**QUESTIONS? Call 866.852.2337**