

CATALYST


A & AC





### INNOVATIVE AND PATENTED FRAME DESIGN

Patented frame concept uses intelligent geometry to create a significantly stronger, lighter and better performing chair.



### OFFSET AXLE PLATE

The Catalyst 4 series vertically adjustable axle plate is mounted forward of the back frame. This small change can significantly improve efficiency.



### SWING AWAY FOOTRESTS

Enjoy the simplicity and strength of Catalyst's patented swing-away footrest hangers. The tapered engagement keeps the footrests snug and rattle free. Comes with a lifetime warranty.

## EASY TO TRANSPORT. EASY TO PUSH.

*The Catalyst 4 & 4C are the lightest, most durable and easiest to use K4 folders on the market.*



### SPECIFICATIONS

Transport Weight:	16 lbs.*	7.25 kg*
Seat Width:	14-22 in.	36-56 cm
Seat Depth:	14-20 in.	36-51 cm
Front Seat Height:	13.5-21 in.	34-53 cm
4 Rear Seat Height:	12.5-19 in.	32-48 cm
4C Rear Seat Height:	12.5-20 in.	32-51 cm
Top Frame Tubing Size:	1 in.	2.54 cm
Bottom Frame Tubing Size:	1.125 in.	2.86 cm
Weight Capacity:	250/350 lbs.	113/158 kg

\* 16"x16" Frame, 5"x1" Poly Casters, Height Adjustable T-Arm Receivers. Footrest and Rear Wheels removed.