



Weighs only 16.8 lbs.



## Fly-Lite Aluminum Transport Chair

Item #	Seat Width	Color		UOM
DFL19-BL	19"	Blue Frame and Blue Plaid Upholstery	<span style="color: blue;">●</span>	1/cs
DFL19-BLK	19"	Black Frame and Black/White/Red Plaid Upholstery	<span style="color: black;">●</span>	1/cs
DFL19-RD	19"	Burgundy Frame and Black Upholstery	<span style="color: red;">●</span>	1/cs

- Weighs only 16.8 lbs. without footrests.
- Aluminum frame is lightweight and strong.
- Comes with carry pocket on backrest.
- Comes with seat belt for added safety.
- Attractive plaid style nylon upholstery. (Figure A)
- Comes with handy carry pocket on back. (Figure B)
- Tool-free, height-adjustable swing-away footrests come standard. Optional leg rests available. (Model# STDELR-TF)
- 8" casters in rear. 6" casters in front.
- Back folds down for storage and transport and features deluxe back release. (Figure C)
- Comes with aluminum rear wheel locks.
- Convenient cup holder can easily attach to chair. Sold separately. (Model# STDS1040S) (Figure D)
- Weight Capacity: 300 lbs.

specifications	
Overall Width	22" (open)   9" (closed)
Overall Height	36"
Seat-To-Floor Height	18"
Back Height	18"
Seat Dimensions	19" (W) x 15.5" (D)
Overall Length w/ Riggings	33"
Weight without Footrests	16.8 lbs.
Carton Shipping Weight	23.8 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

accessories		
STDS1040S	Universal Cup Holder	1/ea