

Lightweight Expedition Aluminum Transport Chair

HCPCS: E1038



FIGURE A




FIGURE B

- Weighs only 19 lbs
- Available in blue or red frames
- Folds flat for easy transport (Figure A)
- Handy lever on chair back makes it easy and convenient to fold back down (Figure B)
- Comes with seat belt
- Companion activated wheel locks
- 12" rear "flat-free" wheels

PRODUCTS

Description

EXP19LTBL	Lightweight Expedition Aluminum Transport Chair, EA
EXP19LTRD	Lightweight Expedition Aluminum Transport Chair, EA

SPECIFICATIONS

Description

Back Height	18"
Carton Shipping Weight	24 lbs.
Overall Height	38"
Overall Length w/ Riggings	38"
Overall Width	22" (open) 10.75" (closed)
Seat Dimensions	19" (W) x 16" (D)
Seat To Floor Height	19"
Seat Width	19"
Warranty	Limited Lifetime
Weight	19 lbs.
Weight Capacity	300 lbs

PARTS

Description

STDS1070	Seat Rail Guide, Removable Arm, 1/ea, EA
STDS871	Hand Grip, Sentra, 1/ea, EA
STDS4020	Tipping Sleeve 1/ea, EA
STDS002001	Arm Pad Screw, 1/ea, EA
STDSFLVB	Universal Wheel Chair Vinyl Full Arm Pad, 1/ea, EA
STDS5D0914	Folding Cane Hinge, ATC 5D, 1 ea, EA
STDS5D5505	Wheel Brake, Left, EXP19LT, 1 ea, EA
STDS5D5506	Wheel Brake, Right, EXP19LT, 1 ea, EA
STDS5D5507	Rear Axle Bearings, EXP19LT, 1 pair, PR
STDS5D5508	Rear Hub Cap, EXP19LT, 1 ea, EA

PARTS	Description
STDS5D5509	Rear Axle Bolt with Sleeve,EXP19LT,1ea, EA
STDS5D5510	Rear Wheel,EXP19LT,1ea, EA
STDS5D5513	Front Axle Bearings,EXP19LT,1pair, PR
STDS5D5514	Front Axle,EXP19LT,1ea, EA
STDS5D5515	Front Caster,EXP19LT,1ea, EA
STDS5D5516	Front Fork,Blue,EXP19LT,1ea, EA
STDS5D5517	Front Fork,Red,EXP19LT,1ea, EA
STDS5D5518	Fork Stem Bearings,EXP19LT,1pair, PR
STDS5D1069	Hanger Pins for ATC, SSP, Cruiser,4/pk, BAG
STDS3J0929	Seat Upholstery,19"ATC 3J,5D,9H,1/ea, EA