

Panther Heavy Duty 4-Wheel Scooter



HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Front and rear suspension on 14" pneumatic tires for the smoothest ride possible
- One-touch, infinite tiller adjustment allows the user to easily position the tiller into the ideal driving position
- Extra-bright LED lighting for superb visibility
- Ergonomic delta control handles for ease of steering and throttle control
- Large front plastic basket with carry handle
- 20" to 22" wide Captain's Seat is height adjustable and features a fold-down, reclining backrest with adjustable headrest
- Hazard lights, turn signals and reflectors for superb visibility
- Dual rear view mirrors
- Positioning belt included
- Anti-tip wheels
- Easy, freewheel mode
- USB port for charging your mobile devices

PRODUCTS

Description

PANTHER20CS	Panther Heavy Duty 4-Wheel Scooter, EA
PANTHER22CS	Panther Heavy Duty 4-Wheel Scooter, EA

SPECIFICATIONS

Description

Anti-Tip Wheels	2.5"
Base Weight	145 lbs.
Batteries (pair)	50AH
Battery Charger	5A
Battery Weight (pair)	64 lbs.
Brakes	Electromagnetic
Controller	Dynamic Rhino 2 120A
Estimated Range	25 miles
Freewheel Mode	Yes
Front Wheels	14" x 4"
Ground Clearance	4"
Max Climbing Angle	7.5°
Motor	800W x 4500RPM
Overall Dimensions	56" (L) x 26" (W) x 53" (H)
Rear Wheels	14" x 4"
Seat Dimensions	20-22" (W) x 20" (D) x 22" (H)
Seat Weight	55 lbs.
Top Speed	8 mph

SPECIFICATIONS**Description**

Total Weight (excl. hand control)	264 lbs.
Turning Radius	53"
Warranty	24 Months on Electronics 12 Months on Batteries
Warranty on Frame	Limited Lifetime
Weight Capacity	425 lbs
Weight Capacity (13046 13046-HD)	425 lbs.