

## Bariatric Sentra EC Heavy-Duty Wheelchair

HCPCS: K0007,K0007+K0195 OR E0990



- Weight Capacity: 450 lbs
- Carbon Steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Double embossed vinyl upholstery is durable, lightweight, attractive and easy-to-clean
- Standard dual axle provides easy transition of seat to hemi-level
- Urethane tires mounted on composite wheels provide durability and low maintenance
- Chrome hand rims standard
- 8" front casters, adjustable in three positions
- No maintenance, sealed, precision bearings
- Padded armrests provide additional comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Aluminum foot plates
- Comes with push-to-lock wheel locks

### PRODUCTS

### Description

STD20ECDDAHD-ELR	Bariatric Sentra EC Heavy-Duty Wheelchair, BX
STD20ECDDAHD-SF	Bariatric Sentra EC Heavy-Duty Wheelchair, CS
STD20ECDFAHDA-ELR	Bariatric Sentra EC Heavy-Duty Wheelchair, BX
STD20ECDFAHDA-SF	Bariatric Sentra EC Heavy-Duty Wheelchair, BX
STD22ECDDA-ELR	Bariatric Sentra EC Heavy-Duty Wheelchair, CS
STD22ECDDA-SF	Bariatric Sentra EC Heavy-Duty Wheelchair, CS
STD22ECDFDA-ELR	Bariatric Sentra EC Heavy-Duty Wheelchair, CS
STD22ECDFDA-SF	Bariatric Sentra EC Heavy-Duty Wheelchair, CS
STD24ECDDA-ELR	Bariatric Sentra EC Heavy-Duty Wheelchair, CS
STD24ECDDA-SF	Bariatric Sentra EC Heavy-Duty Wheelchair, CS
STD24ECDFDA-ELR	Bariatric Sentra EC Heavy-Duty Wheelchair, CS
STD24ECDFDA-SF	Bariatric Sentra EC Heavy-Duty Wheelchair, CS

### SPECIFICATIONS

### Description

Back Height	16"   18"
Carton Shipping Weight	62 lbs.   63 lbs.   64 lbs.
Height	36"
Length with Rigging	42"
Product Weight with 20 Seat	56 lbs.
Product Weight with 22 Seat	57 lbs.
Product Weight with 24 Seat	58 lbs.
Seat Depth	18"

SPECIFICATIONS	Description
Seat To Floor Height	17.5" - 19.5"
Seat Width	20"   22"   24"
Warranty	Limited Lifetime
Weight Capacity	450 lbs
Width Closed	12.5"
Width Open	28"   30"   32"
Arm Rest Type (s)	Detachable Desk Arm   Detachable Full Arm

ACCESSORIES	Description
STDS818	Anti-Tippers,rear, SentraRecliner,1 pr/bx, BX
STDS801	Brake Handle Ext.,8", 1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS821	Overhead AntiTheft w/IV Hks,1/bx, BX
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS803	Oxygn Tank Holdr for W/C,1/cs, CS
STDS831	Heel Loops, 1pr/bx, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
WASR	Limb Support, Right,1/cs, CS
WASL	Limb Support,Left,1/cs, CS

PARTS	Description
STDS1S2650	SeatUpholstery20",SentraEC,1S&Y,1/ea, EA
STDS4S262L-N	SideFrame,NewStyle, Sentra EC, Left,1/ea, EA
STDSY2622	CrossFrame22"SentraEC,1ea, EA
STDSY2620	CrossFrame20"SentraEC,1ea, EA
STDS4S262R-N	SideFrame,NewStyle, Sentra EC, Right,1ea, EA
STDSY2624	CrossFrame24"SentraEC,1ea, EA
STDS1016-8	Spacers & Screws, Sentra Handrim,8/pk, BAG
STDS1S3045	Front Caster Fork,Sentra EC,1/ea, EA
STDS1072-8	Seat Upholstery Screws, 5mm, 8/pkg, EA
STDS5D5504	BackUpholsteryScrews10Pack,EXP19LT,1pkg, PKG
STDS00004A	Footplate,Universal,Alum./1ea, EA
STDS1011HD	Caster bearing HD, 11mm X 23mm, 1/ea, EA
STDSADFASPR	Full Arms,Adj.Ht.Rt.1/ea, EA
STDSADFASPL	Full Arm,Adjt.Ht.Left,1/ea, EA
STDSADDASPL	Desk Arm,Adj.Ht.Left,1/ea, EA
STDSADDASPR	Desk Arm,Adj.Ht.Rt,1/ea, EA
STDS004010	Handrim,24"Wheel,Chrome/1ea, EA
STDSFLVB	UniversalWheelChairVinylFullArmPad,1/ea, EA
STDS8M4603	BrakeTipSleeve,Cruiser,SSP, 1/ea, EA
STDS8B	Spring Btn for Arm Rest,1/ea., EA