

WORLD'S LIGHTEST TRANSPORT WHEELCHAIR!

Feather TRANSPORT

Easily visit offices, restaurants and more!



The Feather Transport sets a new standard for transport wheelchairs. **Anyone can lift** this 13 pound wheelchair into a car, truck or SUV. It is super easy to propel and push. The back folds down and the **frame collapses** into a compact package for travel. It has swing-away footrests and is available with or without caregiver handbrakes.

FEATURES:

- just 13 lbs
- 18" wide seat
- folding frame
- optional handbrakes


SPECIFICATIONS

Weight Capacity	300 lbs.		
Product Weight	13 lbs.	Front Wheel Size	6"
Seat Width	18"	Rear Wheel Size	8.5"
Seat Depth	17"	Armrest Type	Fixed
Seat-to-Floor Height	19.5"	Armrest Height	26"
Overall Width	23"	Back Height	16.5"
Overall Length	24"	Upholstery Type	Nylon
Overall Height	36"	Upholstery Color	Black w/ Orange Trim
Folded Depth	24"	Footrests	Swing-away
Folded Height	29"	Warranty on Frame	One Year Limited
Folded Width	15"	Warranty on Wearables	Six Months Limited

Lightweight



Foldable



Locking Brakes



Arrives Assembled

