

NEW VALIANT



Standard 5 Function Package

- ❖ Head Adjustable
- ❖ Foot Adjustable
- ❖ Hi-Low Adjustable
- ❖ Trendelenburg
- ❖ Reverse Trendelenburg
- ❖ 8 Button Sleek Hand Control



Quick Facts

- ❖ Less than 12 inches from deck to the floor
- ❖ Height range- 12 inches
- ❖ 5 Function Model
- ❖ 300 pound weight capacity
- ❖ Robust, modular control box design
- ❖ Tool-free, easy setup split bed
- ❖ Single or Dual hand controls
- ❖ ECU Compatible**
- ❖ Compatible with most Patient Lift systems**

Deluxe 5 Function Package

- ❖ Head Adjustable
- ❖ Foot Adjustable
- ❖ Hi-Low Adjustable
- ❖ Trendelenburg
- ❖ Reverse Trendelenburg
- ❖ 8 Button Sleek Hand Control
- ❖ Deluxe Illuminated Hand Control
- ❖ Single Variable Massage
(operated only by Deluxe Control)
- ❖ Soft Illuminated Floor Light Feature



Available Sizes**

Twin	38"x80"
Twin SL	38"x84"
Full	53"x80"
Full SL	53"x84"
Queen	60"x80"
Also available in Dual King	



Demonstrating easy split technology for easy delivery and setup