



KARMAN HEALTHCARE INC.
 19255 San Jose Ave., City of Industry, CA 91748
 Tel: 1-800-805-2762 / 626-581-2235 Fax: 626 581-2335
 www.karmanhealthcare.com Email: info001@karmanhealthcare.com

6061 T-6 Aircraft-grade Aluminum

Tilt Foldable Series

VIP 515 Revolutionary New Design Tilt Foldable
- Excellent Pressure Distribution



Foldable & Light Weight Tilt

- * The Only Foldable Light Weight Tilt-in-Space Wheelchair
- * Intelligent 0-35 Center-of-Gravity weight shifting
- * **aegis** Microbe Shield permanently bonded Anti-Microbial barrier upholstery protects from odor, staining and deterioration from bacteria, fungus and other micro-organisms.
- * Exclusive Perfectly-Matched Gas Cylinders
- * Bio-technology: Paint coating uses special Anti-Bacterial materials.
- * Special New Design Ergonomic Pillow
- * Flip-Back & Height Adjustable Armrests
- * Swing-Away Footrests & Anti-tippers
- * Weight Capacity: 253 lbs.
- * **Weight Only 34 lbs** (w/o Footrests & Headrest)



The Lightest Tilt Wheelchair in The World

▲ Pressure is centralized on normal wheelchairs	▲ Pressure is distributed when using tilt-in-space	▲ Center-of-gravity weight shifting tilt-in-space	▲ Height adjustable/ flip-back & armrests	▲ Inward/outward swing-away footrests