

Ventura 3 DLX

HCPCS:



FIGURE A




FIGURE B




FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue (Figure A)
- 33 AH batteries for extended range
- Floorboard with 23" of foot space
- Ergonomic throttle control
- Height-adjustable swivel seat features fold-down backrest and adjustable padded armrests
- Seat position can be adjusted forward or backward
- Rear suspension and Anti-tip wheels
- Programmable controller
- Seat back reclines and features adjustable headrest
- Available with 18" or 20" width Captains Seat
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust Delta tiller
- Headlight, Turn signals/hazard lights and Rearview mirrors
- Easy free wheel operation

PRODUCTS

PRODUCTS	Description
VENTURA320CS	VENTURA 3-WHL DELUXE 20"CAPT.SEAT,RED/BL

SPECIFICATIONS

SPECIFICATIONS	Description
Base Weight	88 lbs.
Batteries (pair)	12V x 33AH
Battery Weight (pair)	46 lbs.
Brakes	Electromagnetic
Charger	4A Offboard
Climbing Angle	7.5°
Dimensions	48" (L) x 25" (W)
Floor to Seat Height	20.5" - 23.5"
Front Wheels	10" x 3" Flat-free
Ground Clearance	2.5"
HCPCS	N/A
Heaviest Piece Weight	60 lbs.
Maximum Range	16 miles
Motor	24V x 350W x 5100 rpm

SPECIFICATIONS	Description
Rear Wheels	10" x 3" Flat-free
Seat Dimensions	18" (W) x 18" (D) 20" (W) x 18" (D)
Seat Type	Captain's
Seat Weight	40 lbs.
Top Speed	6 mph
Total Weight (excl. hand control)	174 lbs.
Turning Radius	53"
Warranty on Batteries	12 Months
Warranty on Electronic Controller/Drive Train Components	24 Months
Warranty on Frame	Lifetime
Weight Capacity	400 lbs
Weight Capacity (13046 13046-HD)	400 lbs.

ACCESSORIES	Description
SF8020	Rear Basket Set, LType, 1/ea, EA
SF8050	Sun Shade for Scooter,1/ea, EA
SF8010	Oxygen Holder Set, L Type,1/ea, EA
SF8030	Crutch Holder Set, SType, 1/ea, EA

PARTS	Description
VENT-74	Seat Plate, Capt. Seat Ventura 3&4, 1/ea, EA
VENT-73	Seat Plate, Folding Seat Ventura 3&4, EA
SEC2010N	HeadRestNew18"SunfireEC,1ea, EA
VENT-65	Seat Slide Set, Ventura 3 & 4, 2/pkg, PKG
GT8072W32	Hazard Button,GladGT,Ventura,1/ea, EA
VENT-53	Mirror,Right, Ventura,Panther,1/ea, EA
VENT-54	Mirror,Left,Ventura,Panther,1/ea, EA
GT8072W35	FrontSignalLeftGladGT,Ventura,1/ea, EA
GT8072W36	FrontSignalRightGladGT,Ventura,1/ea, EA
GT8072W53	DirectionBulbGladGT,Ventura,1/ea, EA
VENT-55	Basket Holder, Ventura, 1/ea, EA
VENT-61	Front Decoration, Red, Ventura3&4, 1/ea, EA
VENT-62	Front Decoration, Blue, Ventura3&4, 1/ea, EA
VENT-46	Coupler, Ventura 3, 1/ea, EA
VENT-49	FrontWheelAssemblyVentura3,1/ea, EA
VENT-45	FrontShroudVentura3,1/ea, EA
VENT-50	FrontWheelAxle,Ventura 3,1/ea, EA
VENT-48	FrontForkAssemblyVentura3,1/ea, EA
VENT-58	FrontWheelW/PneumaticTire,Ventura3,1/ea, EA
VENT-59	OuterTireFront&RearWheelPneumTireVent3&4, EA