

4-Wheel

Model Number	S710LXW
Weight Capacity	400 lbs.
Maximum Speed ¹	Up to 8 mph
Ground Clearance ³	2.2" at motor, 3.15" at mid-frame
Turning Radius ³	61.25"
Overall Length ³	47"
Overall Width ³	25.62"
Seat-to-Ground Height	23" - 25"
Seat-to-Deck Height	16.5" -18.5"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Suspension	Front and rear
Range Per Charge ^{1,2}	Up to 19.7 miles at 200 lbs. Up to 13 miles at 400 lbs.
Total Weight Without Batteries ⁵	160 lbs.
Heaviest Piece When Disassembled ³	77 lbs. (front section)
Standard Seating	Type: Reclining high-back Weight: 43 lbs. Material: Black Vinyl Dimensions: 18" width, 18" depth
Drive System	Rear-wheel drive, sealed transaxle with 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements ^{4,6}	Size: (2) 50 AH Weight: 30.2 lbs. each
Battery Charger	5 amp, off-board
Warranty	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-months

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC
OR DOES NOT MEET CRITERIA

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Includes standard seat weight.
- 6 AGM or gel cell type required.



(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

VICLXSPORT_020119

VICTORY[®]LX SPORT

ENJOY A SMOOTH AND SPORTY RIDE



LIVE YOUR BEST[®]

pridemobility.com |    

SMOOTH, COMFORTABLE RIDE

18" x 18"
Vinyl reclining
high-back seat
with sliders,
height adjustable
headrest and
flip-up, height
and width-
adjustable
armrests

Puddle light
under tiller

Standard
front basket

Easy-to-grip
tiller adjustment
knob

Full LED
lighting
package
includes
curb light,
rear lights
and
directional
signals

Storage
compartment
under the tiller

CTS Suspension includes adjustable
shocks for greater comfort



Tiller Console



CTS Suspension

Maintain your active lifestyle with the **Victory® LX Sport**. Equipped with **Comfort-Trac Suspension (CTS)**, the Victory LX Sport easily handles more varied terrain while providing a smooth and comfortable ride. Enjoy convenient storage and a courtesy light under the tiller. A long-lasting LED headlight ensures optimal pathway illumination. Experience great features and a sporty style with the Victory LX Sport.

- CTS Suspension includes adjustable shocks for greater comfort
- Exclusive one-hand feather touch disassembly
- Charger port located on the tiller lets you conveniently charge your scooter
- Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand and rest your wrist
- Easy-to-read LED voltmeter
- Rearview mirror

Available Colors



True Blue



Candy Apple Red

Optional Seating

Limited recline, black vinyl, high back with height adjustable headrest, flip-up height and width adjustable armrests

- 22" x 18"
- 20" x 20"

Accessories

- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Additional rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Lap belts: 50", 60" and 70"

Military Patches

