

	EW 14
Warranty	3 Year Limited Warranty
Colors	Red, Blue
Weight Capacity	350 lbs
Fully Assembled	Yes
3 or 4 Wheel	4 Wheels
FWD or RWD	Rear Wheel Drive
Max Passengers	2
Airline Approved	No
Seat Dimensions (LxWxH)	14" x 23"
Type of Seat	Extended Saddle Seat
Type of Battery	Deep Cycle Sealed Lead Acid
Overall length	70"
Overall width	31"
Overall height	44"
Armrests	No
Scooter Weight	250 lbs
Shipping Method	LTL Freight
Top speed (Up To)	15
Distance (Up To)	40
Light(s)	Turn Signals, Brake Lights, High/Low Beam, Headlight, Rear Light
Turning Radius	126"
Braking Type	Hand Brakes
Left Hand Control Option	Yes
Horn	Yes
Incline capability	Up to 12 Degrees
Ground Clearance	6.5"
Suspension	Front and Rear
Type of Tires	Tubeless
Battery Charger	48 volt offboard
Speed Control	Low-Med-High
Basket	No
Forward/Reverse	Yes
Throttle Type	Twist
Alarm	Yes
Front Tire Size	15" x 3"
Rear Tire Size	15" x 3"
Cell Phone Holder	No
USB Port	No
Voltage	48 volt offboard
Motor Wattage	500 Watt
Motor Type	Brushless Motor with Differential
Storage	Lockable Rear and Under Seat
Rear View Mirror	yes
Extras	Bag hook, rear reflectors