

	<b><i>EW-54</i></b>
<b><i>Warranty</i></b>	3 Year Limited Warranty
<b><i>Colors</i></b>	Red, orange, green
<b><i>Weight Capacity</i></b>	500 lbs
<b><i>Fully Assembled</i></b>	Yes
<b><i>3 or 4 Wheel</i></b>	4 Wheels
<b><i>FWD or RWD</i></b>	Rear Wheel Drive
<b><i>Max Passengers</i></b>	1
<b><i>Airline Approved</i></b>	No
<b><i>Seat Dimensions (LxWxH)</i></b>	22.5" x 17.5" x 23"
<b><i>Type of Seat</i></b>	Dual adjustment comfort seat
<b><i>Type of Battery</i></b>	Deep Cycle Sealed Lead Acid
<b><i>Overall length</i></b>	72"
<b><i>Overall width</i></b>	27"
<b><i>Overall height</i></b>	65"
<b><i>Armrests</i></b>	Yes
<b><i>Scooter Weight</i></b>	333 lbs
<b><i>Shipping Method</i></b>	LTL Freight
<b><i>Top speed (Up To)</i></b>	15 mph
<b><i>Distance (Up To)</i></b>	43 miles
<b><i>Light(s)</i></b>	Front headlight with high/low beam, rear running/brake lights, turn signals
<b><i>Turning Radius</i></b>	90"
<b><i>Braking Type</i></b>	Foot Brake
<b><i>Left Hand Control Option</i></b>	No
<b><i>Horn</i></b>	Yes
<b><i>Incline capability</i></b>	Up to 12 Degrees
<b><i>Ground Clearance</i></b>	5"
<b><i>Suspension</i></b>	Front and rear shock absorbers
<b><i>Type of Tires</i></b>	Tubeless
<b><i>Battery Charger</i></b>	60 volt offboard
<b><i>Speed Control</i></b>	Low-Med-High
<b><i>Basket</i></b>	Rear Cargo Box
<b><i>Forward/Reverse</i></b>	Yes
<b><i>Throttle Type</i></b>	Twist
<b><i>Alarm</i></b>	Yes
<b><i>Front Tire Size</i></b>	13" x 3"
<b><i>Rear Tire Size</i></b>	15" x 3"
<b><i>Cell Phone Holder</i></b>	No
<b><i>USB Port</i></b>	No
<b><i>Voltage</i></b>	60 volts
<b><i>Motor Wattage</i></b>	700 watt
<b><i>Motor Type</i></b>	Brushless Motor with Differential
<b><i>Storage</i></b>	Lockable Rear and Under Seat
<b><i>Rear View Mirror</i></b>	Yes
<b><i>Extras</i></b>	radio, cruise control, windshield with canopy