

## Cruiser III Wheelchair

HCPCS: K0003,K0003



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E

- Built-in seat rail extensions and extendable upholstery easily adjust seat depth from 16" to 18" (Figure A - C)
- Weighs under 36 lbs (excluding front riggings)
- Carbon steel frame with silver vein finish
- Removable flip-back arms allow for easy transfer (Figure D)
- New frame style eliminates seat guides and allows for custom back inserts and accessories
- Nylon upholstery is durable, lightweight, attractive and easy-to-clean
- Composite, Mag-style wheels are lightweight and maintenance free
- 8" front casters are adjustable in three positions
- Padded armrests provide additional comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings (Figure E)
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Dual axle provides easy transition of seat height to hemi-level
- Comes with push-to-lock wheel locks

### PRODUCTS

### Description

PRODUCTS	Description
K316ADDA-ELR	Cruiser III, 16"Lt. Wt, FlipBkDskArm, 1/cs
K316ADDA-SF	Cruiser III, 16"Lt. Wt, FlipBk, RemvArm, 1/cs
K316ADFA-ELR	Cruiser III, 16"Lt. Wt, FlipBkFullArm, 1/c
K316ADFA-SF	Cruiser III, 16"Lt. Wt, FlipBkFullArm, 1/c
K316DDA-ELR	Cruiser III, 16"Lt. Wt, FlipBk, DeskArm, 1/cs
K316DDA-SF	Cruiser III, 16"Lt. Wt, FlipBk, RemArm, 1/cs
K316DFA-ELR	Cruiser III, 16"Lt. Wt, FlipBk, FullArm, 1/cs
K316DFA-SF	Cruiser III, 16"Lt. Wt, FlipBk, FullArm, 1/cs
K318ADDA-ELR	Cruiser III, 18"Lt. Wt, FlipBk, DeskArm, 1/cs
K318ADDA-SF	Cruiser III, 18" Wt Flip Back
K318ADFA-ELR	Cruiser III, 18"Lt. Wt, FlipBkFullArm, 1/c
K318ADFA-SF	Cruiser III, 18"Lt. Wt, FlipBkFullArm, 1/c
K318DDA-ELR	Cruiser III, 18Lt. Wt, FlipBk, ElvLgRst, 1/cs
K318DDA-SF	Cruiser III, 18" Lt. Wt, FlipBk, SwgFt, 1/cs
K318DFA-ELR	Cruiser III, 18"Lt. Wt, FlipBk, FullArm, 1/cs
K318DFA-SF	Cruiser III, 18"Lt. Wt, FlipBk, 1/cs
K320ADDA-ELR	Cruiser III, 20"Lt. Wt, FlipBk, 1/cs
K320ADDA-SF	Cruiser III, 20"Lt. Wt, FlipBk, RemArm, 1/cs
K320ADFA-ELR	Cruiser III, 20"Lt. Wt, FlipBkFullArm, 1/c
K320ADFA-SF	Cruiser III, 20"Lt. Wt, FlipBkFullArm, 1/c

PRODUCTS	Description
K320DDA-ELR	Cruiser III,20"Lt.Wt,FlipBk,DskArm,1/cs
K320DDA-SF	CruiserIII,20"Lt.Wt,FlipBk,RemDskArm,1/cs
K320DFA-ELR	Cruiser III,20"Lt.Wt,FlipBk,FullArm,1/cs
K320DFA-SF	Cruiser III,20"Lt.Wt,FlipBk,FullArm,1/cs

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	35 lbs 36 lbs 38 lbs
Weight Capacity	300 lbs 350 lbs
Height	35"
Seat To Floor Height	17.5" - 19.5"
Carton Shipping Weight	43 lbs or less
Back Height	16"

ACCESSORIES	Description
STDS860	Forwrd Stablzrs/Univ.Fit , 1pr/bx, BX
STDS4S2432	Caster with 11mm Bearings, 1/ea, EA
STDS819	Anti Tippers,Cruiser,1pr/bx, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS801	Brake Handle Ext.,8", 1/ea, EA
WASR	Limb Support, Right,1/cs, CS
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDSSHK-K3	Super Hemi Kit, Cruiser 3 & X4, 1/ea, EA
STDSK234-PTL	Pull-To-Lock Brake, Cruiser, Left, 1/ea, EA
STDS803	Oxygn Tank Holdr for W/C,1/cs, CS

PARTS	Description
STDS8M462R	BrakeAsm.DetachableArmRT,SS2,1/ea, EA
STDS8M4518	SeatUpholstery18",SS1,1/ea, EA
STDS8M4651	ForkStemCap,SS2,BS,Cruiser,1/ea, EA
STDS3J2450	Footrest Bumper,Cruiser3,1/ea, EA
STDS8M4601	ForkStemBearing,SSP2,Cruiser,1/ea, EA
STDS3J2424	Wheel24in Blk.Handrl,Cruiser3,SSP,BS,1/e, EA
STDS3J2400	Handrim,Cruiser3,1/ea, EA
STDS3J2413	Rear Axle Bearing,Cruiser3,1/ea, EA
STDS3J2421	HubCap,RearWheel,Cruiser3,SSP,BS,1/ea, EA
STDS3J2410R	Footplate&HeelLoop,Right,Cruiser,SSP,CS, EA
STDS3J2432	Caster with 8mm Bearings, 1/ea, EA
STDS001000	Axle Cap 8mm, 1ea, EA
STDS3J2408	BackUphol.18"Cruiser3,Blk.1/ea, EA
STDS3J24SF	Footrest Set,Cruiser3,SSP, 1pr/set, EA
STDS8M462L	BrakeAsm.DetachableArmLt,SS2,1/ea, EA
LK3JELR	ELR Set,Black,1Pair eq EA, PR
STDS831	Heel Loops, 1pr, BX
STDSDSKVB	UniversalWheelChairVinylDeskArmPad,1/ea, EA
STDSFLVB	UniversalWheelChairVinylFullArmPad,1/ea, EA

PARTS

Description

STDS004609-9

Hand Rim Screw Set of 9,1Bag, BAG