

AL003

TILT-N-TOTE WHEELCHAIR LIFT



INNOVATIVE AND LIGHTWEIGHT WHEELCHAIR LIFT

Lightweight, non-powered solution for manual folding wheelchairs. Innovative tilting joint helps the wheelchair roll on/off easily and the adjustable hold-down arm safely secures it. Compact design folds down for tail gate/lift gate access, folds up when not needed.

AL003 PRODUCT SPECIFICATIONS

LIFTING CAPACITY	100 lb	ACCOMMODATES	Manual folding wheelchairs
INSTALLED WEIGHT	25 lb	PLATFORM SIZE	14" x 30"
HITCH CLASS	I or greater	POWERED	Manual
WARRANTY	3-year	INCLUDED	Ball mount
AVAILABLE OPTION(S)	Manual folding wheelchairs	SECUREMENT OPTION(S)	- Manual hold-down arm - Wheel cradle



Harmar.com
800.833.0478

02JULY2020 | MKT-000130 REV B

