

Buzzaround EX

Bumper Conversion Kit &
Battery Docking Station



BUZZ™
AROUND EX
EXTREME



 **GOLDEN**



Bumper Conversion Kit & Battery Docking Station

Features

- **Control Panel:**
Easy to see LED battery gauge; easy accessible speed control; key switch and horn button.
- **Delta Tiller handle:**
Offers comfort and easy right or left hand use.
- **Headlight On/Off Switch:**
Located on the right side of the tiller for easy access.
- **Front Basket:**
Easily store your items.
- **LED Headlight:**
Angle-adjustable headlight, easier to see and be seen.
- **Tiller Adjustment and Lock Knobs:**
Lower and lock the tiller in place for more compact storage without removing the battery packs.
- **Optional Battery Docking Station**



Full Suspension:
Front and rear springs offer an all-day comfort ride.




Storage:
Offers vertical storage with a compact footprint of 22 inches wide by 23 inches deep.



Front Bumper and Side Rails: Protect the scooter completely from impacts.



Large Rear Bumper:
Offers added protection from mishaps. The caster wheels offer ease of storage by allowing the scooter to be rolled into place while standing vertically.



Camouflage & other optional color shrouds available.
*Upcharge applies.



Battery Docking Station for Models GB118 & GB148

The new Battery Docking Station is portable and charges the interlocking U1 batteries. Perfect for keeping a set of charged batteries on hand.

Item Description	Item Number
Camouflage Shroud Kit	MBB-BZ-SHRD-CMO
Pink Shroud Kit	MBB-BZ-SHRD-PNK
Silver Shroud Kit	MBB-BZ-SHRD-SIL
Yellow Shroud Kit	MBB-BZ-SHRD-YEL
3-Wheel Bumper Kit	MBB-118R-KIT
4-Wheel Bumper Kit	MBB-148R-KIT
Battery Docking Station	MBA-BAT-DOCK