



# Invacare® IVC™ 9000 Jymni™ Wheelchair

## Thank you for using Invacare eForms.

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

For additional information please visit [www.invacare.com/eForms](http://www.invacare.com/eForms)



## Adobe Acrobat Reader DC

\* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

### Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

### Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

### Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

**800-678-4682**

SUBMIT

CLEAR FORM

PRINT FORM



Yes, you can.



# Account Information

Request Type: Quote                      Order

Date: \_\_\_\_\_

Account #: \_\_\_\_\_

Company: \_\_\_\_\_

SHIP TO

Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Country: \_\_\_\_\_

Purchase Order #: \_\_\_\_\_

CONTACT

Name: \_\_\_\_\_

Back Up Contact: \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

Comments: \_\_\_\_\_

## Special Note

The Invacare IVC 9000 Jymni pediatric wheelchair offers seat widths and depths to accommodate the pediatric client for short-term rentals.

## FRAMES

### FRAME TYPE

9JYLT 9000 Jymni Light Weight Chair .....\$1,252.00 E1236

### FRAME WIDTH AND DEPTH

WD00 10" x 10" .....No Charge  
Not available with 17" seat-to-floor height, unless chair is a recliner.

WD20 12" x 10" .....No Charge  
Not available with 17" seat-to-floor height, unless chair is a recliner.

WD22 12" x 12" .....No Charge  
Not available with 17" seat-to-floor height, unless chair is a recliner.

WD42 14" x 12" .....No Charge  
Not available with 17" seat-to-floor height, unless chair is a recliner.

WD44 14" x 14" .....No Charge

WD64 16" x 14" .....No Charge  
16" Seat width changes base code to K0003

### SEAT TO FLOOR HEIGHT

15STF 15" Seat to Floor Height .....No Charge  
Only available with 20" rear wheels.

16STF 16" Seat to Floor Height .....Standard  
Only available with 20" rear wheels.

17STF 17" Seat to Floor Height .....No Charge  
Only available with 22" rear wheels.



Yes, you can.



### ARM TYPE

All above selections come standard with fixed height arms 7 1/2", adjustable height arms 8" - 10" in 1" increments. Removable flip-back arms.

- 18JY Full Length Fixed Height Conventional .....  
.....No Charge  
Recliner and angle back must use conventional arms. Arms do not flip back on the recliner.
- 19JY Full Length Adjustable Ht Conventional .....  
.....\$139.23 E0973  
Recliner and angle back must use conventional arms. Arms do not flip back on the recliner.
- 28JY Desk Length Fixed Height Conventional.....  
.....No Charge  
Recliner and angle back must use conventional arms. Arms do not flip back on the recliner.

- 29JY Desk Length Adjustable Ht Conventional . \$139.23  
Recliner and angle back must use conventional arms. Arms do not flip back on the recliner.
- 78JY Full Length Fixed Height Space-Saver .....  
.....No Charge
- 79JY Full Length Adjustable Height Space-Saver .....  
..... \$139.23
- 88JY Desk Length Fixed Height Space-Saver .....  
.....No Charge
- 89JY Desk Length Adjustable Height Space-Saver .....  
..... \$139.23

### BACK TYPE

All below selections come standard with fixed height arms 7.5", adjustable height arms 8" - 10" in 1" increments. Removable flip-back arms.

- ADJBH Adjustable Ht Back with 10° Bend ..... Standard  
Adjustable back height 14"-15". Factory set at 14".
- R Recliner Frame ..... \$417.69 E1014  
Recliner is only available in U240 black nylon upholstery.  
Recliner back height 16" plus 11" with headrest extension.  
This selection comes standard with rear anti-tippers.  
This selection must be made when WD42 and I7STF are selected.

- ANGLE Adjustable Angle Back with 10° Bend .... \$107.10  
Recliner and angle back must use conventional arms.  
Arms do not flip back on recliner.  
This selection comes standard with rear anti-tippers.

## POSITIONING STRAP OPTIONS

### SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap ...\$53.55 E0978
- 1515 Push Button Style Seat Positioning Strap .....  
.....\$53.55 E0978

- 1321BK Hook and Loop Seat Positioning Strap .....  
..... \$40.00 E0978

### CHEST POSITIONING STRAPS

- 9515 Push Button Style Chest Positioning Strap .....  
..... \$107.10 E0960

## CASTER AND WHEEL OPTIONS

### CASTER OPTIONS

- 1687 6" x 1" Dark Grey Comp Urethane Casters .....  
.....Standard
- 1686 6" x 1" Composite Urethane Casters .....No Charge

- AT902 6" Composite Urethane 3 Spoke .....No Charge



Yes, you can.



### WHEEL OPTIONS

- 1763BK 20" Composite Pneumatic Tires .....\$125.31 E2211+E2212
- 4763BK 20" Composite Pneumatic/Flat Free Inserts ..... \$248.47 E2211+E2213
- U1763BK 20" Composite Urethane Tires ..... Standard

- 1764BK 22" Composite Pneumatic Tires .....\$125.31 E2211+E2212
- 4764BK 22" Composite Pneumatic/Flat Free Inserts .....\$248.47 E2211+E2223
- U1764BK 22" Composite Urethane Tires ..... \$121.02

### HANDRIM OPTIONS

- BHDRM Black Handrims .....Standard

## ACCESSORIES

### SEAT ACCESSORIES

- L8651 Omit Seat Upholstery ..... Less \$10.00

- 9120 Solid Seat Insert .....\$104.96

### BACK ACCESSORIES

- 1483JY Stroller Handle .....\$97.46  
Not available on Recliner or Angle back type.
- 9130 Solid Back Insert .....\$114.60  
Only available in 14" back height.

- L8601 Omit Back Upholstery ..... Less \$10.00
- 1483R Headrest w/Push Handle Extensions ..... \$101.75 E0966

### MISCELLANEOUS ACCESSORIES

- A1576 Wheel Lock Extensions .....\$41.77 E0961
- 9691 Push to Lock Wheel Locks .....No Charge  
Pull-to-Lock wheel locks are standard.
- 1577QR Quick Release Axles .....\$92.11
- 1487 IV/Oxygen Holder Package Assembly .....\$169.22 E2208+K0105  
Not available on Recliner or Angle back type.

- 1488 Oxygen Holder Package Assembly ..... \$161.72 E2208  
Not available on Recliner or Angle back type.
- 1360 Adjustable Rear Anti-Tipper ..... \$87.82 E0971  
Anti-Tippers for 22" rear wheels.
- 1360B Adjustable Rear Anti-Tippers- 20" Wheels ..... \$87.82 E0971  
Anti-Tippers for 20" rear wheels.

## COLOR OPTIONS

### FRAME FINISH

- 104P Electric Teal .....No Charge
- 105P Midnight Blue .....No Charge
- 08P Deep Purple .....No Charge
- 24P Wet Black .....No Charge
- 34P Satin (Flat) Black .....No Charge

- 61P Electric Red .....No Charge
- 62P Electric Blue .....No Charge
- P Electric Purple .....No Charge
- 71P Silver Vein .....No Charge



### UPHOLSTERY COLORS

U240 Black Nylon Upholstery .....Standard  
Recliner is only available with U240 black nylon upholstery. Recliner back height 16" plus 11" with headrest extension.

U550 Midnight Blue Vinyl Upholstery .....No Charge

U67 Black Vinyl Upholstery .....No Charge

## FRONT RIGGING OPTIONS

### FRONT RIGGINGS

J932A 70° Footrests with Aluminum Footplates .....  
.....\$202.52

J932A-14 80° Footrests with Alum Ftplts - 12" and  
14" deep .....\$202.52  
Not Available on 10" deep frames.

J933A 90° Footrests with Aluminum Footplates .....  
.....\$202.52

J934A Elevating Legrests with Aluminum Footplates .....  
..... \$287.79 K0195  
Knee-to-Heel range is required with J934A

J934A-14 Elevating Legrests with Alum Ftplts - 14" deep  
..... \$270.00 K0195  
Not available on 10" and 12" deep frames.

### FOOTREST/LEGREST ACCESSORIES

1820 Non-Removable Device ..... \$12.79

1346 Heel Loops ..... \$41.57 E0951  
Standard on Swingaway Footrests.

### KNEE-TO-HEEL RANGE

Knee-to-Heel Range must be selected if J934A is chosen.

SHEXT 7" - 10" (Short Extended) .....No Charge

LGEXT 11" - 14" (Long Extended) .....No Charge

MDEXT 9" - 12" (Medium Extended) .....No Charge

### Seat to Floor Matrix 6" Casters

SEAT TO FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
15"	Top	20"	Top
16"	Middle	20"	Bottom
17"	Bottom	22"	Bottom

#### Weight Limitation 125 Pounds

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.



## Adobe Acrobat Reader DC

\* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

### Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

### Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

### Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

**800-678-4682**

SUBMIT

CLEAR FORM

PRINT FORM