

SB2 Order Form

For assistance call (877) 769-2833



Company	<input type="text"/>	Name	<input type="text"/>
Address 1	<input type="text"/>	Phone	<input type="text"/>
Address 2	<input type="text"/>	Email *	<input type="text"/>
City	<input type="text"/>	PO Number	<input type="text"/>
State	<input type="text"/>	Reference	<input type="text"/>
Postal Code	<input type="text"/>	Notes	<input type="text"/>

SolutionBased BathLyft (White)

The BathLyft is a battery-operated bath chair that lowers and lifts the user in and out of the bathtub, reducing the strain on the user and preventing the risk of slip and falls. The stable frame and base plate ensure the user's stability and safety, while the padding and optimized dimensions ensure a comfortable bathing experience. The Bathlyft's superior design allows for a backrest recline of up to 52° and lowers into the tub with a minimum sitting position of 2.5" from the bottom. The side flaps move into and out of place to create a transfer surface, allowing the user even more peace of mind. The water-tight hand control will provide the user with approximately 35 lifts on a single charge. The built-in safety switch will prevent use if the battery is too low to perform a complete lift. The two-piece construction makes this lift easy to assemble, disassemble, store, and transport. Assembly is tool-free and ensures easy installation and operation. Available in Blue, Orange, Pink, and White.

The BathLyft weighs 26 lbs with a weight capacity of 300 lbs.

1. Primary	Description	MSRP	Quantity
BLFT-W	SolutionBased BathLyft (White)	\$465.70	1

2. Other	Description	MSRP	Quantity
BLFT-1	SolutionBased BathLyft remote	\$120.00	
BLFT-3	SolutionBased BathLyft remote charger	\$85.00	
BLFT2-PT6	SolutionBased Upholstery White Bathlyft (Second Gen)	\$45.00	