


CSL500 HELIX

AS UNIQUE AS YOU AND YOUR HOME



ONE OF A KIND

A Helix Curved Stair Lift is one of a kind. Harmar designs are crafted to fit your needs, staircase, and decor.

The Helix is built to last. Every staircase that is fitted with a Helix is measured with a calibrated camera kit for digital accuracy. Constructed with highly durable steel, custom precision-bent rails, and finished by experienced craftsmen, this true curved rail is highly engineered to give you an extraordinarily smooth ride.


Your unique path to independence will allow you to enjoy time with your family and your home on every level giving you independence, freedom, and the quality of life you deserve.

TIGHTER CURVES


The Helix's narrow profile mounts closer to the wall, hugging the stairway with tighter bends and radius curves.

This means more room for you, and a comfortable smooth ride.


CSL500 MODELS




CSL590: 90° first turn



CSL580: 180° first turn



CSL570: Sweep



CSL560: Spiral



My Helix was made just for me. Most importantly, it was made for my independence, my confidence, my quality of life. It provides a solid and stable ride that has given me full access to my home.

CSL500 PRODUCT SPECIFICATIONS

WEIGHT CAPACITY	350 lb	SPEED	20 FPM
TRACK LENGTH	- 20' standard - Extra rail up to 164'	INCLINE LIMIT(S)	up to 60°
SEAT SIZE	19"W x 15"D	FOOTREST SIZE	11"W x 12"L
FOLDED WIDTH	15.25"	DRIVE	Steel rack and pinion
CERTIFIED	ETL listed	WARRANTY	- Limited lifetime - 2-year parts - 1-year battery
SAFETY FEATURE(S)	- Direction limit switches - Final limit switch - Footrest obstruction sensors - Seat swivel interlock - Constant pressure controls - Retractable seat belt - Chassis obstruction sensors	OPTION(S)	- Emergency stop - Key switch - Additional call/send - Additional charging station - Black or beige rail, chassis, upholstery