

RevVA2

VISION CF

FULL SIZE



merits

8/12/2020

FEATURES

- 16" width with 16"-18" depth adjustability
- 18" or 20" seat width with 18" - 20" depth adjustability
- Premium, semi-reclining ergonomic seat; up to 135°
- 9" knobby drive tires for superior traction
- Dynamic LiNX controller for convenient Bluetooth programming and diagnostics
- True U1 batteries with 34 Ah

SPECIFICATIONS

Weight Capacity	300 lbs
Max Speed	5 mph
Per-Range Charge	18 miles
Batteries	U1-34Ah
Ground Clearance	2.4"
Turning Radius	21"
Overall Length	33"
Overall Width	24"
Drive Wheels	9" foam filled tires
Rear Casters	6" foam filled tires
Seat Width	16", 18" or 20"
Seat Depth	16", 17", 18", 19" or 20"
Back Height	18" or 20"
Controller	Dynamic LiNX
Warranty	18 months : Motor & Electronics 5 Year/ Limited Lifetime Frame

OPTIONS

- Swing-away joystick
- USB charger
- Oxygen tank holder
- Foldable rear basket
- Cane & crutch holder
- Walker holder
- Weather cover
- Elevating legrests

MODELS

Vision CF - G2, Standard, Captain	P322-2S C
Vision CF - G2, Standard, Pan	P322-2S S
Vision CF - G2, Portable, Captain	P322-2SP C
Vision CF - G2, Portable, Pan	P322-2SP S
Vision CF - G1, Standard, Captain	P322-1S C
Vision CF - G1, Standard, Pan	P322-1S S



Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

Phone: 1-800-963-7487 Fax: 239-772-32

merits