

SOPHISTICATED AND STYLISH

VICTORY^{® 1}

Versatility knows no bounds with Victory[®] Scooters.¹ Enjoy simple disassembly with one hand and advanced features designed for convenience and peace of mind. You can feel confident that the Victory's long-lasting batteries will keep you going. Anything is possible with a Victory[®] Scooter!

All models feature...

- Feather-touch disassembly
- Standard front basket
- Delta tiller with port to conveniently charge your scooter



FDA Class II Medical Device¹



SCAN THE QR CODE
TO LEARN MORE ABOUT
VICTORY[®] SCOOTERS!



FOLLOW US ON SOCIAL MEDIA

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

K0801 Group 1 Heavy Duty

VICTORY® 10.2¹



FDA
Class II
Medical
Device¹

SPEED²

Up to
5.2 mph

WEIGHT CAPACITY

400 lbs.

COLORS

FEATURES

- Frosted LED headlights and low-mounted LED pathway light brightens your route
- 18"x17" seat with swivel and sliders
- Integrated storage/cup holder in tiller
- Maximum range up to 13 miles with 3-wheel model and up to 13.5 miles with 4-wheel model (40AH batteries)^{3,4}

WARRANTY

Lifetime limited on frame; 2-year limited on electronics
2-year limited on drivetrain

K0801 Group 1 Heavy Duty

VICTORY® 10¹



FDA
Class II
Medical
Device¹

SPEED²

Up to
5.3 mph

WEIGHT CAPACITY

400 lbs.

COLORS

FEATURES

- LED headlight
- 18"x17" black vinyl seat with 8-point swivel and slider
- Ambient voltmeter
- Maximum range up to 16 miles³

WARRANTY

Lifetime limited on frame; 2-year limited on electronics
2-year limited on drivetrain

K0800 Group 1 Standard

VICTORY® 9¹



FDA
Class II
Medical
Device¹

SPEED²

Up to
5.3 mph

WEIGHT CAPACITY

300 lbs.

COLORS

FEATURES

- LED headlight
- 18"x17" black vinyl seat with 8-point swivel and slider
- Ambient voltmeter
- Maximum range up to 13 miles³

WARRANTY

Lifetime limited on frame; 2-year limited on electronics
2-year limited on drivetrain

K0899 - POWER MOBILITY DEVICE,
NOT CODED BY DME PDAC OR DOES
NOT MEET CRITERIA

VICTORY® LX SPORT¹



FDA
Class II
Medical
Device¹

SPEED²

Up to
8 mph

WEIGHT CAPACITY

400 lbs.

COLORS

FEATURES

- CTS Suspension
- Full LED lighting package includes curb light, rear lights and directional signals
- Puddle light under tiller
- 18"x18" high back seat with sliders
- Maximum range up to 19.7 miles with 200 lb. user weight and up to 13 miles with 400 lb. user weight³

WARRANTY

Lifetime limited on frame; 2-year limited on electronics
2-year limited on drivetrain

1) Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2) Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This can be subject to a variance of (+ or -) 10%.
3) Range varies based on user weight, terrain and battery type.
4) Maximum range up to 11 miles with U-1 batteries.

AUTHORIZED PROVIDER



pridemobility.com

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]