

P326A Vision Sport

Technical Specs

Max. speed up to	8 kph / 5 mph
Range up to	28 km / 18 mi.
Turning radius	50 cm / 20"
Ground clearance	5 cm / 2"
Motor	DC24V, 160W (Rating) 520W (Max)
Controller	Dynamic Shark 50A
Battery	12V / U1 x 2 pcs
Charger	4 A off-board
Gradient	6°
Caster	F : 6" / R : 6" solid tire
Drive wheel	10" foam filled tire
Brake	intelligent, regenerative, electromagnetic brakes



* Actual product may not be exactly as pictured

Features

- Compact mid-wheel drive design for tight space maneuverability
- Dual in-line motors for enhanced efficiency, torque, range and performance
- Low center of gravity design delivers a smooth and stable ride
- New PG VR-2 or Dynamic shark fully programmable controller
- Swing back control arm standard
- Height and width adjustable armrest
- Semi reclining backrest
- Height adjustable seat
- Colors available Red and Blue
- Full suspension



Patented obstacle climb suspension on front castor



No tipping when standing on footplate



Dimensions

	Seat		Backrest	Overall		
	Width	Depth	Height	Length	Height	Width
cm	51	46	55-60	86	123	61
in	20	18	21.5-23.5	34	48.5	24

Weight

	kg	lbs
Base	41	90
Battery/ea	12	26
Seat	19	42
Total w/batt.	84	184
Capacity	136	300