

# Vision Super

## Mid-Wheel Drive -Bariatric

The Vision Super mid-wheel drive delivers outstanding maneuverability in tight compact areas and a smooth and stable ride both indoors and outdoors. Its high performance suspension, together with the advantages of in-line motor technology for enhanced efficiency, torque and range make the Vision Super best in class. Now Available with optional 10" seat elevator.



# Vision Super

## Features

- Mid-wheel drive design delivers outstanding maneuverability in tight, compact spaces and a smooth stable ride both indoors and outdoors even over rough terrain
- Height and width adjustable armrest, semi reclining backrest with height adjustable seat
- Footplate is rated to hold 250lbs and rear casters will not tip
- Dual in-line motors for enhanced efficiency, torque, range and performance
- Lap Belt (Standard with chair)
- Full Suspension

## Specification

Weight Capacity	450 lbs, 400 lbs seat lift
Maximum Speed	5 mph
Per-Range Charge	20 miles
Batteries	22NF (50AH)
Ground Clearance	2.2"
Turning Radius	21"
Overall Length	35"
Overall Width	27"
Drive Wheels	12" flat-free
Casters	6" front, 6" rear
Suspension	Yes full
Seat Width	22", 24"
Seat Depth	22"
Back Height	22"
Warranty	18 months: Motor & Electronics 5 Year/Limited Lifetime Frame
Controller	LinX

## Options

- Solid Seat Pan
- Cup Holder
- Weather Cover
- Rear basket
- Oxygen Tank holder
- Can/Crutch Holder
- Swing-Away Elevating Legrests
- 10" Seat Elevator



## Models

Vision Super	P327
Vision Super with Seat Lift	P3274

