

Dimensions and Specifications



TRAVEL BUGGY 
See *Your* World Again

DASH Ultra-Lite Foldable Power Chair
2nd Generation

Performance Specifications

Driving Range	9.5 mi / 15 km	Turning Radius	24"
Maximum Speed	4.5 mph / 7 kph	Stopping Distance	5 ft (at full speed)
Drive System	Two Motor Drive System	Obstacle Clearance	2"
Motor Type	170 Watt Brushless	Braking System	Electromagnetic
Incline Capability	9.0° / 15.9%	Max Weight Capacity	265 lbs

Dimensional Specifications

Unfolded Dimensions	38" L x 24" W x 34" H	Folded Dimensions	30" L x 24" W x 13.5" H
Chair Gross Weight	43 lbs	Front Wheel Diameter	8"
Chair Net Weight	39 lbs (no battery)	Rear Wheel Diameter	12"
Seating Width	18.5"	Floor to Seat Height	20"
Seat Cushion Size	17" W x 16" D x 2" H	Back Support Height	15"

Other Specifications




Battery Type	Lithium-ion (24V 10Ah)	Charging Time	4 to 6 hours
Batteries Provided	One (1)	Charger Input	100-240V, 50/60Hz ^(international)
Battery Watt Hour	240 Wh (Airline Approved)	Charger Output	29.4V 2.0A
Battery Charge Cycles	~1,000	Lap Belt	Provided
Tires	Non-Puncture	Horn	Built-In



**Live Freely with
Travel Buggy
Electric Wheelchairs**

Convenient to use. Quick to fold. Easy to transport.

TRAVEL BUGGY 
See *Your* World Again

 1-855-56-BUGGY (28449)
 us.travelbuggy.com
 contact@travelbuggy.com



So Light Anyone Can Lift It

Introducing DASH Ultra Lite by Travel Buggy. Weighing only 39 lbs, DASH is our lightest power chair ever. DASH can easily be lifted into the back of any vehicle without the use of ramps or specialty lifts, and is even lighter than most manual wheelchairs! But don't let its lightweight fool you, DASH can still hold up to 265 lbs and offers powerful motors to tackle inclines and moderate outdoor terrains.

Easy to Fold and Transport

Folds quickly like a lawn chair to the size of a suitcase for convenient transport in any vehicle or for storage.



1-Step Fold Sequence



Net weight
39 lbs



Fits any vehicle

Will Travel!

Travel Buggy electric wheelchairs do not require any additional travel batteries or approvals like other "Travel Chairs". Our LiPo batteries are rated for safe airline travel by TCCA (CA), FAA (US) and various other aviation authorities.

AIR CANADA

WESTJET

Southwest

American Airlines

HAWAIIAN AIRLINES

Alaska

AIRFRANCE

BRITISH AIRWAYS

DELTA

UNITED

Air transat

sunwing

jetBlue

easyJet

KLM

RYANAIR



Intelligently Designed

Our Travel Buggy electric wheelchairs are intuitive to drive and convenient to manage.



Joystick can be adjusted in length and right or left position for user comfort



Flip up arm rests for a barrier-free side entry




When powered off, freewheel (push) mode can be utilized and the wheels turn freely

Just Load and Go

DASH is perfect for outings to the park, coffee shop or meet-ups with friends and family. DASH is the ideal power chair for a day out and about.



Use Indoors & Outdoors



Top Speed



Long Driving Range



Bold Colors to Match Your Personality

Our new DASH Ultra Lite power chair offers 3 vibrant colors so your new wheels stand out from the pack. Pick a color as unique as you and make DASH ULTRA-LITE your own.



The information contained herein is correct as at the time of printing, we reserve the right to alter specifications without prior notice