

	EW-66
Warranty	3 Year Limited Warranty
Colors	Red
Weight Capacity	600 lbs
Fully Assembled	Yes
3 or 4 Wheel	4 Wheels
FWD or RWD	Rear Wheel Drive
Max Passengers	2
Airline Approved	No
Seat Dimensions (LxWxH)	(Front) 12" x 16" x 14"
Type of Seat	(Front) Stationary Seat with Armrests
Type of Battery	Deep Cycle Sealed Lead Acid
Overall length	74"
Overall width	26"
Overall height	43.5"
Armrests	Yes
Scooter Weight	261 lbs
Shipping Method	LTL Freight
Top speed (Up To)	15 mph
Distance (Up To)	35 miles
Light(s)	Front headlight with high/low beam, rear running/brake lights, turn signals
Turning Radius	76"
Braking Type	Hand Brake with Electric Brake Assist
Left Hand Control Option	Yes
Horn	Yes
Incline capability	Up to 12 Degrees
Ground Clearance	4.5"
Suspension	Front fork shocks, rear shock absorbers
Type of Tires	Tubeless
Battery Charger	48/60 volt offboard
Speed Control	Low-Med-High
Basket	Rear
Forward/Reverse	Yes
Throttle Type	Twist
Alarm	Yes
Front Tire Size	16" x 3"
Rear Tire Size	16" x 3"
Cell Phone Holder	No
USB Port	No
Voltage	60 volts
Motor Wattage	700 watt
Motor Type	Brushless Motor with Differential
Storage	Basket & Lockable Under Seat
Rear View Mirror	Yes
Extras	Key Fob, Parking Brake