

Jazzy®

1450

The Jazzy® 1450 offers heavy-duty construction and performance components. This bariatric power chair features a 600 lb. weight capacity, front-wheel drive for superb indoor maneuverability in the tightest of spaces and large front drive wheels for extra absorption over obstacles. The Jazzy 1450 combines superb climbing capabilities with front anti-tips for transitioning over grade changes and handling various terrains easily.



ULTIMATE TRACTION DRIVE TIRES



**OVERSIZED REAR CASTERS
OPTIMIZE PERFORMANCE
AND COMFORT**



**EXTENDED SEAT SIZES
UP TO 32" x 28"**

**AMERICA'S FAVORITE
POWER CHAIR**

Learn more online.



gojazzy.com

Specifications

Weight Capacity	600 lbs.
Maximum Speed	Up to 4.5 mph
Ground Clearance	1.5"
Turning Radius	26.5"
Overall Length	38.25"
Overall Width	29.25"
Drive Wheels	14" solid
Front Wheels	4" solid anti-tips
Rear Wheels	9" solid casters
Suspension	Swing Arm
Drivetrain	Two-motor, front-wheel drive
Braking System	Intelligent braking (electronic, regenerative disc brakes)
Standard Electronics	90 amp, PG VR2
Battery Charger	Off-board, 8 amp
Per-Charge Range	Up to 15.7 miles
Battery Requirements	(2) 12 volt, Group 24
Battery Weight	52 lbs. each
Base Weight	195 lbs.
Standard Seat Weight	52 lbs. (high-back)
Standard Seat	24" x 24"
Max. Seat Dimensions	32" x 28"
Warranty (K0826/K0827)	Lifetime limited on frame; 13-month limited on batteries; 13-month limited on electronics; 13-month limited on drive motors
Warranty (K0839)	Lifetime limited on frame; 2-year limited on electronics; 18-month limited on drive motors



Candy Apple Red



Viper Blue

Colors

*Additional charges may apply.

Features and Benefits

- Depth adjustable, fold flat, limited recline high back seat
- Flat-free drive tires
- 60" lap belt
- Large 17" x 13" foot platform
- Dual-layer powder coated frame for increased durability

Seating Options*

- Limited recline high back 20" x 22"
- Pinchless hinge, high back, available sizes: 20" x 24" - 24" x 24"
- Heavy duty contoured captains seat (includes lumbar cutout) available sizes: 26" x 20" or 32" x 28"
- Limited recline, high back with solid seat pan available sizes: 20" x 20" or 20" x 22"
- Pinchless hinge, high back with solid seat pan available sizes: 20" W x 24" or 24" x 22"
- Heavy duty, high back with solid seat pan (includes lumbar cutout) available sizes: 26" x 20" or 32" x 28"

Options

- 17" x 16" foot platform
- Heavy duty drop-in leg rests
- Heavy duty elevating leg rests
- Swing-away joystick

Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger

Models

Jazzy 1450-2VHD-SS	Code: K0826
Jazzy 1450-2VHD-C	Code: K0827
Jazzy 1450-2SPVHD-SS	Code: K0839



(US) 800-800-8586 • (Canada) 888-570-1113 • pridemobility.com •



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.