

## Prime Care P503

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- 450 lbs. weight capacity, Resident Weight 425 lbs (Figure A)
- Best in class height range 7.8" to 30" (Figure B)
- Integrated 76" and 80" length extension
- Eight (8) Function Handheld Pendant Control includes Auto Contour (Figure C)
- 3 DC Motors keeps bed level with true vertical lift
- Four sealed bearings for smooth/quiet lifting
- Eight 3" Heavy Duty casters roll at any height, 4 locking, 2 caster guide locks (Figure D)
- Fixed wall bumper (Figure E)
- Options:
- Assist Rails, Assist Bars
- Bed ends, Six styles and twelve wood grain laminate colors

### PRODUCTS

### Description

P503	Prime Care Model P503
------	-----------------------

### SPECIFICATIONS

### Description

Bed Height Range	7.8" - 30"
Bed Length Range	76" or 80"
Bed Width	36"
Warranty	10 Years on Frame, 15 Years on Welds, 2 Years on Electrical
Weight	210 lbs
Weight Capacity	450 lbs
Bed Height Min w casters	7.8"
Bed Height Max w casters	30"

### ACCESSORIES

### Description

PLTCAB	LTC Assist Bar, 1PR, PR
PLTCAR	LTC Assist Rail, 1PR, PR
35BWBE-B	35" BRNWOOD BEDEND, OKOAK, EA, EA
35BWBE-A	35" BRNWOOD BEDEND, CHERRY, EA, EA
35MCBE-B	35" MILL CRK BEDEND, OKOAK, EA, EA
35MCBE-A	35" MILL CRK BEDEND, CHERRY, EA, EA
42BWBE-B	42" BRNWOOD BEDEND, OKOAK, EA, EA
42BWBE-A	42" BRNWOOD BEDEND, CHERRY, EA, EA
42MCBE-A	42" MILL CRK BEDEND, CHERRY, EA, EA
42MCBE-FM	42" MILL CRK BED END, FIG MAHOG, EA, EA

**ACCESSORIES****Description**

42MCBE-B	42" MILL CRK BEDEND,OKOAK,EA, EA
35MCBE-FM	35" MILL CRK BED END,FIG MAHOG,EA, EA
35BWBE-FM	35" BRNWOOD BED END,FIG MAHOG,EA, EA
42BWBE-FM	42" BRNWOOD BED END,FIG MAHOG,EA, EA