

Prime Care Bed Model P703



HCPCS:



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E

- 550 lbs. weight capacity, Resident Weight 500 lbs
- 7" to 30" height range
- Optional integrated width adjusts to 39" - 42" (Figure A)
- Integrated 76", 80" and 84" sleep surface length (Figure B)
- Two way Red/Green color coded floor lock system (Figure C)
- Eight (8) Function Handheld Pendant Control
- Auto Contour
- 4 DC Motors Synchronized Hi/Lo control system keeps bed level
- Four sealed bearings secured in "G" track
- IV pole holders
- Extended Heel Lift
- Eight 3" Heavy Duty casters roll at any height
- Fixed adjustable wall bumper (Figure D)
- Options:
 - > Staff Control with Trend/Reverse Trend and Cardiac Chair (Figure E)
 - > 88" length Extension (Figure F)
 - > Deluxe Assist Rails, Assist Bars
 - > Bed ends, two styles and three wood grain laminate colors
 - > Battery Back-up, Bed Transporter, Trapeze - Battery Back-up (P703-BBK), Bed Transporter (P903TRP)

PRODUCTS	Description
P703	Prime Care Model P703, 1 ea

SPECIFICATIONS	Description
Bed Height Range	7" - 30"
Bed Length Range	76", 80" or 84"
Bed Width	36"
Unit Dimensions	38.5" (w) x 85.25" (L) x 8" (H)
Warranty	Lifetime on Welds; 15 Year on Frame; 5 years on Electrical
Weight	205 lbs.
Weight Capacity	Accommodate up to 550 lbs. Maximum Resident Weight 500 lbs.
Bed Height Min w casters	7"
Bed Height Max w casters	30"

ACCESSORIES	Description
PLTCDR	LTC Deluxe Rail, 1PR, PR
PLTCAR	LTC Assist Rail, 1PR, PR

ACCESSORIES**Description**

PLTCAB	LTC Assist Bar,1PR, PR
P703LE	P703 Length Extension 88",1ea, EA
P703-BBK	P703 Battery Backup,1ea, EA
P703SC04	Staff Control Panel f/P703 ,1ea, EA
P703WE	P703 Width Extension 39-42",1pr, PR
35BWBE-B-SC04	35" BRNWOOD BEW/SC04,OKOAK,EA, EA

PARTS**Description**

SP02-TH21060003	HandPendant,8-Button,703/1750/15904C, EA
SP02-PPB1750-3402E	Power Cord, 503/1750/15904C/200, 1/ea, EA
SP01-PCB900-24BS	BrakeStop,900/601/901/902/903/15902C/3, EA
SP01-UF16108061E	Hd/FtBrdBolt,P703/503,1ea, PKG
SP01-PCB300-01-57E	FtEdRchtAssy,900/01/02/03/TL500/703,1/PR, PR
SP01-TL500-29E	NonLckCastr,900/901/2/3/601/301/500/P200, EA