



NEW & IMPROVED CRADLE

13244
with Sling
(Sold Separately)

CAN FIT UNDER A LOW BED!

BARIATRIC
600 LB
WEIGHT CAPACITY

REFER TO PAGE 182-190 FOR LIFT ACCESSORIES AND SLINGS.

Bariatric Battery Powered Lift with Six Point Cradle

13244 1/cs

FEATURES:

- 600 lb. weight capacity.
- To charge, simply plug in unit. Battery does not need to be removed to be charged. Lift will not operate while plugged in.
- Jumbo actuator pump with emergency button can switch to manual mode lowers patient down safely.
- Base legs adjust easily and lock securely into open position with padded shifter handle.
- **Ideal for homecare use.**
- Extra wide coated cradle with 360° rotation with 6 point hook up.
- High lift range allows transfer to high surfaces.
- Dual, 3" heavy duty front casters.
- Warranty: Lifetime on frame, 1 year on components and electrical parts.

NEW & IMPROVED CRADLE

BARIATRIC
600 LB
WEIGHT CAPACITY



CAN FIT UNDER A LOW BED!

13245
with Sling
(Sold separately)

Bariatric Battery Powered Patient Lift with Rechargeable, Removable Battery

13245 1/cs

FEATURES:

- Same features as 13244
- Includes removable, rechargeable battery, (item #13258) external charger (item #13255) and wall mounted bracket (item #13260).
- **Ideal for institutional use.**

To avoid lift downtime - order an additional re-chargeable battery, item # 13258. One for the lift and one for the charging station.

REPLACEMENT PARTS:



External Charger with Power Cord

13255 1/ea



Wall Mounted Bracket

13260 1/ea



Rechargeable Battery

13258 1/ea

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*

SPECIFICATIONS

Item #	Boom Maximum Height	Boom Minimum Height	Base Width (open)	Base Width (closed)	Base Height	Base Clearance	Casters	Product Weight	Carton Shipping Weight	Weight Limit
13244	74.5"	32"	45.75"	25.75"	4.25"	2"	3" Front/4" Rear	127 lbs.	139 lbs.	600 lbs.
13245	74.5"	32"	45.75"	25.75"	4.25"	2"	3" Front/4" Rear	127 lbs.	139 lbs.	600 lbs.