

Invacare® IVC™ Tracer® SX5 Wheelchair

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

IVC™ TRACER® SX5 WHEELCHAIRS

Quantity	Model Number	Description	List Price
<input type="checkbox"/>	TRSXS56FBP	16" x 16" Frame with Flip-Back Fixed Height Desk Length Arm (9153653563)	\$1,035.20 K0003
<input type="checkbox"/>	TRSXS58FBP	18" x 16" Frame with Flip-Back Fixed Height Desk Length Arm (9153653564)	\$1,035.20 K0003
<input type="checkbox"/>	TRSXS58FBFP	18" x 16" Frame with Flip-Back Fixed Height Full Length Arm (9153653565)	\$1,035.20 K0003
<input type="checkbox"/>	TRSXS50FBP	20" x 16" Frame with Flip-Back Fixed Height Desk Length Arm (9153653560)	\$1,035.20 K0006 E2201
<input type="checkbox"/>	TRSXS50FBFP	20" x 16" Frame with Flip-Back Fixed Height Full Length Arm (9153653561)	\$1,035.20 K0006 E2201
<input type="checkbox"/>	TRSXS52FBP	22" x 16" Frame with Flip-Back Fixed Height Desk Length Arm (9153653562)	\$1,035.20 K0006

STOCK FRONT RIGGING OPTIONS

Quantity	Model Number	Description	List Price
<input type="checkbox"/>	T93HAP	Hemi Footrests with Aluminum Footplates and Heel Loops (9153653624)	\$99.13 E0951
<input type="checkbox"/>	T93HCP	Hemi Footrests with Composite Footplates and Heel Loops (9153653625)	\$99.13 E0951
<input type="checkbox"/>	T94HAP	Hemi Elevating Legrests with Aluminum Footplates and Padded Calf Pads (9153653627)	\$191.86 K0195 E0990
<input type="checkbox"/>	T94HCP	Hemi Elevating Legrests with Composite Footplates and Padded Calf Pads (9153653628)	\$191.86 K0195 E0990
<input type="checkbox"/>	AHL4A_19911	Hemi Smartleg® Articulating Legrests with Aluminum Footplates and Padded Calf Pads (9153619911)	\$341.09 K0053
<input type="checkbox"/>	T93HEP	Hemi Footrests without Heel Loops, Composite Footplate (9153653626)	\$69.28
<input type="checkbox"/>	T94HEP	Hemi Elevating Legrests with Non-padded Calf Pads and Composite Footplate (9153653629)	\$117.25

All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare® IVC™ Tracer® SX5 Wheelchair

FRAME AND SEAT SELECTIONS

CHAIR TYPE

- TRSX5 Tracer SX5 Builder Wheelchair \$1,055.50 K0003

FRAME TYPE

The below selections are based on 24" wheels and 8" casters, see seat-to-floor matrix.

- ADULT Adult 19 1/2" No Charge
 HEMI Hemi 17 1/2" No Charge

FRAME WIDTH AND DEPTH

- WD46 14" Wide x 16" Deep No Charge A990
 WD66 16" Wide x 16" Deep No Charge
 WD86 18" Wide x 16" Deep No Charge
 WD06 20" Wide x 16" Deep No Charge E2201
 20" and 22" widths have 300 lbs. weight limit requiring 1228 casters (upcharge necessary, see "Casters" section), and 24" rear wheels (no upcharge).
 WD26 22" Wide x 16" Deep \$152.08 K0006 E2201
 20" and 22" widths have 300 lbs. weight limit requiring 1228 casters (upcharge necessary, see "Casters" section), and 24" rear wheels (no upcharge).

ARM AND BACK OPTIONS

ARM TYPE

The below arm types are flip-back, non-removable.

- 18FB Full Length Fixed Height Flip-Back 10" No Charge A990
 19FB Full Length Adj. Height Flip-Back 10"-14" \$139.23 E0973
 28FB Desk Length Fixed Height Flip-Back 10" No Charge A990
 29FB Desk Length Adj. Height Flip-Back 10"-14" \$139.23 E0973

BACK TYPE

- BH16 Back Height 16" Standard
 BH18 Back Height 18" No Charge

COLOR OPTIONS

UPHOLSTERY COLORS

- U240 Black Nylon Upholstery Standard
 U550 Midnight Blue Vinyl Upholstery No Charge
 U67 Black Vinyl No Charge
 U86 Ruby Red Vinyl No Charge

POSITIONING STRAP OPTIONS

SEAT POSITIONING STRAPS

- 1515 Push Button Style Positioning Strap \$53.55 E0978 A990
 1321BK Hook and Loop Seat Positioning Strap \$42.84 E0978

CASTERS

CASTER OPTIONS

- AT902 6" Composite Urethane 3 Spoke No Charge
 AT903 8" x 1 1/4" Composite Urethane No Charge
 1255 8" x 1" Solid Rubber Dark Grey Standard
 1228 8" x 1 3/4" Semi Pneumatic Casters Lt Grey \$59.98 A990

WHEELS

WHEEL OPTIONS

- U1763BK 20" Composite Urethane Tires \$121.02 A990
 U1764BK 22" Composite Urethane Tires \$121.02 A990
 2244BK 24" Composite Pneumatic Tires \$148.51 A990
 4244BK 24" Comp Pneumatic w/FF Inserts \$238.68 E2213
 U2222C 24" Composite Urethane Corded Tires Standard
 U2244BK 24" Composite Urethane Treaded Tires \$128.52 A990

WHEEL OPTIONS (ONE ARM DRIVE)

The below selections are not available with the amputee attachment.

- 9456L One Arm Drive with U2222 LH \$642.84 E0958
 9456R One Arm Drive with U2222 RH \$642.84 E0958

HANDRIM OPTIONS

- COM Composite Handrims Standard
 Available with 24" rear wheels only.
 AHDRM Aluminum Handrims No Charge
 Not available with U2222C rear wheels.
 CHDRM Chrome Plated Handrims No Charge

WHEEL LOCKS

- 9690B Pull-To-Lock Wheel Locks (Std on Hemi) No Charge
 9691B Push-To-Lock Wheel Locks (Std on Adult) No Charge

ACCESSORIES

BACK ACCESSORIES

- STC Stencil Back Upholstery \$31.06 A9270
 1160 Carrying Pocket \$57.83 A9270

SEAT ACCESSORIES

- L8651 Omit Seat Upholstery Less \$10.61

MISCELLANEOUS ACCESSORIES

- 1360 Adjustable Rear Anti-Tipper For 8" Casters \$87.82 E0971
 1360A Adj Rear Anti-Tipper For 6" Casters \$87.82 E0971
 9357L Amputee Attachment \$100.78 E0959
 1329BK Toggle Wheel Lock Extension (8") \$41.77 E0961
 1335 Non-Folding Device \$42.85 K0108
 Not available on 20" and 22" widths.
 Not available with One Arm Drive.
 1514 Overhead Anti-Theft Device \$49.27
 Not available on 14" and 22" widths.
 1497 Oxygen Holder Package Assembly \$161.72 E2208
 Not available with One Arm Drive.
 1496 IV/Oxygen Holder Package Assembly \$169.22 K0105 E2208
 Not available with One Arm Drive.
 EXT2 Seat Extension Kit 2" No Charge

All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare® IVC™ Tracer® SX5 Wheelchair

FRONT RIGGING OPTIONS

FRONT RIGGINGS

- T93HE Hemi Footrests without heel loop, Composite Ftpl\$69.28 A9900
Not available with 14" wide.
- T93HA Hemi Footrests Aluminum Footplates\$99.13
Not available with 14" wide.
- T93HC Hemi Footrests \$99.13
Not available with 14" wide.
- T94HE Hemi Elevating Legrest Composite Ftpl/Blk Plastic Calf Pad.....\$117.25 A9900
Not available with 14" wide.
- T94HA Hemi Elev. Lgrt Aluminum Ftpl\$191.86 K0195
Not available with 14" wide.
All calf pads are black.
- T94HC Hemi Elevating Legrests\$191.86 K0195
Not available with 14" wide.
All calf pads are black.
- AHL4A Hemi Smartleg Articulating Lgrt Alum Ftpl\$341.09 K0053
Not available with 14" wide.
All calf pads are black.
- AHL4C Hemi Smartleg Art. Lgrt Comp Ftpl\$319.77 K0053
Not available with 14" wide.
All calf pads are black.
- AHS4A Hemi Smartleg Art. Lgrt-14" W Chairs\$341.09 K0053
All calf pads are black.
- T93HS Hemi Footrests Aluminum Footplates 14" W.....\$95.93 A9900
- T94HS Hemi Legrests Aluminum Footplates 14" W.....\$186.53 A9900
All calf pads are black.

FOOTPLATE TYPE

- 1350 Extra Large Footplates\$84.21 K0041
Not available with 14" wide.
- AT5543 Adjustable Angle Flip-Up Footplates\$71.42

FOOTREST/LEGREST ACCESSORIES

- 1820 Non-Removable Device\$12.79 A9900
- 1346 Heel Loops\$41.57 E0951
Already included on swingaway footrests T93HA, T93HC.


Specifications	
Seat Widths	14", 16", 18", 20" or 22"
Overall Width	Seat Width plus 8 1/2"
Seat Depth	16"
Overall Depth	3 1/8" (18 x 16 without Front Riggings)
Back Heights	16" or 18"
Seat-to-Floor	17 1/2" or 19 1/2"
Weight Limit	250 lbs. (300 lbs. on 20" or 22" widths)
Chair Weight	34 3/4" (18 x 16 without Front Riggings)
Shipping Weight	44 1/4 lbs.

20-22 Wide (300 lb. weight limitation) Wheelchairs only: 8 x 1 3/4 inch semi-pneumatic front casters MUST be installed. DO NOT replace with any other size/style front caster. Otherwise, injury or damage may occur.


6-INCH FRONT CASTERS			
SEAT-TO-FLOOR HEIGHT (IN INCHES)	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
15 1/2	TOP	20 inch	TOP
16 1/2	MIDDLE	22 inch	TOP
17 1/2	BOTTOM	24 inch	TOP
17 1/2	BOTTOM	20 inch	BOTTOM

8-INCH FRONT CASTERS			
SEAT-TO-FLOOR HEIGHT (IN INCHES)	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
17 1/2	TOP	20 inch	BOTTOM
17 1/2	TOP	24 inch	TOP
18 1/2	MIDDLE	22 inch	BOTTOM
19 1/2	BOTTOM	24 inch	BOTTOM

Rear Wheel Mounting Positions - Adult/Hemi Frame Style



6 or 8-Inch Fork Mounting Positions



All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare Corporation

www.invacare.com

USA

One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

©2020 Invacare Corporation
All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.



Yes, you can.®

Revised January 1, 2020