

## Bariatric Sentra Extra-Heavy-Duty Wheelchair

HCPCS: K0007,K0007+K0195 OR E0990



- Weight Capacity: 500 lbs
- Carbon steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Reinforced steel gussets at all weight bearing points provide additional strength
- Durable, heavy-gauge reinforced naugahyde upholstery
- Comes standard with carry pocket on backrest
- Reinforced side frame, caster journals and dual cross bars
- Dual-position front fork
- Dual-position front frame
- Composite Mag-style wheels with chrome hand rim are lightweight and maintenance free
- Hemi-adjustable caster is optional (Part #STDS1074)
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Padded armrests provide added patient comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Comes with push-to-lock wheel locks

### PRODUCTS

### Description

STD20DDA-ELR	Bariatric Sentra Extra-Heavy-Duty Wheelchair, CS
STD20DDA-SF	Bariatric Sentra Extra-Heavy-Duty Wheelchair, CS
STD22DDA-ELR	Bariatric Sentra Extra-Heavy-Duty Wheelchair, CS
STD22DDA-SF	Bariatric Sentra Extra-Heavy-Duty Wheelchair, CS
STD24DDA-ELR	Bariatric Sentra Extra-Heavy-Duty Wheelchair, CS
STD24DDA-SF	Bariatric Sentra Extra-Heavy-Duty Wheelchair, CS
STDS004010	Bariatric Sentra Extra-Heavy-Duty Wheelchair, EA
STDS1016-8	Bariatric Sentra Extra-Heavy-Duty Wheelchair, BAG
STDS1070	Bariatric Sentra Extra-Heavy-Duty Wheelchair, EA
STDS1072-8	Bariatric Sentra Extra-Heavy-Duty Wheelchair, EA
STDS1S2650	Bariatric Sentra Extra-Heavy-Duty Wheelchair, EA
STDS1S2652	Bariatric Sentra Extra-Heavy-Duty Wheelchair, EA
STDS1S2654	Bariatric Sentra Extra-Heavy-Duty Wheelchair, EA
STDS4019L	Bariatric Sentra Extra-Heavy-Duty Wheelchair, EA
STDS4019R	Bariatric Sentra Extra-Heavy-Duty Wheelchair, EA
STDS4020	Bariatric Sentra Extra-Heavy-Duty Wheelchair, EA
STDS8M4603	Bariatric Sentra Extra-Heavy-Duty Wheelchair, EA
STDSADDASPL	Bariatric Sentra Extra-Heavy-Duty Wheelchair, EA

### SPECIFICATIONS

### Description

Back Height	18"
-------------	-----

SPECIFICATIONS	Description
Carton Shipping Weight	60 lbs.   62 lbs.   72 lbs.
Height	37"
Length with Riggings	42"
Product Weight with 20 Seat	54 lbs.
Product Weight with 22 Seat	56 lbs.
Product Weight with 24 Seat	66.25 lbs.
Seat Depth	18"
Seat To Floor Height	17.5" - 19.5"
Seat Width	20"   22"   24"
Warranty	Limited Lifetime
Weight Capacity	500 lbs
Width Closed	13"
Width Open	29"   31"   32.5"
Arm Rest Type (s)	Detachable Desk Arm   Detachable Full Arm   Detachable and Adjustable Height Desk Arm   Detachable and Adjustable Height Full Arm

ACCESSORIES	Description
WASL	Limb Support,Left,1/cs, CS
WASR	Limb Support, Right,1/cs, CS
STDS831	Heel Loops, 1pr/bx, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS803	Oxygn Tank Holdr for W/C,1/cs, CS
STDS821	Overhead AntiTheft w/IV Hks,1/bx, BX
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDELR-AL	ArticulatingElev.Legrest, Std,1pr/cs, CS
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS801	Brake Handle Ext.,8", 1/ea, EA
STDS860	Forwrđ Stablzrs/Univ.Fit , 1pr/bx, BX
STDS2420	Seat Rail Ext.Kit for Sentra,24x20,1/bx, BX
STDS802	Anti-Tipprs,w/Wheels,Adj,1pr/bx, BX

PARTS	Description
STDS1051	Dust Cap, fork stem, 1/ea, EA
STDS871	Hand Grip,Sentra,1/ea, EA
STDS1026	Front Bumper for Ft.Rst/LegRest,1/ea, EA
STDS1002-1	Axle for 8" caster, 1/ea, EA
STDS1060	Wheelchair Fork Bearing,1/ea, EA
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, EA
STDS8M4601	ForkStemBearing,SSP2,Cruiser,1/ea, EA
STDS1070	Seat Rail Guide, Removable Arm, 1/ea, EA
STDS4020	Tipping Sleeve 1/ea, EA
STDS002001	Arm Pad Screw,1/ea, EA
STDS005200	SlideTube,Sentra,1/ea, EA
STDS4011	BrakeTip, 1/ea, EA
STDS2329BLK	BackUphol black,SentraHD 20inx 18in,1/bx, BX
STDS2331BLK	Back Uphol 22"x18",Black,1/bx, BX
STDS2333BLK	Back Uphol 24"x18", Black,1/bx, BX
STDS1052	Rear AxleBolt,SentraHD,1/ea, EA
STDS2A5222	Crossframe,22"SentraHD/1ea, EA

**PARTS****Description**

STDS4010R	Brake Asmb.Right Push to Lock,1/ea, EA
STDS1074	6"Hemi Caster,SentraHD,1ea, EA
STDS1000HD	Wheel,24"Hvy.Dty.Hand Rl.Sentra,1/ea, EA