



Invacare® IVC™ 9000 XDT Wheelchair

The Invacare IVC 9000 XDT wheelchair is designed for individuals who demand an extra wide, heavy-duty or tall chair with seat-to-floor heights ranging from hemi through tall. The high-performance 9000 XDT provides a 350-lb. weight capacity and components interchangeable with the other 9000 series wheelchairs.

FEATURES

- Durable, low-maintenance, triple chrome-plated, carbon steel frame is long-lasting
- Tall frames are designed to accommodate taller individuals and provide a more comfortable ride
- Multiple frame color variations for customizing
- Triple-axle positions allow wide variation of seat-to-floor heights
- Urethane rear tires mounted on “no flex” wheels offer superior performance
- Urethane casters for better ride and lighter weight
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Meets RESNA standards for manual wheelchairs in product performance

Specifications

Dimensions	
Seat sizes	
Widths	16", 17", 18", 19", 20", 22" tall to adult to hemi frame
Depths	16", 17", 18", 19", 20"
Seat-to-floor heights	17.5", 19.5" 21.5" with different wheel and caster packages
Back heights - adjustable	15", 16", 17", 18", 19"
Front riggings	
Hemi swingaway footrest	17"-21"
Hemi swingaway padded elevating legrest	17.5"-21"
Super-hemi swingaway footrest	11"-17"
Super-hemi swingaway padded elevating legrest	11"-17"
Articulating elevating legrest	17"-25"
(Footplates available in aluminum or composite)	
Arm styles	
Full-length fixed height, padded, flip back, removable	
Desk-length fixed height, padded, flip back, removable	
Full-length adjustable height, padded, flip back, removable	
Desk-length adjustable height, padded, flip back, removable	
(All arm styles available conventional or space saver)	
Weight capacity	350 lb.
Product weight (9XDT/WD08/28)	37 lb.
Warranty	lifetime on frame and crossbraces