

	EW-72
Warranty	3 Year Limited Warranty
Colors	Red, Blue, Black, Yellow, Orange
Weight Capacity	500 lbs
Fully Assembled	Yes
3 or 4 Wheel	4 Wheels
FWD or RWD	Rear Wheel Drive
Max Passengers	1
Airline Approved	No
Seat Dimensions (LxWxH)	18.5" x 18.5" x 30"
Type of Seat	Adjustable Captains seat with head rest
Type of Battery	Deep Cycle Sealed Lead Acid
Overall length	56"
Overall width	27"
Overall height	52"
Armrests	Yes
Scooter Weight	275 lbs
Shipping Method	LTL Freight
Top speed (Up To)	18 mph
Distance (Up To)	43 miles
Light(s)	Front headlight with high/low beam, rear running/brake lights, turn signals
Turning Radius	92"
Braking Type	Hand Brake with Electric Brake Assist
Left Hand Control Option	Yes
Horn	Yes
Incline capability	Up to 12 Degrees
Ground Clearance	6"
Suspension	Front leaf springs, rear shock absorbers
Type of Tires	Tubeless
Battery Charger	48 volt offboard
Speed Control	Variable Speed Control Knob
Basket	Rear
Forward/Reverse	Yes
Throttle Type	Twist
Alarm	Yes
Front Tire Size	14" x 3"
Rear Tire Size	16" x 3"
Cell Phone Holder	Yes
USB Port	Yes
Voltage	48 volts
Motor Wattage	700 watt
Motor Type	Brushless Motor with Differential
Storage	Basket & Lockable Under Seat
Rear View Mirror	Yes
Extras	Key Fob, Parking Brake