



LC-358M shown in  
Cloud 9 Stone fabric

MADE IN  
**USA**  
OF FOREIGN & DOMESTIC COMPONENTS

# HERITAGE COLLECTION

Comfort and style complement one another perfectly in the **Heritage Collection** from Pride® Power Lift Recliners. With an array of plush fabrics to fit any décor and a superior design that provides exceptional strength, Heritage Collection delivers style and durability that will last for years to come. Pride's Heritage Collection offers numerous standard features, superb comfort and quiet and smooth performance, all at a great value.

AMERICA'S #1  
POWER LIFT RECLINER

POWER  
LIFT  
RECLINERS

US: 1-800-800-8586 • CAN: 1-888-570-1113  
pridemobility.com |    

**PRIDE**  
MOBILITY  
LIVE YOUR BEST™

# HERITAGE COLLECTION

AMERICA'S #1  
POWER LIFT RECLINER

POWER LIFT  
RECLINERS

**FIRST UP**  
TECHNOLOGY

Available on the Heritage Collection Series (except 358XL & XXL)



| SPECIFICATIONS      | Petite Wide                            | Small                                  | Medium                             | Large                                | Heavy-Duty                           | Very Heavy-Duty                      |
|---------------------|--|--|------------------------------------|--------------------------------------|--------------------------------------|--------------------------------------|
|                     | LC-358PW <sup>2,5</sup>                | LC-358S <sup>5</sup>                   | LC-358M <sup>2,5</sup>             | LC-358L <sup>2,5</sup>               | LC-358XL <sup>4,6</sup>              | LC-358XXL <sup>4,6,7</sup>           |
|                     | 3-Position Chaise Lounger 5'3" & Below | 3-Position Chaise Lounger 5'3" & Below | 3-Position Chaise Lounger 5'4"- 6' | 3-Position Chaise Lounger 5'8"- 6'2" | 3-Position Chaise Lounger 5'6"- 6'1" | 2-Position Chaise Lounger 5'6"- 6'1" |
| Weight Capacity     | 375 lbs.                               | 375 lbs.                               | 375 lbs.                           | 375 lbs.                             | 500 lbs.                             | 600 lbs.                             |
| Seat-to-Floor       | 19"                                    | 19"                                    | 20"                                | 21"                                  | 20"                                  | 20"                                  |
| Seat Depth          | 19"                                    | 19"                                    | 20"                                | 22"                                  | 22"                                  | 22"                                  |
| Seat Width          | 24"                                    | 18"                                    | 20"                                | 22"                                  | 26"                                  | 30"                                  |
| Top of Back to Seat | 23.5"                                  | 24"                                    | 24.5"                              | 25.5"                                | 29"                                  | 29"                                  |
| Back Style          | Blown Fiber Button Back                | Blown Fiber Button Back                | Blown Fiber Button Back            | Blown Fiber Button Back              | Blown Fiber Button Back              | Blown Fiber Button Back              |
| Fully Padded Chaise | Yes                                    | Yes                                    | Yes                                | Yes                                  | Yes                                  | Yes                                  |
| Heat & Massage      | Optional                               | Optional                               | Optional                           | Optional                             | Optional                             | NA                                   |
| Seating Type        | Total Comfort Seating                  | Total Comfort Seating                  | Total Comfort Seating              | Total Comfort Seating                | Sinusuous Spring                     | Sinusuous Spring                     |
| Distance From Wall  | 18"                                    | 18"                                    | 17"                                | 19"                                  | 19"                                  | 14"                                  |
| Overall Chair Width | 40"                                    | 35"                                    | 35"                                | 37"                                  | 39.5"                                | 43"                                  |
| Chair Weight        | 130 lbs.                               | 118 lbs.                               | 117 lbs.                           | 120 lbs.                             | 159 lbs.                             | 168 lbs.                             |
| Warranty            | 7-Year Prorated*                       | 7-Year Prorated*                       | 7-Year Prorated*                   | 7-Year Prorated*                     | 7-Year Prorated                      | 7-Year Prorated                      |

## HERITAGE FABRIC SELECTION

### CLOUD 9



BLACK CHERRY



PACIFIC



STONE



WALNUT

### LEXIS STA-KLEEN®

sta\*kleen.  
PERFORMANCE FABRIC



BLACK



MUSHROOM



BURGUNDY



CHESTNUT

### ULTRALEATHER™



BUFF



FUDGE



GARNET



PECAN



CHARCOAL

### ARIA<sup>3</sup>

CRYPTON.  
super fabric



ESPRESSO



COOL GREY



RED



LAZULI



SAND

- Batteries not included.
- Standard footrest extension.
- Crypton® Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.
- Removable, fully padded chaise.
- Unique spring design and a one-piece soft foam seat for outstanding relaxation.
- Features Dual Motor Sync System.
- Only available in Stone, Walnut, Lexis Sta-Kleen, Chestnut and Mushroom.
- Standard with footrest extension and lumbar pillow.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.  
Lifetime warranty on steel lift frame.  
The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



Look for the Quick Ship Indicator for the shortest lead times available.

NEW!

## FEATURES AND BENEFITS

- Standard footrest extension on 358PW, M and L for additional leg support while reclined
- Easily adjustable leg levelers provides additional support for uneven surfaces
- Easy-to-use hand control with USB charger port (not available when adding optional heat and massage)
- Standard head and arm covers
- Large dual pockets for extra storage
- Industry-first Total Comfort Seating features a unique spring design and one-piece soft foam seat for outstanding relaxation (except LC-358XL and LC-358XXL)
- Quiet and smooth lift and scissor mechanism system
- Integrated battery backup ensures chair is operational when power goes out
- Energy-efficient power box equipped with self-testing electronics and 9' cord
- Premium furniture grade laminate/hardwood frame for added strength and durability

## OPTIONS\*

- Heat and massage (except XXL)
- Cleaning kit
- Head pillow
- Lumbar pillow
- Wireless remote (not available on LC-358XL and LC-358XXL)
- FirstUp® technology (except LC-358XL and LC-358XXL)
- Non-skid leg levelers

Additional lead time may apply for some options. \*Additional charges apply for all options.

(US) 800-800-8586 | (CAN) 888-570-1113  
pridemobility.com | prideliftchairs.com

**PRIDE**  
MOBILITY  
LIVE YOUR BEST®