

bellavita DIVE

Seat height adjusts from 2.6" to 18.9" to fit almost any tub

Easy-to-use floating remote control

New removable backrest with 50° recline

Installs easily in almost any bathtub — no plumbing or electrical work required


Improved dual column design provides optimal stability for up to 310 lbs.

The lightest and lowest bath lift on the market now helps users dive deeper and get cleaner!

Enjoy relaxing, immersive baths at the touch of a button with the Bellavita Dive, the latest innovation in our beloved line of Bellavita bath lifts. This next generation automated lift lowers users to a minimum height of 2.6 inches, reclines to 50 degrees, and now features a removable backrest for our most immersive bathing experience ever. Its redesigned, larger seat with hygiene cutout makes it easy to get clean with minimal reaching and shifting. An easy-to-use floating remote control keeps users in the driver's seat as they get in, out, and find the perfect height and angle for their bath.

With a new, patented dual column design, the Dive supports up to 310 pounds with unparalleled stability. And with an ergonomic, slimline design and maximum height of 18.9 inches, the Bellavita Dive can be used in almost any bathtub — no plumbing or electrical work required.

Why should they settle for a dip when they can take a dive? The new Bellavita Dive from Drive Medical lets users enjoy the comfort of fully immersive baths safely and independently.



Larger seat with hygiene cutout

ITEM #	PRODUCT DESCRIPTION	UPC	UOM
477400252	Bellavita Dive Reclining Bathlift with classic white cover	822383019604/ea	1/ea, 1/cv

PRODUCT SPECS

Color	White
Product Dimensions	27.6" (W) x 40.2" (L) x 2.6" – 18.9" (H)
Product Weight	10.89 lbs.
Product Weight without Hand Control	21 lbs.
Product Weight with Hand Control	21.8 lbs.
Seat Dimensions	15" (W) x 20.7" (D)
Backrest Dimensions	11.8–12.6" (W) x 23.6" (H)
Backrest recline	50°
Baseplate Dimensions	11.8" (W) x 23.6" (H)
Weight Capacity	310 lbs.
Warranty	5 years on the bath lift, 2 years on the handset

