

Invacare® Top End® Terminator™ Everyday/Titanium Upgrade Ultra Lightweight Active Performance Rigid Wheelchair

Price List and Order Form

Effective January 1, 2006 Revised January 1, 2015

For accuracy of ordering, please fax order to Top End at 727-522-1007 or e-mail to TopEndOrders@invacare.com

For questions call 800-532-8677 or 727-522-8677

To order Literature: 800-828-6282

www.topendwheelchair.com



Provider Name _____ Date _____ Top End Job Number _____
Account # _____ P.O. # _____ Phone # _____ Fax # _____
Address _____ City _____ State _____ Zip Code _____
Provider Primary Contact Name _____ Email _____ Phone No _____
Customer Name _____ Email _____ Phone No _____
Customer Address _____ City _____ State _____ Zip Code _____

Indicate: ☐ Order or ☐ CAD or ☐ Quote Ship order to: ☐ Provider ☐ Customer Order will ship to Provider unless otherwise noted.

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

☐ RACE* 101,102,103 CAD (see definition on page 3) .. 125.00-225.00

PO required _____ E-mail address: _____

FRAME TYPE:

- ☐ TED Top End Terminator Everyday (K0005)..... 2,490.00 ¹
☐ TEDTiU Titanium Upgrade (K0108) must select TED 1000.00 ²
☐ TED61 Heavy Duty Upgrade (K0005) 500.00 ³

FRONT FRAME STYLE: (CHOOSE ONE)

- Tapered Front End Standard
☐ Tapered Seat Width No Charge ⁴
Taper Difference rear to front ☐ -1" ☐ -2" ☐ -3" ☐ -4" Other _____
☐ Non Tapered V Front End No Charge ⁵

FOOTREST OPTIONS FOR TAPERED FRONT END OR TAPERED SEAT WIDTH:

- One Piece Tubular Aluminum Footrest Standard ⁶
☐ TED96 One Piece Tubular Titanium Footrest Upgrade ... 200.00
☐ TED10 Plastic Footrest Cover (tubular footrest only).....75.00
☐ TED89 Carbon Fiber Footrest Cover(tubular footrest only)125.00
☐ Adjustable Angle One Piece Footplate..... No Charge
☐ TED191 Adjustable Angle Footplate w/Risers.....95.00 ⁶
Footrest Height from Floor _____

FOOTREST OPTIONS FOR NON TAPERED V FRONT END

- Adjustable Angle One Piece Footplate Standard ⁷
☐ TED191 Adjustable Angle Footplate w/Risers.....95.00 ⁶
Footrest Height from Floor _____

LOWER REAR FRAME

- Curved Standard
☐ TEDTi5 Square 100.00

FOOTNOTES

- ¹ Terminator Everyday built of 6061T6 Aluminum in 1" tubing only.
² Titanium frame upgrade built of 1" 3-2.5 seamless titanium. Must select TED.
³ 350 lb weight limit for TED / 400 lb weight limit for TEDTiU. Restrictions: Only available with folding back, square rear frame, single camber system (A). Solid tires, swing-away arms, one-sided forks, anodized forks and suspension not available. TED Heavy Duty will be built of 4130 Chrome-moly steel.
⁴ Tapered seat width tapers rear to front up to a maximum difference of 4". The front seat width measures the same as the outside of the footrest or straight so there is no tapered or V front end. Not available with Ergo seat or screw-on seat upholstery.
⁵ Not available w/tapered seat width or tubular footrest.
⁶ One piece tubular footrest is open and does not have a cover. Either a plastic or carbon fiber footrest cover is available. Tubular footrest not available with V style front end. For short leg lengths see knee-to-heel range chart on page 2 and consider option TED191.
⁷ Non-tapered V style front end requires adjustable angle one-piece footplate.

SUSPENSION OPTIONS:

- ☐ TEDTi4 Rear Suspension Option (E1015)..... 350.00 ⁸
Weight of User _____

SEAT WIDTH: ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18"
☐ Other _____ (20" or wider require HD Option)

SEAT DEPTH: ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18"

FRONT SEAT-TO-FLOOR HEIGHT:

☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20" ☐ 20.5" ☐ 21"


REAR SEAT-TO-FLOOR HEIGHT: ⁹

☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5"
☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20" ☐ 20.5" ☐ 21"

Knee-to-Heel Range

Front Seat-to-Floor

17"	12 ¼" – 16 ¼"	20"	15 ¼" – 19 ¼"
18"	13 ¼" – 17 ¼"	21"	16 ¼" – 20 ¼"
19"	14 ¼" – 18 ¼"		



FOOTNOTES

- ⁸ Front suspension forks not available at this time but can be ordered from Frog Legs and installed by provider. Not compatible with 85 or 90 degree front frame angle, seat depths less than 16" or one-sided forks. Maximum weight capacity is 250 lbs. Minimum rear seat height is 16" and for 26" or 25" knobby wheels is 17" with rear suspension option.
⁹ All rear seat heights are without suspension. 14.5" is lowest rear seat height for 24" wheels with Primo, Kenda or solid tires, 15" for 24" with Schwalbe and 1 3/8" tires and 25" wheels except w/knobby, 15.5" for 24" w/knobby tires and 16" for 25" knobby or 26" wheels.

Invacare®Top End® Terminator™ Everyday/ Titanium Upgrade

Ultra Lightweight Active Performance Rigid Wheelchair

ERGO SEAT OPTION


- ☐ TED104 Ergo Seat w/tapered front end 250.00 ¹⁰
☐ TED110 Ergo Seat w/V style front end 350.00 ¹⁰

REAR TO FRONT CASTER DISTANCE:

STD= Seat Depth plus 3" ☐ Other _____

FRONT FRAME ANGLE: ¹¹

- ☐ 70° ☐ 75° ☐ 80° ☐ 85° ☐ 90°



FOOT SPACE WIDTH: standard as per diagram above

Exception: Seat widths of 12", 13", and 14" will have footspace width of 8" ☐ Other _____

BACK HEIGHT: (CHOOSE ONE) ¹²

- ☐ Adjustable Height Back without Push Handles Standard
☐ 9" – 11" ☐ 10" – 14" ☐ 12" – 16" ☐ 14" – 18" ☐ 16" – 20"
☐ Fixed Ht. Back without Push Handles (8" – 16") No Charge
☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16"

BACK TYPE: (FOLDING OR OPTIONAL FIXED)

FOLDING BACK: Fold-down back, Locking, Adjustable height/ angle
..... Standard ¹³

SIDE GUARD OPTIONS FOR FOLDING BACK:

- ☐ TED18 Removable Side Guards (with hole) 125.00
☐ 6" tall (sm) ☐ 7.5" tall (lg)
☐ TED79 Removable Side Guards/Fender 175.00
☐ TED121 Fold-in Style Black Aluminum Side Guards 125.00 ¹⁴
☐ TED142 Fold-in Style Carbon Fiber Side Guards 325.00

☐ **OPTIONAL FIXED BACK** No Charge

Fixed Angle Welded at: ☐ 0° ☐ 3° ☐ 6° ☐ 9°
Backrest Crossbar Opening (3.5" std.) Standard

SIDE GUARD OPTIONS FOR FIXED BACK:

- ☐ TED20 Fixed Aluminum Side Guards 175.00 ¹⁵
☐ Even w/Wheels ☐ 1" Above Wheels
☐ TED113 Fixed Aluminum Side Guards with/Fender 225.00 ¹⁶
☐ TED83 Fixed Carbon Fiber Side Guards 250.00 ¹⁷
☐ Even w/Wheels ☐ 1" Above Wheels
☐ TED119 Fixed Carbon Fiber Side Guards w/Fender 300.00

FOOTNOTES

¹⁰ Ergo seat option has the most benefit when seat height differentials are between 2" and 4". Not available with tapered seat width. A minimum rear seat height of 15" is required with 24" rear wheels.

¹¹ An 85 degree front frame angle not available with 6" front casters. 90 degree front frame angle not available with 5" or 6" front casters.

¹² 3/4" tube telescoping in 1" tube. Back height limited to 16" on fixed back style. Fixed height not available on folding back.

¹³ Folding lock-down back is great for transferring the chair in/out of the car.

¹⁴ Fold-in Style side guards flip back for easier transfers and are adjustable in height. Compatible with both high (push or pull to lock) or undermount (scissor) wheels locks. Not available with adjustable height push handles if back height is 10-14" or less.

¹⁵ Fixed back only.

¹⁶ Fixed back only. Not available w/tapered seat width.

¹⁷ Fixed back only. Not available w/tapered seat width.

BACK UPHOLSTERY AND ACCESSORIES:

Ballistic Nylon Adjustable Tension Upholstery with /Straps &

- Standard**
☐ MESH Breathable Mesh Adjustable Tension Upholstery with/Straps
& Anchors No Charge
☐ TED16 Fold Over Back Upholstery 100.00 ¹⁸
☐ TED30 Back Pack 50.00
☐ TED82 Push Handles Adjustable Height Bolt-On 200.00 ¹⁹
☐ TED13 Push Handles, Non-Adjustable Height, Bolt-On ... 75.00
☐ TED90 Welded Integral Push Handles 200.00 ²⁰
☐ TED97 Bolt-Through Integral Push Handles 100.00 ²⁰

FRONT WHEELS/FORK UPGRADE:

- Semi-Pneumatic** ☐ 4" ☐ 5" ☐ 6" Standard ²¹
☐ 3" Micro No Charge
☐ TED165 3" Aluminum Soft Roll 90.00
☐ TED65 Aluminum Soft Roll ☐ 4" ☐ 5" ☐ 6" 90.00 ²¹
☐ TED133 One Sided Forks 150.00 ²²

Wheel/Tire Chart – Applies to Standard & Spinergy® Wheels				
Tire Choices	22"	24"	25"	26"
Primo: 1" high pressure grey clincher	X	X	X	X
Kenda: 1" high pressure clincher		X	X	X
Schwalbe 1" Black flat resistant		X	X	X
1 3/8" grey treaded	X	X		X
1 3/8" w/flat free	X	X		
Solid 1" no-flat		X	X	
Knobby® 2" wide grey		X	X	

REAR WHEEL OPTIONS: (SEE WHEEL/TIRE CHART ABOVE)

- ☐ 22" Sun® Spoke No Charge
☐ 24" Sun Spoke No Charge
☐ 25" Sun Spoke No Charge
☐ 26" Sun Spoke No Charge
☐ TED67 22" Spinergy® Spox Wheels 695.00 ²³
☐ TED67 24" Spinergy Spox Wheels 695.00 ²³
☐ TED134 24" Spinergy LX Wheels *(12 spoke) 895.00
☐ TED67 25" Spinergy Spox Wheels 695.00
☐ TED135 25" Spinergy LX Wheels*(12 spoke) 895.00
☐ TED67 26" Spinergy Spox Wheels 695.00

***Spinergy Spoke Color: (If color is not chosen, spokes will be black.)**

- ☐ SPCBK – Black ☐ SPCY – Yellow ☐ SPCBW - White
☐ SPCR -- Red ☐ SPCBL – Blue

TIRE OPTIONS: (SEE WHEEL/TIRE CHART ABOVE)

- Primo** Standard
☐ Pneumatic – 1 3/8" threaded No Charge
☐ TED123 Schwalbe Flat resistant: black 90.00
☐ TED73 Pneumatic – FF for 22 or 24 x 1 3/8" Tires 125.00 ²⁴
☐ TED60 Knobby Tires (K0067) 75.00
☐ TED7 Solid/Flat Free Tires ☐ Grey ☐ Black 75.00 ²⁵
☐ TED101 Kenda® Tires ☐ Red ☐ Yellow 75.00

FOOTNOTES

¹⁸ TER195, TER197 and UPH1/UPH2 upholstery not available with fold over back upholstery.

¹⁹ Not available with Swingaway tubular padded arms or back heights less than 11" (fixed height TED13 push handles or any of the integral push handles will accommodate Swingaway arms). Push handles are 16" tall and bolted to the back post and have 10" of adjustment up or down.

²⁰ Welded integral push handles only available on Fixed Back with fixed height. Not available with foldover back upholstery.

²¹ 6" casters not available with seat depths of 14" or less, 85 degree front frame angle, 5" and 6" casters not available with 90 degree front frame angle.

Aluminum soft roll casters are 1.4" wide with soft roll tires. For suspension fork option, see page 1. The caster housing is flared 2" on each side for 6" casters and 1.5" for 4" and 5" casters.

²² Not available with 6" casters or suspension. Black anodized color only.

²³ Flat free inserts not available on Spinergy wheels.

²⁴ Not available with 25" wheels.

²⁵ Not available w/heavy duty option.

Invacare®Top End® Terminator™ Everyday/ Titanium Upgrade

Ultra Lightweight Active Performance Rigid Wheelchair

AXLES:

Quick Release.....	Standard
☐ TED27 Quad Release Axles.....	125.00
☐ TED76 Titanium Axles	150.00

HANDRIM OPTIONS:

Aluminum Anodized Handrims	Standard
☐ TED75 Titanium Handrims.....	600.00 ²⁶
☐ TED8 Coated handrims (K0059)	100.00
☐ TED99 Natural Fit® Standard Grip Handrims	295.00 26
☐ TED98 Natural Fit Super Grip Handrims.....	340.00 26
☐ TED199 Natural Fit LT Standard Grip Handrims	325.00 21
☐ TED198 Natural Fit LT Super Grip Handrims.....	370.00 21

MOUNTING HANDRIMS:

Standard Mount Handrim	Standard
☐ Close Mount Handrim.....	No Charge
☐ TER40 Cut Tabs Aluminum handrim	40.00
☐ TERTi40 Cut Tabs Titanium handrim.....	80.00

Overall Width by Camber Angle					
Wheel Type	0°	3°	6°	9°	12°
Spoke	6 ¾"	8"	10"	12 ¼"	14 ¼"
Add the amount above to the seat width to obtain overall width at its widest point					

CAMBER SYSTEM: (CHOOSE STYLE A OR B/ CAMBER) ²⁷

☐ A	Camber Tube with single camber.....	Standard
	Camber: ☐ 0° ☐ 3° ☐ 6° ☐ 9°	
☐ Ted64	Extra Camber Tube	200.00
	☐ 0° ☐ 3° ☐ 6° ☐ 9°	
☐ B	Adjustable Camber System (Choose Two)....	No Charge
	Primary Angle	No Charge
	☐ 0° ☐ 3° ☐ 6° ☐ 9° ☐ 12°	
	Secondary Angle	No Charge
	☐ 0° ☐ 3° ☐ 6° ☐ 9° ☐ 12°	
☐ TED46	Additional Camber Inserts (Choose Two)	125.00
	☐ 0° ☐ 3° ☐ 6° ☐ 9° ☐ 12°	

CAMBER INSERT RELEASE TYPE:

Nut & Bolt	Standard
☐ TED48 Quick Release Lever	30.00

FOOTREST/LEGREST ACCESSORIES:

Calf Strap	Standard
☐ TED28 Frame Protector Set (2 front, 1 rear)	60.00
☐ TED136 Chair Armor Black Neoprene Frame Protectors.....	65.00
☐ TED137 Chair Armor Blue Neoprene Frame Protectors.....	65.00
☐ TED138 Chair Armor Red Neoprene Frame Protectors.....	65.00
☐ TED103 Luggage Carriers	100.00

SEAT ACCESSORIES/ CUSHION:

☐ TED29 Seat Pouch	35.00
☐ TED144 Heavy Duty Seat Pouch.....	100.00
☐ TED72 Positioning Strap – Hook & Loop (K0978)	40.00
☐ TED74 Positioning Strap – Auto Style (K0978)	50.00
☐ TED14 2" Soft Cushion	75.00 ²⁸

SEAT RIGIDIZER/SOLID SEAT:

☐ TEDTi3 Welded titanium rigidizer bars w/ABS solid seat	425.00
☐ TED31 Welded aluminum rigidizer bars w/ABS solid seat	175.00
☐ TED112 Screw-On Seat Terminator Everyday	75.00 ²⁹
☐ TEDTi112 Screw-On Seat Titanium	300.00 29

FOOTNOTES

²⁶ Not available with 22" wheels.

²⁷ Single camber (choice A) is recommended for those who do not change camber. Adjustable camber (choice B) is not available w/heavy duty option.

²⁸ Cushions are filled with 2" soft foam which offer minimum protection. If a therapeutic cushion is needed please order from Invacare Customer Service 1-800-333-6900.

²⁹ Screw-on seat upholstery not available with tapered seat width.

ANODIZING ACCENTS: STD. FORKS AND BACK ANGLE BRACKET

ANOBLK	Bold Black.....	Standard
☐ ANOSLV	Classic Silver	No Charge
☐ TED140	Shadow Blue	50.00
☐ TED141	Victory Red.....	50.00

Arm Heights		
(Top of Seat Rail to Top of Armrest)		
Swingaway	10.5" – 12.5"	(Adjustable in 1" Increments)
T-Arm	9.5"-11.5"	(Adjustable in 1" Increments)

ARM TYPE

☐ TED78	Swingaway Arms (E0973): not available with HD option	175.00 ³⁰
☐ TED17	T-Arm Adjustable Height. Desk (E0973)	175.00 ³¹

WHEEL LOCKS (CHOOSE ONE)

Hideaway Undermount Wheel Locks		Standard
☐	Push or Pull To Lock Wheel Locks.....	No Charge
☐ TED81	Wheel Lock Extensions	39.00
☐ TED22	Hill Holders (E0974)	150.00
☐	Omit Wheel Locks	No Charge

MISCELLANEOUS ACCESSORIES

☐ TED21	Spoke Guard ☐ Clear ☐ Black (K0065).....	75.00
☐ TED23	Rear Anti Tippers (E0971)	100.00
☐ TED26	Crutch Holder	125.00
☐ TED52	Anti Tipper w/3" Swivel Wheel	150.00
☐ TED94	Alignment Gauge.....	75.00
☐ TED145	Water bottle and cage	50.00
☐ RACE22	Top End T-Shirt ☐ L ☐ XL	10.00 ³³

FRAME COLOR AND UPHOLSTERY

UPHOLSTERY: Black.....	Standard
☐ UPH1 Woodland Camouflage	100.00 ³⁴
☐ UPH2 Digital Desert Camouflage	100.00 34

FRAME COLOR:

☐ **Bead Blast Finish for Titanium only**

CONTEMPORARY COLORS

☐ Wet Black	☐ Black Prism	☐ White Snow
☐ Silver Metallic	☐ Stealth Black	

WAY COOL COLORS

☐ Electric Blue	☐ Midnight Blue	☐ Grasshopper Green
☐ Electric Blueberry		

TOO HOT COLORS

☐ Tangerine Red	☐ Electric Red	☐ Grape Madness
☐ Safety Yellow	☐ Sunkissed Orange	☐ Bubblegum Pink

NOTE: CAD drawings are computer aided design that translate the measurements to a 3D preview. Our experienced draftsman completes each CAD request and will send it back for final approval. The charge for RACE101 is \$125 with standard order form selections. A separate purchase order is required to begin a CAD and another purchase order is needed to order the product.

FOOTNOTES

³⁰ Swingaway tubular padded arms interfere with adj. height push handles. Fixed height, bolt-on or integral push handles are compatible with Swingaway Arms.

³¹ Rigid side guards are included with the T-Arms.

³² Push or Pull to Lock Wheel Locks only.

³³ Dealer discount does not apply.

³⁴ Camouflage upholstery not available for screw-on seat, mesh or fold over back upholstery. Camouflage upholstery is ballistic nylon adjustable tension with anchors and straps for the back upholstery. If cushion is ordered, it will also be camouflage.

Top End® Terminator™ Everyday/Titanium Upgrade Ultra Lightweight Active Performance Rigid Wheelchair

	Width & Depth Matrix									
Width (across) Depth (down)	12"	12"	13"	14"	15"	16"	17"	18"	19"	20"
12"	18.43	18.53	18.63	18.73	18.83	18.93	19.03	19.13	19.23	19.33
13"	18.53	18.63	18.73	18.83	18.93	19.03	19.13	19.23	19.33	19.43
14"	18.63	18.73	18.83	18.93	19.03	19.13	19.23	19.33	19.43	19.53
15"	18.73	18.83	18.93	19.03	19.13	19.23	19.33	19.43	19.53	19.63
16"	18.83	18.93	19.03	19.13	19.23	19.33	19.43	19.53	19.63	19.73
17"	18.93	19.03	19.13	19.23	19.33	19.43	19.53	19.63	19.73	19.83
18"	19.03	19.13	19.23	19.33	19.43	19.53	19.63	19.73	19.83	19.93
19"	19.13	19.23	19.33	19.43	19.53	19.63	19.73	19.83	19.93	20.03
20"	19.23	19.33	19.43	19.53	19.63	19.73	19.83	19.93	20.03	

Section 1

Standard Chair Includes:

- * Titanium Frame
- * Tubular footrest w/o plate
- * 4" Semi-pneumatic casters/forks
- * 3° Camber
- * 24" High performance wheels
- * 24" High pressure tires
- * Medium height folding back
- * Adjustable tension seating sling
- * Adjustable tension back upholstery

Top End Terminator Component Weights – Section 2

Seating

Omit adjustable tension seat upholstery	-.85
No cushion	0
2" cushion	+1.10
Screw-on seat upholstery	+.05

Back Type

Folding	0
Fixed	-.80

Push Handles

Adjustable height bolt-on	+1.60
Fixed height bolt-on	+1.70
Integral: welded	+.75
Integral: bolt-through	+.75

Back Upholstery

Adjustable tension w/straps	0
Fold over	+.09
Omit back upholstery	-.97

Armrests

Swing away padded tubular	+2.64
T-Arm, Adj. Height Desk	+6.25

Footrest

Tubular	0
Adjustable angle one piece footplate	+.60
Plastic footrest cover	+.35
Carbon fiber footrest cover	-.30

Front Casters (pair)

3" Micro	-.13
3" Softroll	+.19

4" Semi pneumatic	0
5" Semi pneumatic	+.31
4" Aluminum softroll	+.44
5" Aluminum softroll	+.92

One-sided forks	-.05
-----------------	------

Wheel Locks

Push or pull-to-lock	+1.00
Undermount scissor	+1.25
Wheel lock extensions	+.40
Hill holders/grade aids	+.51

Rear Suspension	+.40
------------------------	------

Sideguards

Fold-in Style carbon fiber	+.95
Fold-in Style black aluminum	+1.45
Folding	+1.45
Removable	+1.45
Removable side guards with fender	+1.94
Fixed aluminum	+1.45
Fixed carbon fiber	+1.00

Rear Tires (pair)

High pressure, Primo or Kenda	0
Schwalbe	+.75
Pneumatic 1 3/8"	+1.40
Pneumatic Flat Free	+1.45
Knobby	+3.10
Solid 1" flat free	+2.15

Rear Wheels (pair)

22" spoke	-.60
24" Sun high performance spoke	0
25" Sun high performance spoke	+.30
26" Sun high performance spoke	+.65
24" Spinergy Spox	-.80
24" Spinergy LX	-1.20
25" Spinergy Spox	-.60
25" Spinergy LX	-1.00
26" Spinergy Spox	-.25

Handrim

Aluminum	0
Titanium	+.20
Plastic coated	+1.40
Natural fit standard or super grip	+1.70
Natural fit LT	+1.25

Axles

Standard quick release	0
Titanium quick release	-.17
Quad quick release	+.22

Accessories

Frame protectors/Chair armor	+.10
Luggage carriers	+.75
Calf strap	+.10
Seat pouch	+.35
Positioning strap – hoop and loop	+.10
Positioning strap – auto style	+.50
Spoke guards – black	+1.40
Spoke guards – clear	+1.60
Anti-tip bars	+1.63
Crutch holder	+1.45

Terminator Everyday will be approximately 2 pounds heavier than Terminator Titanium Upgrade.
Section 1 _____ + Section 2 _____ = Total Weight _____

Invacare Corporation

www.invacare.com



TOP END

Sports and Recreation Products
4501 63rd Circle North
Pinellas Park FL USA
33781-5914

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Sun® is a registered trademark of Sun Components Inc. Natural Fit is a registered trademark of Three Rivers Holding, LLC. Spinergy™ is a registered trademark of Spinergy Corp. Kenda is a registered trademark of Kenda Rubber Industrial Co. Ltd. Schwalbe is a registered trademark of RALF BOHLE GmbH in Reichshof Germany. Specifications and prices are subject to change without notification.
©2015 Invacare Corporation