



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Top End® Terminator™ Titanium/Terminator Everyday Ultra Lightweight Active Performance Rigid Wheelchair

Price List and Order Form

www.topendwheelchair.com

Please fax order to TopEnd at 727-522-1007
For questions call: 800-532-8677 or 727-522-8677

Or order by email: TopEndOrders@invacare.com
To order literature: 800-828-6282



TRANSACTION TYPE:

☐ Order ☐ Quote

Provider Name _____ Date _____ Email _____

Account # _____ P.O. # _____ Phone # _____ Fax # _____

Address _____ City _____ State _____ Zip Code _____

Customer Name _____ Email _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____

Ship order to: ☐ Provider ☐ Customer

Sales Order Number _____

Order will ship to Provider unless otherwise noted.

Internal Use Only

FRAME TYPE

- ☐ TEDTi Top End Terminator™ Titanium \$3,218.00 K0009
Terminator Titanium built of 1" 3-2.5 seamless titanium.
☐ TED Top End Terminator Everyday \$2,490.00 K0009
Terminator Everyday built of 6061T6 Aluminum in 1" tubing only.

HEAVY DUTY OPTION

- ☐ TED61 Heavy Duty Option \$500.00 K0009
350 lbs. weight limit for TED/ 400 lb weight limit for TEDTi.
Restrictions: Only available with folding back, single camber
system (A). TED Heavy Duty will be built of 4130 Chrome-moly
steel. Solid tires, swingaway arms, one-sided forks, anodized
forks and suspension not available.

FRAME MODIFICATIONS

Tubing Size: 1" Standard

- ☐ Tapered Front End No Charge
☐ Non Tapered "V" Front End No Charge
Not available w/tapered seat width or tubular footrest.
☐ Tapered Seat Width No Charge
Tapered seat width tapers near rear to front up to a maximum
difference of 4". The front seat width measures the same as the
outside of the footrest or straight so there is no tapered or V front end.
Not available with Ergo seat or screw-on seat upholstery.
☐ Taper Difference rear-to-front: -1" No Charge
☐ -2" No Charge
☐ -3" No Charge
☐ -4" No Charge

LOWER REAR FRAME

- ☐ Curved Standard
☐ TEDTi5 Square \$100.00

SUSPENSION OPTIONS

Maximum weight capacity is 250 lbs. Minimum rear seat height is 16" and for 26"
or 25" knobby wheels is 17" with rear suspension option.

- ☐ TEDTi4 Rear Suspension Option \$350.00

E1015

SEAT WIDTH

- ☐ W12 12" Seat width No Charge
☐ W13 13" Seat width No Charge
☐ W14 14" Seat width No Charge
☐ W15 15" Seat width No Charge
☐ W16 16" Seat width No Charge
☐ W17 17" Seat width No Charge
☐ W18 18" Seat width No Charge
☐ Other
20" width and wider require HD option.

SEAT DEPTH

- ☐ D14 14" Seat depth No Charge
☐ D15 15" Seat depth No Charge
☐ D16 16" Seat depth No Charge
☐ D17 17" Seat depth No Charge
☐ D18 18" Seat depth No Charge

FRONT SEAT-TO-FLOOR HEIGHT

- ☐ FSTF17 17" Front seat-to-floor height No Charge
☐ FSTF175 17.5" Front seat-to-floor height No Charge
☐ FSTF18 18" Front seat-to-floor height No Charge
☐ FSTF185 18.5" Front seat-to-floor height No Charge
☐ FSTF19 19" Front seat-to-floor height No Charge
☐ FSTF195 19.5" Front seat-to-floor height No Charge
☐ FSTF20 20" Front seat-to-floor height No Charge
☐ FSTF205 20.5" Front seat-to-floor height No Charge
☐ FSTF21 21" Front seat-to-floor height No Charge


Knee-to-Heel Range		
Front Seat-to-Floor	Minimum	Maximum
17"	12 1/4"	16 1/4"
18"	13 1/4"	17 1/4"
19"	14 1/4"	18 1/4"
20"	15 1/4"	19 1/4"
21"	16 1/4"	20 1/4"

Comments

All prices subject to change without notice. All prices in U.S. dollars. All orders F.O.B. Pinellas Park, Florida USA

For the most current pricing information go to www.invacare.com.

Top End® Terminator™ Titanium/Terminator Everyday Ultra Lightweight Active Performance Rigid Wheelchair



REAR-TO-FRONT CASTER DISTANCE

- ☐ Seat Depth plus 3" Standard
☐ Other _____

FRONT FRAME ANGLE

An 85 degree front frame angle not available with 6" front casters. 90 degree front angle not available with 5" or 6" front casters.

- ☐ FFA70 70° No Charge
☐ FFA75 75° No Charge
☐ FFA80 80° No Charge
☐ FFA85 85° No Charge
☐ FFA90 90° No Charge

REAR SEAT-TO-FLOOR HEIGHT


All rear seat heights are without suspension. 14.5" is lowest rear seat height for 24" wheels with Primo, Kenda, or solid tires, 15" for 24" with Schwalbe and 1 3/8" tires and 25" wheels except with knobby, 15.5" for 24" with knobby tires and 16" for 25" knobby or 26" wheels.

- ☐ RSTF145 14.5" Rear Seat-to-Floor Height No Charge
☐ RSTF15 15" Rear Seat-to-Floor Height No Charge
☐ RSTF155 15.5" Rear Seat-to-Floor Height No Charge
☐ RSTF16 16" Rear Seat-to-Floor Height No Charge
☐ RSTF165 16.5" Rear Seat-to-Floor Height No Charge
☐ RSTF17 17" Rear Seat-to-Floor Height No Charge
☐ RSTF175 17.5" Rear Seat-to-Floor Height No Charge
☐ RSTF18 18" Rear Seat-to-Floor Height No Charge
☐ RSTF185 18.5" Rear Seat-to-Floor Height No Charge
☐ RSTF19 19" Rear Seat-to-Floor Height No Charge
☐ RSTF195 19.5" Rear Seat-to-Floor Height No Charge
☐ RSTF20 20" Rear Seat-to-Floor Height No Charge
☐ RSTF205 20.5" Rear Seat-to-Floor Height No Charge
☐ RSTF21 21" Rear Seat-to-Floor Height No Charge

ERGO SEAT OPTION

Ergo seat option has the most benefit when seat height differentials are between 2" and 4". Not available with tapered seat width. A minimum rear seat height of 15" is required with 24" rear wheels. See drawing on pg. 2

- ☐ TED104 Ergo Seat with tapered front end \$170.00
☐ TED110 Ergo Seat with V style front end \$330.00



Foot Space Width

Standard as per diagram above. Exception: seat widths of 12", 13" and 14" will have foot space width of 8".

BACK TYPE

- ☐ LKDOWN Lock-Down, Folding Ht/angle Adj. Back Standard
☐ Fixed Back Angle welded at:
Back height limited to 16".
☐ 0° No Charge
☐ 3° No Charge
☐ 6° No Charge
☐ 9° No Charge
☐ Backrest Crossbar Opening (3.5" std)..... Standard

Comments

All prices subject to change without notice. All prices in U.S. dollars. All orders F.O.B. Pinellas Park, Florida USA

For the most current pricing information go to www.invacare.com.

Top End® Terminator™ Titanium/Terminator Everyday Ultra Lightweight Active Performance Rigid Wheelchair

BACK HEIGHT ADJUSTMENTS

3/4" tubing telescoping in 1" tube. Back height limited to 16" on fixed back style. Fixed height not available on folding back.

Adjustable Height Back without Push Handles:

- ☐ 9" - 11" No Charge
☐ 10" - 14" No Charge
☐ 12" - 16" No Charge
☐ 14" - 18" No Charge
☐ 16" - 20" No Charge

Fixed Height Back without Push Handles (8" - 16"):

- ☐ 8" No Charge
☐ 9" No Charge
☐ 10" No Charge
☐ 11" No Charge
☐ 12" No Charge
☐ 13" No Charge
☐ 14" No Charge
☐ 15" No Charge
☐ 16" No Charge

BACK UPHOLSTERY AND ACCESSORIES

- ☐ CRFBAT Ballistic Nylon Adj. Tnsn Uph w/Straps & Anchors Standard
☐ MESH Breathable Mesh Adj. Tnsn Uph w/Anchors No Charge
☐ TED16 Fold Over Back Upholstery \$100.00
☐ TED30 Back Pack \$50.00
☐ TED82 Push Handles Adjustable Height Bolt-On \$200.00
Not available with Swingaway tubular padded arms or back heights less than 11" (fixed height TED13 push handles or integral push handles will accommodate Swingaway arms). Push handles are 16" tall and bolted to the back post and have 10" of adjustment up or down.
☐ TED13 Push Handles, Non-Adjustable Height, Bolt-On \$75.00

INTEGRAL PUSH HANDLES

- ☐ TED90 Welded Integral push handles \$200.00
Welded integral push handles only available on Fixed Back, with fixed height. Not available with foldover back upholstery.
☐ TED97 Bolt-Through Integral Push Handles \$100.00

FRONT WHEELS/FORKS

- ☐ TED133 Frog Legs® One Sided Forks \$150.00
Not available with 6" casters or suspension. Black anodized color only.
☐ 3" Micro No Charge
☐ TED165 3" Aluminum Soft Roll \$90.00
☐ Semi-Pneumatic No Charge
6" casters not available with seat depths of 14" or less, 85 degree front frame angle, 5" and 6" casters not available with 90 degrees front frame angle. For suspension fork option see page 1. The caster housing is flared 2" on each side for 6" casters and 1.5" for 4" and 5" casters.
☐ 4" See Above
☐ 5" See Above
☐ 6" See Above
☐ TED65 Aluminum Soft Roll \$90.00
6" casters not available with seat depths of 14" or less, 85 degree front frame angle. Aluminum softroll casters are 1.4" wide with soft roll tires. For suspension fork option see page 1. The caster housing is flared 2" on each side for 6" casters and 1.5" for 4" and 5" casters.
☐ 4" See Above
☐ 5" See Above
☐ 6" See Above
☐ TED88 Lighted Front Wheel \$70.00
☐ 3" See Above
☐ 4" See Above
☐ 5" See Above

REAR WHEEL OPTIONS

Sun Performance Rear Wheels Standard

- ☐ 22" Sun® Spoke No Charge
☐ 24" Sun Spoke No Charge
☐ 25" Sun Spoke No Charge
☐ 26" Sun Spoke No Charge

Spinergy™ Spox \$695.00

- ☐ TED67 22" Spinergy Spox Wheels See above
☐ TED67 24" Spinergy Spox Wheels See above
☐ TED67 25" Spinergy Spox Wheels See above
☐ TED67 26" Spinergy Spox Wheels See above

Spinergy LX \$895.00

- ☐ TED134 24" Spinergy LX Wheels (12 spoke) See above
☐ TED135 25" Spinergy LX Wheels (12 spoke) See above

☐ SPC Choose Spinergy Spoke Color No Charge

- ☐ Black Standard
☐ Yellow No Charge
☐ White No Charge
☐ Red No Charge
☐ Blue No Charge

AXLES

- ☐ Quick Release axles Standard
☐ TRA27 Quad Release axles \$125.00
☐ TRA77 Titanium axles \$150.00

WHEELS: EXTRA SET FOR OFF ROAD

24" knobby match up with 25" daily use wheels with Primo tires and 25" knobby match up best with 26". Extra set includes standard wheels, anodized pushrims and knobby tires.

- ☐ TED91 24" Knobby Tires (540) \$550.00
☐ TED92 25" Knobby Tires (559) \$550.00
☐ TED93 Extra Set Of Axles \$55.00

Wheel/Tire Chart (Applies to Standard and Spinergy Wheels)

Tire Choices	22"	24"	25"	26"
Primo: 1" high pressure grey clincher	X	X	X	X
Kenda: 1" high pressure clincher		X	X	X
Schwalbe 1" Black flat resistant		X	X	X
1 - 3/8" grey treaded	X	X		X
1 - 3/8" w/ flat free	X	X		
Solid 1" no-flat		X	X	
Knobby: 2" wide grey		X	X	

Comments

All prices subject to change without notice. All prices in U.S. dollars. All orders F.O.B. Pinellas Park, Florida USA

For the most current pricing information go to www.invacare.com.

Top End® Terminator™ Titanium/Terminator Everyday Ultra Lightweight Active Performance Rigid Wheelchair

TIRE OPTIONS (SEE WHEEL/TIRE CHART ON PAGE 2)

- ☐ Primo No Charge
- ☐ TED123 Schwalbe Flat resistant: black \$90.00
- ☐ Pneumatic – 1 3/8" threaded No Charge
- ☐ TED73 Pneumatic – FF for 22 or 24 x 1 3/8" Tires \$125.00 K0064
Not available with Spinergy wheels.
- ☐ TED60 Knobby Tires \$75.00 K0067
- ☐ TED7 Solid/Flat Free Tires \$75.00
Not available with Heavy-Duty option.
- ☐ ☐ Grey See above
- ☐ ☐ Black See above
- ☐ TED101 Kenda® Tires \$75.00
- ☐ ☐ Red See above
- ☐ ☐ Yellow See above

MOUNTING HANDRIMS

- ☐ OUT Standard Mount Handrim Standard
- ☐ IN Close Mount Handrim No Charge
- ☐ TER40 Cut Tabs Aluminum pushrim \$40.00
- ☐ TERT140 Cut Tabs Titanium pushrim \$80.00

HANDRIM OPTIONS

- ☐ Aluminum Anodized Handrims Standard
- ☐ TED75 Titanium Handrims \$600.00
Not available with 22" wheels.
- ☐ TED8 Coated handrims \$100.00 K0059
- ☐ TED99 Natural Fit® Standard Grip Handrims \$295.00
Not available with 22" wheels.
- ☐ TED98 Natural Fit Super Grip Handrims \$340.00
Not available with 22" wheels.
- ☐ TED199 Natural Fit LT Standard Grip Handrims \$325.00
Not available with 22" wheels.
- ☐ TED198 Natural Fit LT Super Grip Handrims \$370.00
Not available with 22" wheels.

CAMBER SYSTEM (CHOOSE A OR B)

Single camber is recommended for those who do not change camber.
Adjustable camber is not available with heavy-duty option.

Single Camber with Camber Tube Standard

- ☐ 0° See above
- ☐ 3° See above
- ☐ 6° See above
- ☐ 9° See above

☐ Ted64 Extra Camberbe \$200.00

- ☐ 0° See above
- ☐ 3° See above
- ☐ 6° See above
- ☐ 9° See above

Adj. Camber System (Choose 1 Primary and 1 Secondary) No Charge

Primary Angle:

- ☐ 0° See above
- ☐ 3° See above
- ☐ 6° See above
- ☐ 9° See above
- ☐ 12° See above

Secondary Angle:

- ☐ 0° No Charge
- ☐ 3° See above
- ☐ 6° See above
- ☐ 9° See above
- ☐ 12° See above

☐ TED46 Additional Camber Inserts (Choose Two) \$125.00

- ☐ 0° See above
- ☐ 3° See above
- ☐ 6° See above
- ☐ 9° See above
- ☐ 12° See above

Overall Width by Camber Angle					
Degrees	0°	3°	6°	9°	12°
Seat width plus	6 3/4"	8"	10"	12 1/4"	14 1/4"
Add the amount above to the seat width to obtain overall width at its widest point.					

CAMBER INSERT RELEASE TYPE

- ☐ Nut and Bolt Standard
- ☐ TED48 Quick Release Lever \$30.00

FOOTREST OPTIONS

- ☐ 1 Piece Tubular Aluminum Footrest Standard
One piece tubular footrest is aluminum and does not have cover. Either a plastic or carbon fiber footrest cover is available. Tubular footrest not available with V style front end. For customers with short leg lengths, please indicate top of footplate to floor measurement. See Knee-to-heel range length chart page 1 and consider option TED191.
- ☐ TED96 1 Piece Tubular Titanium Footrest Upgrade \$200.00
Top of Footplate to Floor _____
- ☐ TED10 Tubular Footrest with Plastic Footrest Cover \$75.00
- ☐ TED89 Tubular Footrest with Carbon Fiber Footrest Cover \$125.00
- ☐ Adjustable Angle One Piece Footplate No Charge
Non-tapered V-style front end requires adjustable angle one-piece footplate.
- ☐ TED191 Adj Angle Footplate with 2" Risers \$95.00
Top of Footplate to Floor _____

FOOTREST/LEGREST ACCESSORIES

Frame Protector Set includes 2 front, 1 rear and Chair Armor includes 2 front.

- ☐ TED28 Frame Protector Set (2 front, 1 rear) \$60.00
- ☐ TED136 Chair Armor Black Neoprene Frame Protectors \$65.00
- ☐ TED137 Chair Armor Blue Neoprene Frame Protectors \$65.00
- ☐ TED138 Chair Armor Red Neoprene Frame Protectors \$65.00
- ☐ TED103 Luggage Carriers \$100.00
- ☐ Calf Strap Standard

SEAT ACCESSORIES

- ☐ TED29 Seat Pouch \$35.00
- ☐ TED144 Heavy Duty Seat Pouch \$100.00
- ☐ TED72 Positioning Strap – Hook and Loop \$40.00 K0978
- ☐ TED74 Positioning Strap – Auto Style \$50.00 K0978

CUSHIONS

- ☐ TED14 2" Soft Cushion \$75.00
Cushions are made with 2" soft foam which offer minimum protection. If a therapeutic cushion is needed please order from Invacare Customer service 1-800-333-6900.

Comments

All prices subject to change without notice. All prices in U.S. dollars. All orders F.O.B. Pinellas Park, Florida USA

For the most current pricing information go to www.invacare.com.

Top End® Terminator™ Titanium/Terminator Everyday Ultra Lightweight Active Performance Rigid Wheelchair

SEAT RIGIDIZER/SOLID SEAT

- ☐ TEDT13 Welded titanium rigidizer bars with ABS solid seat .. \$425.00
☐ TED31 Welded aluminum rigidizer bars with ABS solid seat \$175.00
☐ TED112 Screw-On Seat Terminator Everyday \$75.00
Screw-on seat upholstery not available w/tapered seat width.
☐ TEDTi112 Screw-On Seat Titanium..... \$300.00

ANODIZING ACCENTS (STD. FORKS AND BACK ANGLE BRACKET)

- ☐ ANOBLK Bold Black Standard
☐ ANOSLV Classic Silver No Charge
☐ TED140 Shadow Blue \$50.00
☐ TED141 Victory Red \$50.00

ARM TYPE

- ☐ TED78 Swingaway Arms \$175.00 E0973
Swingaway tubular padded arms interfere w/adj. height push handles. Fixed height, bolt-on or integral push handles are compatible with Swingaway arms.
☐ TED17 T-Arm Adj. Hgt. Desk \$175.00 E0973
Rigid side guards included with the T-arms.

Arm Heights (Top of Seat Rail to Top of Armrest)		
Swingaway	10 1/2" – 12 1/2"	Adjustable in 1" Increments
T-Arm	9 1/2" - 11 1/2"	Adjustable in 1" Increments

ARM ACCESSORIES

Not available with T-Arms.

- ☐ TED18 Removable Side Guards: 6" tall (small) \$125.00
☐ TED18 Removable Side Guards: 7.5" tall (large) \$125.00
☐ TED79 Removable Side Guards/Fender \$175.00
☐ TED20 Fixed Aluminum Side Guards: Even with Wheels \$175.00
Fixed back only.
☐ TED20 Fixed Aluminum Side Guards: 1" Above Wheels \$175.00
Fixed back only.
☐ TED113 Fixed Aluminum Side Guards w/Fender \$225.00
Fixed back only.
☐ TED83 Fixed Carbon Fiber Side Guards: Even with Wheels . \$250.00
Fixed back only. Not available w/tapered seat width.
☐ TED83 Fixed Carbon Fiber Side Guards: 1" Above Wheels .. \$250.00
Fixed back only. Not available w/tapered seat width.
☐ TED119 Fixed Carbon Fiber Side Guards with Fender \$300.00
Fixed back only. Not available w/tapered seat width.
☐ TED121 Fold-in Style Black Aluminum Side Guards \$125.00
Fold-in Style sideguards flip back for easier transfers and are adjustable in height. Compatible with both high (push or pull to lock) or undermount (scissor) wheels locks. Not available with adjustable height push handles if back height is 10-14" or less.
☐ TED142 Fold-in Style Carbon Fiber Side Guards \$325.00
Fold-in Style sideguards flip back for easier transfers and are adjustable in height. Compatible with both high (push or pull to lock) or undermount (scissor) wheels locks. Not available with adjustable height push handles if back height is 10-14" or less.

WHEEL LOCKS (Choose One)

- ☐ Hideaway Undermount Wheel Locks No Charge
☐ Push or Pull To Lock Wheel Locks No Charge
☐ TED81 Wheel Lock Extensions \$39.00
Push or Pull to Lock Wheel locks only.
☐ TED22 Hill Holders \$150.00 E0974
☐ Omit Wheel Locks No Charge

MISCELLANEOUS ACCESSORIES

- ☐ TED21 Spoke Guard \$75.00 K0065
☐ Clear See Above
☐ Black See Above
☐ TED23 Rear Anti Tippers \$100.00 E0971
☐ TED26 Crutch Holder \$125.00
☐ TED52 Anti Tipper with 3" Swivel Wheel \$150.00
☐ TED94 Alignment Gauge \$75.00
☐ TED145 Water bottle and cage \$50.00
☐ RACE22 Top End T-Shirt \$10.00
Dealer discount does not apply.
☐ Large See Above
☐ Extra Large See Above

FRAME COLOR AND UPHOLTERY

UPHOLSTERY

- ☐ Black Standard
☐ UPHI Woodland Camouflage \$100.00
☐ UPH2 Digital Desert Camouflage \$100.00
Camouflage upholstery not available for screw-on seat, mesh or fold over back upholstery. Camouflage upholstery is ballistic nylon adjustable tension with anchors and straps for the back upholstery. If cushion is ordered it will also be camouflage.

FRAME FINISH

- ☐ Bead Blast Finish for Titanium No Charge

CONTEMPORARY COLORS

- ☐ 24P Wet Black No Charge
☐ 115P Black Prism No Charge
☐ 29P White Snow No Charge
☐ 60P Silver Metallic No Charge
☐ 146P Stealth Black No Charge

WAY COOL COLORS

- ☐ 62P Electric Blue No Charge
☐ 105P Midnight Blue No Charge
☐ 127P Grasshopper Green No Charge

TOO HOT COLORS

- ☐ 117P Tangerine Red No Charge
☐ 120P Grape Madness No Charge
☐ Raspberry Red No Charge
☐ 132P Safety Yellow No Charge
☐ 116P Sunkissed Orange No Charge
☐ 119P Bubblegum Pink No Charge

Comments

All prices subject to change without notice. All prices in U.S. dollars. All orders F.O.B. Pinellas Park, Florida USA

For the most current pricing information go to www.invacare.com.

**Top End® Terminator™ Titanium/Terminator Everyday
Ultra Lightweight Active Performance Rigid Wheelchair**

Component Weights				
Seating			Fold-in Style black aluminum	+1.45
Omit adjustable tension seat upholstery	-.85		Folding	+1.45
No cushion	0		Removable	+1.45
2" cushion	+1.10		Removable side guards with fender	+1.94
Screw-on seat upholstery	+.05		Fixed aluminum	+1.45
Back Type			Fixed carbon fiber	+1.00
Folding	0		Rear Tires (pair)	
Fixed	-.80		High pressure, Primo or Kenda	0
Push Handles			Schwalbe	+.75
Adjustable height bolt-on	+1.60		Pneumatic 1 3/8"	+1.40
Fixed height bolt-on	+1.70		Pneumatic Flat Free	+1.45
Integral: welded	+.75		Knobby	+3.10
Integral: bolt-through	+.75		Solid 1" Flat Free	+2.15
Back Upholstery			Rear Wheels (pair)	
Adjustable tension w/ straps	0		22" spoke	-.60
Fold over	+.09		24" Sun high performance spoke	0
Omit back upholstery	-.97		25" Sun high performance spoke	+.30
Armrests			26" Sun high performance spoke	+.65
Swingaway padded tubular	+2.64		24" Spinergy Spox	-.80
T-Arm, Adjustable Height Desk	+6.25		24" Spinergy LX	-1.20
Footrest			25" Spinergy Spox	-.60
Tubular	0		25" Spinergy LX	-1.00
Adjustable angle one piece footplate	+.60		26" Spinergy Spox	-.25
Plastic footrest cover	+.35		Handrim	
Carbon fiber footrest cover	-.30		Aluminum	0
Front Casters (pair)			Titanium	+.20
3" Micro	-.13		Plastic coated	+1.40
3" Softroll	+.19		Natural fit standard or super grip	+1.70
3" Micro lighted	-.13		Natural fit LT	+1.25
4" Semi pneumatic	0		Axles	
5" Semi pneumatic	+.31	Standard quick release	0	
4" Aluminum softroll	+.44	Titanium quick release	+.17	
5" Aluminum softroll	+.92	Quad quick release	+.22	
4" Lighted	+.22	Accessories		
5" Lighted	+.53	Frame protectors/Chair armor	+.10	
		Luggage carriers	+.75	
		Calf strap	+.10	
Wheel Locks		Seat pouch	+.35	
Push or pull-to-lock	+1.00	Positioning Strap – hook and loop	+.10	
Undermount scissor	+1.25	Positioning Strap – auto style	+.50	
Wheel Lock extensions	+.40	Spoke guards – black	+1.40	
Hill holders/grade aids	+.51	Spoke guards – clear	+1.60	
Rear Suspension	+.40	Anti-tip bars	+1.63	
Sideguards		Crutch holder	+1.45	
Fold-in Style carbon fiber	+.95	Weights are approximate +/- .50 lbs.		

All prices subject to change without notice. All prices in U.S. dollars. All orders F.O.B. Pinellas Park, Florida USA

For the most current pricing information go to www.invacare.com.

Top End® Terminator™ Titanium/Terminator Everyday Ultra Lightweight Active Performance Rigid Wheelchair

Width and Depth Weight Chart									
	D12	D13	D14	D15	D16	D17	D18	D19	D20
W12	18.43	18.53	18.63	18.73	18.83	18.93	19.03	19.13	19.23
W13	18.53	18.63	18.73	18.83	18.93	19.03	19.13	19.23	19.33
W14	18.63	18.73	18.83	18.93	19.03	19.13	19.23	19.53	19.43
W15	18.73	18.83	18.93	19.03	19.13	19.23	19.53	19.63	19.53
W16	18.83	18.93	19.03	19.13	19.23	19.53	19.63	19.73	19.63
W17	18.93	19.03	19.13	19.23	19.53	19.63	19.73	19.83	19.73
W18	19.03	19.13	19.23	19.53	19.63	19.73	19.83	19.93	19.83
W19	19.13	19.23	19.53	19.63	19.73	19.83	19.93	20.03	19.93
W20	19.23	19.53	19.63	19.73	19.83	19.93	20.03	20.13	20.03

Total Weight	
Weight from Width/Depth Chart	
Total Weight from Components	
Total Weight	
Note: Terminator Everyday will be approximately 2 lbs. heavier than Terminator Titanium.	

Standard Chair Includes

- Titanium Frame
- Tubular footrest without plate
- 4" Semi-pneumatic casters/forks
- 3° Camber
- 24" High performance wheels
- 24" High pressure tires
- Medium height folding back
- Adjustable tension seating sling
- Adjustable tension back upholstery

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation

TOP END

Sports and Recreation Products
4501 63rd Circle North
Pinellas Park FL USA
33781-5914

www.invacare.com

©20010 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols TM and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Sun® is a registered trademark of Sun Components Inc. Natural Fit is a registered trademark of Three Rivers Holding, LLC. SpinergerTM is a registered trademark of Spinerger Corp. Frog Legs is a registered trademark of Frog Legs, Inc. Kenda is a registered trademark of Kenda Rubber Industrial Co. Ltd. Schwalbe is a registered trademark of RALF BOHLE GmbH in Reichshof Germany. Specifications and prices are subject to change without notification.