INVACARE® 9000 XT SERIES WHEELCHAIRS





The Invacare® 9000 XT Wheelchair is a high performance, low maintenance and low total lifetime cost product. The 9000 XT wheelchair is simply the most versatile lightweight wheelchair available from Invacare, boasting an extensive offering of frame styles, seat-to-floor heights, back styles, widths, depths, frame colors and a wide range of options and accessories. The 9000 XT wheelchair components are interchangeable with all other 9000 series and Invacare® Tracer® series chairs.

The **Invacare® 9000 XT Recliner Wheelchair** is the number one choice for individuals who demand a truly versatile, yet moderately priced, lightweight wheelchair. The 9000 XT Recliner wheelchair can be configured into an endless number of seat widths, depths and heights, and has multiple frame and upholstery options. Components and accessories are interchangeable with the other 9000 series chairs.

Invacare® 9000 XT Wheelchair Features

Dual axle positions allow variation of seat-to-floor heights

Offset axle position allows adjustment of center-of-gravity to increase maneuverability or stability to meet each user's skill level

Urethane rear tires mounted on "no flex" wheels offer superior performance

Urethane casters for better ride and lighter weight

Invacare® 9000 XT Recliner Wheelchair Features

Dynamic recline range from 95° to 175°, +/- 5°

Rear axle base moves, providing stable center-of-gravity in recline position

Multiple wheel offerings allow variation of seat-to-floor heights

Urethane rear tires mounted on "no flex" wheels offer superior performance

Urethane casters for better ride and lighter weight Low-maintenance, dual-sealed precision bearings

Heavy-duty inner liners keep seat and back from stretching

Stocked with attendant (standard) or patient-operated wheel locks

HCPCS Code: K0004 & E1226 (recliner)











- 1. 10" Removable headrest
- 2. Dual axle positions
- 3. Urethane casters

Model#

9XT

- 4. Seating options available
- 5. Builder options available



Invacare Corporation

www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668 - 5324

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 11-227 Rev. 04/14 131068



Overall

Width

21", 22", 23"-29"



Overall

Height

36"

52"



Overall Depth

w/o riggings

28"

32"



Folded

Width

12"

17"









	H			
Seat-to-	Back			
loor Height	Height			
13 5"-19 5"	15"-19"			













Seat

Dimensions

14"-20", 22"

x 16"-18"

14"-22" (2" inc)



17.5"-19.5"



34"

	Caster Height	Wheel Height	Weight w/o Rigging	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
9XT	6", 8"	20", 22", 24"	31 lb.	250 lb.	39 lb.	32"x 13"x 37"	5 Years
9RC	8"	20", 22", 24"	42 lb.	250 lb.	50 lb.	32"x 13"x 37"	5 Years



All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.