




VECTOR

HCPCS CODING: E2624/E2625
ADJUSTABLE SKIN PROTECTION & POSITIONING



PATENTED VICAIR® TECHNOLOGY

Add or remove "Comfort Cells" as needed to customize the cushion to fit each person's individual needs. Newly designed comfort cells are softer and quieter. This gives the cushion a softer feel and allows the user to immerse into the cushion more.



PRODUCT FEATURES

Adjustable lateral leg contours promote postural alignment.

Allows for deep immersion to redistribute pressure and maintain sitting stability.



Air Cells may be added or removed from medial thigh separator.

Split elevated side compartments running front to back provide independently adjustable lateral pelvic and upper leg support.

- Eight fully adjustable compartments provide custom fit with maintenance free air cells.
- Newly designed Comfort Cells are softer and quieter, allowing the user to immerse into the cushion more than ever before.

COVER OPTIONS



STRETCH-AIR FABRIC: designed for air transmission and heat reduction. This breathable laminate combines a high-stretch polyester top layer with a breathable vapor barrier.



COMFORT-TEK FABRIC: the fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multi-dimensional stretch and helps alleviate pressure on bony prominences.

VECTOR LOW PROFILE WITH VICAIR TECHNOLOGY

stretch-air	comfort-tek	w x l	front	rear	abductor	leg trough	wt. limit
VT-SL-1616	VT-FL-1616	16 x 16	2.5"	2.5"	3"	2"	550 lbs
VT-SL-1618	VT-FL-1618	16 x 18	2.5"	2.5"	3"	2"	550 lbs
VT-SL-1816	VT-FL-1816	18 x 16	2.5"	2.5"	3"	2"	550 lbs
VT-SL-1818	VT-FL-1818	18 x 18	2.5"	2.5"	3"	2"	550 lbs
VT-SL-1820	VT-FL-1820	18 x 20	2.5"	2.5"	3"	2"	550 lbs
VT-SL-2016	VT-FL-2016	20 x 16	2.5"	2.5"	3"	2"	550 lbs
VT-SL-2018	VT-FL-2018	20 x 18	2.5"	2.5"	3"	2"	550 lbs

VECTOR

WITH VICAIR® TECHNOLOGY

stretch-air	comfort-tek	w x l	front	rear	abductor	leg trough	wt. limit
VT-S-1010	VT-F-1010	10 x 10	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-1212	VT-F-1212	12 x 12	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-1414	VT-F-1414	14 x 14	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-1416	VT-F-1416	14 x 16	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-1616	VT-F-1616	16 x 16	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-1618	VT-F-1618	16 x 18	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-1620	VT-F-1620	16 x 20	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-1816	VT-F-1816	18 x 16	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-1818	VT-F-1818	18 x 18	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-1820	VT-F-1820	18 x 20	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-2016	VT-F-2016	20 x 16	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-2018	VT-F-2018	20 x 18	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-2020	VT-F-2020	20 x 20	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-2216	VT-F-2216	22 x 16	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-2218	VT-F-2218	22 x 18	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-2220	VT-F-2220	22 x 20	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-2416	VT-F-2416	24 x 16	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-2418	VT-F-2418	24 x 18	3.5"	3.5"	4.5"	3"	550 lbs
VT-S-2420	VT-F-2420	24 x 20	3.5"	3.5"	4.5"	3"	550 lbs

NOTES

+ CUSHIONS AVAILABLE WITH RIGID INSERTS & EXTRA INCONTINENT PROOF LINER. SEE PAGE 19 FOR DETAILS