ADJUSTER

HCPCS CODING: E2622/E2623 ADJUSTABLE SKIN PROTECTION





- Five fully adjustable compartments provide custom fit with maintenance free air cells.
- Newly designed Comfort Cells are softer and quieter, allowing the user to immerse into the cushion more than ever before.

ADJUSTER WITH VICAIR® TECHNOLOGY								
stretch-air part #	comfort-tek part #	w x l	front	rear	wt. limit			
AJ-S-1414	AJ-F-1414	14 x 14	3.5"	3.5"	550 lbs			
AJ-S-1416	AJ-F-1416	14 x 16	3.5"	3.5"	550 lbs			
AJ-S-1616	AJ-F-1616	16 x 16	3.5"	3.5"	550 lbs			
AJ-S-1618	AJ-F-1618	16 x 18	3.5"	3.5"	550 lbs			
AJ-S-1620	AJ-F-1620	16 x 20	3.5"	3.5"	550 lbs			
AJ-S-1816	AJ-F-1816	18 x 16	3.5"	3.5"	550 lbs			
AJ-S-1818	AJ-F-1818	18 x 18	3.5"	3.5"	550 lbs			
AJ-S-1820	AJ-F-1820	18 x 20	3.5"	3.5"	550 lbs			
AJ-S-2016	AJ-F-2016	20 x 16	3.5"	3.5"	550 lbs			
AJ-S-2018	AJ-F-2018	20 x 18	3.5"	3.5"	550 lbs			
AJ-S-2020	AJ-F-2020	20 x 20	3.5"	3.5"	550 lbs			
AJ-S-2216	AJ-F-2216	22 x 16	3.5"	3.5"	550 lbs			
AJ-S-2218	AJ-F-2218	22 x 18	3.5"	3.5"	550 lbs			
AJ-S-2220	AJ-F-2220	22 x 20	3.5"	3.5"	550 lbs			
AJ-S-2416	AJ-F-2416	24 x 16	3.5"	3.5"	550 lbs			
AJ-S-2418	AJ-F-2418	24 x 18	3.5"	3.5"	550 lbs			
AJ-S-2420	AJ-F-2420	24 x 20	3.5"	3.5"	550 lbs			

NOTES

+ CUSHIONS AVAILABLE WITH RIGID INSERTS & EXTRA INCONTINENT PROOF LINER. SEE PAGE 19 FOR DETAILS





PATENTED VICAIR® TECHNOLOGY_

Add or remove "Comfort Cells" as needed to customize the cushion to fit each person's individual needs. Newly designed comfort cells are softer and quieter. This gives the

cushion a softer feel and allows the user to immerse into the cushion more.

COVER OPTIONS

STRETCH-AIR FABRIC: designed for air transmission and heat reduction. This breathable laminate combines a high-stretch polyester top layer with a breathable vapor barrier.

COMFORT-TEK FABRIC: the fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multidimensional stretch and helps alleviate pressure on bony prominences.

ADJUSTER LOW PROFILE WITH VICAIR TECHNOLOGY									
stretch-air part #	comfort-tek part #	w x l	front	rear	wt. limit				
AJ-SL-1616	AJ-FL-1616	16 x 16	2.5"	2"	550 lbs				
AJ-SL-1618	AJ-FL-1618	16 x 18	2.5"	2"	550 lbs				
AJ-SL-1816	AJ-FL-1816	18 x 16	2.5"	2"	550 lbs				
AJ-SL-1818	AJ-FL-1818	18 x 18	2.5"	2"	550 lbs				
AJ-SL-1820	AJ-FL-1820	18 x 20	2.5	2"	550 lbs				
AJ-SL-2016	AJ-FL-2016	20 x 16	2.5"	2"	550 lbs				
AJ-SL-2018	AJ-FL-2018	20 x 18	2.5"	2"	550 lbs				