

Invacare® Top End® Force™ 3 Handcycle



Product Description:

The Invacare® Top End® Force™-3 Handcycle is in stock and ready to ship. It features a wrap-around frame design and a cable driven parking brake and upper chain ring derailleur all at a great price! This is the perfect handcycle for customers looking to just get started in handcycling.

About New Smart Forms:

Top End is excited to introduce a new way to order Top End products. Our new Smart Forms provide an easy and electronic way to fill out order forms, save copies, and submit them via e-mail. Please see below for the features and useful information.

How to use:

- Use any Adobe Acrobat application to electronically fill out the order form
- Use the default selections to order a product with the most common configuration
- Use the **"Reset"** button located in the forms footer to reset selections to their defaults
- Use the **"Clear All" button** located in the forms footer to clear selections, including defaults. This feature allows for a blank order form to be printed and completed by hand
- Use the **"E-mail to Top End order box"** button located on the last page of the form to e-mail the completed form directly to the topendorders@invacare.com e-mail inbox
- Forms can be completed electronically, printed, and e-mailed or faxed to topendorders@invacare.com or 727-522-1007 if access to Acrobat Reader for Windows is unavailable

(Acrobat Reader for Windows with java-script enabled is required to save and attach to e-mail)

Features:

- Tabulated options for an organized and standardized work-flow
- New more specific option codes which will appear on confirmations to improve customer / manufacturing communication
- Common special requests have been given option codes
- Built in exclusions to help eliminate errors
- Ability to save the form after it is filled out
- Ability to e-mail the form inside Acrobat Reader rather than saving a copy and attaching to an e-mail
- Ability to clear the entire form so it can be printed and completed by hand in the field when a computer is not available
- Ability to reset form to default selections

Invacare® Top End® Force™ 3 Handcycle



PROVIDER / CUSTOMER INFORMATION:			
Provider Information:		End User Information:	
Date:	Name:	Height:	
IVC Acct #:	Ship to Address:	Weight (lbs.):	
Contact:	Address Line 2:	Leg Length (in.):	
Phone #:	City:	Hip Width (in.):	
Fax/Email:	State:	Zip:	Phone #:
<input type="radio"/> Ship to Provider		<input type="radio"/> Ship to End User	
		Email:	

ORDER REQUEST TYPE:	
<input type="radio"/> Quote	<input type="radio"/> Order
PO #	<input type="text"/>

BASE MODEL:		
<input type="radio"/> 10039	Force 3 Handcycle	\$3,499

Assembly: <i>(see definitions)</i>		
<input type="radio"/> 10046	Unassembled	\$0
<input type="radio"/> 10047	Assembled	\$250
Shipping:		
Shipping Fee		\$150
<i>(Continental US Only - See definitions)</i>		

BACKREST TYPE:		
<input type="radio"/> 15172	Straight, Upright position, Adjustable Angle Backrest	\$0

FORK TYPE:		
<input type="radio"/> 75026	Wide Fork	\$0

DRIVE CONTROL:		
CRANK ASSEMBLY AND FORK:		
<input type="radio"/> 65179	Length and Height Adjustable Fork	\$0
CRANK ARM TYPE:		
<input type="radio"/> 65048	STD Crank 170mm Long X 12-1/4" Wide	\$0
HAND PEDALS:		
<input type="radio"/> 65137	Medium, 1-3/4" Shifter Post on Right	\$0
SHIFTING:		
<input type="radio"/> 65230	SRAM GX 10 Speed Shifter	\$0
CHAINRING:		
<input type="radio"/> 65104	48-36-26 Tooth Upper Chainring	\$0
CASSETTE:		
<input type="radio"/> 65161	11-36 Tooth Front Wheel Cassette	\$0
BRAKING:		
<input type="radio"/> 65090	Locking Parking Brake	\$0

FOOTREST:		
<input type="radio"/> 60005	5.5" Footrest Width	\$0

WHEELBASE: <i>(Front to Rear Axle)</i>		
<input type="radio"/> 75078	60" Wheelbase	\$0

CAMBER ANGLE:		
<input type="radio"/> 45025	1 Degree Camber	\$0



WHEEL SETS:		
WHEEL SIZE:		
<input type="radio"/> 50192	26" (559)	\$0
WHEEL TYPE:		
<input type="radio"/> 50196	High Performance	\$0
TIRES:		
<input type="radio"/> 50248	Cruiser Clincher Tires	\$0
<input type="radio"/> 50215	Knobby Tires <i>(Self Installation Required)</i>	\$80
AXLES:		
<input type="radio"/> 50402	Threaded Axles	\$0

FRAME COLOR:			
<input type="radio"/> 55017	Black	Stealth Black	\$0

UPHOLSTERY:		
MATERIAL & COLOR:		
<input type="radio"/> 30001	Nylon, Black Upholstery	\$0
SEAT UPHOLSTERY:		
<input type="radio"/> 25060	Hook and Loop Positioning Strap	\$0
BACK UPHOLSTERY:		
<input type="radio"/> 25064	Back Cushion	\$0
SEAT CUSHION:		
<input type="radio"/> 25063	Seat Cushion	\$0

Invacare® Top End® Force™ 3 Handcycle



ACCESSORIES:		
STANDARD:		
<input type="radio"/> 25061	Safety Flag	\$0
<input type="radio"/> 25072	Chain Guard	\$0
<input type="radio"/> 25067	Drafting Bumper	\$0
UPGRADABLE:		
CRANK UPGRADES:		
<input type="radio"/> 65053	V-Crank Upgrade 170mm X 20" Cranks with a 48-36-26 Chainring	\$400
SAFETY:		
<input type="radio"/> 25030	Extra Steering Dampener	\$15
<input type="radio"/> 25015	Light, Front <i>(White)</i>	\$15
<input type="radio"/> 25016	Light, Rear <i>(Red)</i>	\$15
<input type="radio"/> 25041	Mirror, Footrest Mounted	\$35
<input type="radio"/> 25018	Helmet, Small/Medium	\$30
<input type="radio"/> 25019	Helmet, Medium/Large	\$30
<input type="radio"/> 25029	Tool and Tire Repair Kit	\$95
HYDRATION:		
<input type="radio"/> 25021	Water Bottle and Cage	\$50
POSITIONING:		
<input type="radio"/> 25001	Click Strap, 8" Small	\$75
<input type="radio"/> 25002	Click Strap, 10" Medium	\$80
<input type="radio"/> 25003	Click Strap, 12" Large	\$90
PERFORMANCE ENHANCEMENTS:		
<input type="radio"/> 25023	Crank Width Extensions <i>(+1.5" of crank width per set)</i>	\$20
<input type="radio"/> 25033	Heartrate Monitor	\$95
<input type="radio"/> 25024	Cordless Computer <i>(speed, distance, etc.)</i>	\$100
<input type="radio"/> 25022	Wheel Alignment Gauge	\$215
<input type="radio"/> 25036	Indoor Stationary Training Roller	\$250
TRANSPORTATION:		
<input type="radio"/> 25006	Crutch Holder	\$150
<input type="radio"/> 25025	Tow Bar for Wheelchair	\$225
BIKE RACKS:		
<input type="radio"/> 25007	Bike Rack, 1-1/4" Hitch	\$550
<input type="radio"/> 25008	Bike Rack, 2" Hitch	\$600
APPAREL:		
<input type="radio"/> 25084	Lemon Yellow Top End T-Shirt	\$20
<input type="radio"/> 25085	Grey Top End T-Shirt	\$20
<input type="radio"/> 25086	Grey Top End Racer Back Tank Top	\$20
<i>Please select a size for your shirt selection:</i>		
<input type="radio"/> 25087	Extra Small	
<input type="radio"/> 25088	Small	
<input type="radio"/> 25089	Medium	
<input type="radio"/> 25090	Large	
<input type="radio"/> 25091	Extra Large	
<input type="radio"/> 25092	XX Large	
<input type="radio"/> 25093	XXX Large	

HANDCYCLE SPECIFICATIONS:	
Weight:	45 Pounds
Turning Radius:	18 feet for 360 degree circle, 12 feet with a 3-point turn
Transfer Height:	9" from ground with 2" cushion
Frame Style:	Wrap around ovalized lower frame, 7000 series aluminum
Fork:	Removable, 7000 series aluminum, adjustable height of 22.5" to 24.5
Propulsion / Shifting:	Shimano Crank set and shifter, 10 X 3 (30) speeds
Cranks:	170mm standard bicycle crankset. Optional add-on V crank upgrade available
Back:	Angle adjustable straight back 35-70 degrees
Footrest:	Dual adjustable length, 5.5" wide, includes straps
Hand Pedals:	Ovalized ergonomic with frictionless bearings
Brake:	Caliper locking parking brake with lever mounted on the right handpedal. Includes a locking lever to be used as the parking brake.
Wheelbase:	60" Axle to axle
Camber:	1 Degree
Wheels:	559 metal spoked wheels
Axles:	Threaded
Tires/Tubes	559 Cruiser Tires (100 PSI)
Positioning Straps:	Hook and loop seat restraint included
Cushions:	Seat and back cushions included
Safety / Visibility:	Flag and wheel reflectors included

DEFINITIONS:	
Force-3:	The Force-3 has an updated and improved design. The frame is 17" with 15" seat and back cushions. Included with all Force-3's is a seat wedge which can be used to allow for a tighter fit into the frame.
Assembled vs Unassembled and Shipping:	If TE10046 (Unassembled) is selected the bike will require about 60 minutes of assembly upon delivery (instructions provided). If TE10047 (Assembled) is selected the bike will be fully assembled upon delivery. Shipping rate is \$150 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
V Cranks:	This upgraded crank option is great for customers whom require more clearance due to long legs or larger midsections. V crankset will be shipped separately and required self installation.
Straight Back:	The straight back measures 22.5" tall.
Disc Brake:	The Force-3 is manufactured with the welded mounts on the fork to allow for future disc brake installation. Parts required are available for purchase separately. Parts include a front wheel that is disc brake compatible, disc brake rotor and caliper.
Upgrade Accessories:	All Add-on options and accessories are boxed separately and require self installation.

Warranty Info:

Invacare warrants the frames and forks when purchased new and unused to be free from defects in materials and workmanship for a period of three (3) years from the date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Invacare warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty