

Pinniped™ Pediatric Commode

FAX TO 310-305-1718 OR SUBMIT ONLINE



Dealer: _____ Marked For: _____ Account #: _____

Date: _____ PO #: _____ Ordered By: _____

Phone: _____ Fax: _____ Email: _____

Ship to Address: _____

PRODUCT SPECIFICATIONS

\$373.86

PC 1000 BL

Seat Depth	15"
Seat Width	13"
Seat Height	10.5–14.5"
Armrest Height (from floor)	15.5–19.5"
Armrest Pad	2" W x 9" D
Overall Width	17"
Total Opening	6.5" W x 8" D
Weight	9.5 lb.
Carton Shipping Weight	13 lb.
Weight Limit	250 lb.
Warranty	Limited Lifetime

ACCESSORIES

Poziflex™ Chest Strap, Extra Small	\$96.41	PCB-900XS 3" Wide, 12.5" Padded Length
Poziflex Chest Strap, Small	\$96.41	PCB-900S 3.75" Wide, 15.75" Padded Length