

Ultra Lightweight Transport Chair



Item # 379B (Blue)
UPC # 652308125956
Item # 379R (Red)
UPC # 652308125963
Suggested Retail Price \$229.95

Specs

- **Weight Capacity:** 300 lbs.
- **Weight:** 17 lbs.
- **Seat Height:** 19"
- **Seat Dimensions:** 18"W x 16"D
- **Overall Dimensions:** 37.5"H x 23"W
- **Overall Width Folded:** 8.5"
- **Width Between Arms:** 17.5"

Features

- Lightweight aluminum compact frame **weighs only 17 lbs.** without footrests
- Large 8" wheels provide increased maneuverability
- Provides quick and easy transfer for shopping, appointments and sporting events
- Quick release fold down back allows the chair to be transported and stored in compact areas
- Quick adjust footrests swing away for safe patient transfer
- Padded armrests
- Seat belt included



See all of Nova's products, videos, custom mobility and bathroom safety tools, and so much more!

www.NovaMedicalProducts.com

Los Angeles
Chicago

© Copyright 2011



Phone: 1-800-557-6682

Fax: 1-800-551-1229

www.NovaMedicalProducts.com