

INVACARE®
820DLX
CARROLL DLX SERIES BED



820-3MDLX Bed featured with Invacare® ThinkSoft® Positioning Device and Kendall Bed Ends

The **Invacare® Carroll DLX Series 820DLX Bed** is the workhorse of long-term care, offering proven, dependable performance at an affordable price. This durable bed includes features such as an expandable length deck, removeable bed-end brackets and a built-in hand crank so the bed can be raised or lowered by hand in case of a power outage, or if the battery runs out.

Features

Mobility at Low Height – The 820DLX bed is secured to the floor in all elevated positions, providing reassurance that the bed is locked and secured for maximum resident safety.

Durable Grid Deck – Provides a supportive and comfortable sleep surface for the resident.

Ergonomic Hand Pendant – Waterproof, sealed hand pendant designed to fit comfortably in the resident's hand.


Dual Pendant Connections – Pendant connects to either side of the bed, providing customized convenience for the resident.

Expandable Length – Length expands from 76" to 80" with adjustable bed-end brackets.

Emergency Override – The 8203M-DLX bed features a built-in hand crank in case of a power outage.

820DLX Manual bed available.

Invacare® Carroll DLX Series beds are made in the USA.



1. Optional Pan Deck
2. Adjustable Mattress Deck
3. Optional 6-Function Ergonomic Hand Pendant



Invacare Corporation
www.invacare-cc.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 668-2337

Canada

994 Hargrieve Road
London, Ontario
N6E 1P5 Canada
(800) 668-2337

©2013
Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 12-638
121139

Model #	Height Range	Sleep Surface	Weight Capacity	Regulation Compliance	Limited Warranty
IH820-3MDLX	14" - 27" H	76" - 80" L 36" W	450 lb.	UL 60601-1 IEC 60601-2-38 CAN/CSA C22.2 No. 601.1	10 Years Frame 2 Years Mechanical 1 Years Electronics
IH820DLX	14" - 27" H	76" - 80" L 36" W	450 lb.	UL 60601-1 IEC 60601-2-38 CAN/CSA C22.2 No. 601.1	10 Years Frame 2 Years Mechanical

