



Invacare® IVC™ 9000 Topaz™ Wheelchair

Thank you for using Invacare eForms.

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

For additional information please visit www.invacare.com/eForms



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-678-4682

SUBMIT

CLEAR FORM

PRINT FORM



Yes, you can.®

Invacare® IVC™ 9000 Topaz™ Wheelchair



Account Information

Request Type: Quote

Order

Date: _____

Account #: _____

Company: _____

SHIP TO

Name: _____

Address: _____

City: _____ State: _____ Zip: _____

Country: _____

Purchase Order #: _____

CONTACT

Name: _____

Back Up Contact: _____

Phone: _____

Email: _____

Comments: _____

Special Note

The IVC 9000 Topaz with a 700 pound weight limit.

FRAME TYPE

☐ 9TPZ 9000 Topaz \$2,979.00 K0007

CUSTOM PACKAGES

☐ HDPG Heavy Duty Package \$321.30
Reinforced Frame Width weight capacity of 1,000 lbs. Only available in 28" and 30" Frame Widths.

FRAME WIDTH

☐ W20 20" Seat Width No Charge E2201
This selection is not available with HDPG option.

☐ W22 22" Seat Width No Charge E2201
This selection is not available with HDPG option.

☐ W24 24" Seat Width No Charge E2202
This selection is not available with HDPG option.

☐ W26 26" Seat Width No Charge E2202
This selection is not available with HDPG option.

☐ W28 28" Seat Width \$267.75 K0108

☐ W30 30" Seat Width \$267.75 K0108

FRAME DEPTH

☐ D18 18" Seat Depth No Charge

☐ D20 20" Seat Depth No Charge E2203

SEAT-TO-FLOOR HEIGHT

☐ ADULT Adult No Charge
Based on 24" wheels and 8" casters.

☐ H Hemi No Charge
Based on 24" wheels and 8" casters.

☐ 9TSH Super Hemi 15 1/2" No Charge
Not available with HDPG option.
Only available with 6" casters and 20" rear wheels.



Yes, you can.®

Invacare® IVC™ 9000 Topaz™ Wheelchair



ARM TYPE

Fixed height arms 9 1/2" adjustable height arms.

- ☐ T18 Full Length Fxd Ht ConvNo Charge
- ☐ T19 Full Length Adj Ht Conv\$134.95 E0973

- ☐ T28 Desk Length Fxd Ht ConvNo Charge
- ☐ T29 Desk Length Adj Ht Conv\$134.95 E0973

BACK TYPE

- ☐ ADJH Adj Height Back 16", 17", 18" with 10° bend
.....Standard
Adjustable back height 16"-18". Factory set at 16".

FRAME FINISH

- ☐ 71P Silver VeinNo Charge

- ☐ 24P Wet BlackNo Charge

CASTER OPTIONS

- ☐ 1259 6" x 2" Semi Pneu Casters \$103.89
Only available with super hemi seat-to-floor, NA with HDPG option.
- ☐ 1228BK 8" x 1 3/4" Semi Pneu Casters Black
.....No Charge
Only available with adult and hemi seat-to-floor, NA with HDPG option.

- ☐ 1228HD 8" x 2" Semi Pneu Casters BlackNo Charge
Only available with HDPG option.

WHEEL OPTIONS

- ☐ 5763BK 20" HD Pneu with Flat Free Inserts
.....\$215.27 E2211+E2213
Not available with HDPG option.
Only available with super hemi seat-to-floor.

- ☐ U5222SP Spoked HD Urethane Treaded Tires
.....No Charge

HANDRIM OPTIONS

- ☐ AHDRM Aluminum HandrimsStandard

- ☐ CHDRM Chrome HandrimsNo Charge

UPHOLSTERY COLORS

- ☐ U240 Black Nylon UpholsteryStandard
- ☐ U550 Midnight Blue Vinyl UpholsteryNo Charge

- ☐ U67 Black Vinyl UpholsteryNo Charge

FRONT RIGGINGS

- ☐ TZ33A Sw Footrests with Extra Large Footplates
.....\$205.52 K0041

- ☐ TZ34A Elevating Legrests with Extra Large Footplates ...
.....\$426.36 K0195 + K0041



Yes, you can.®

Invacare® IVC™ 9000 Topaz™ Wheelchair



FOOTREST/LEGREST ACCESSORIES

☐ 1337 Calf Strap \$38.37 K0038

☐ 1346 Heel Loops \$41.57 E0951
Standard on swingaway footrests.

SEAT POSITIONING STRAPS

☐ T1515 Topaz Push Button Style Seat Positioning Strap
..... \$53.55 E0978

BACK ACCESSORIES

☐ T1483 Push Bar \$166.01
This accessory is shipped in a separate carton.

SEAT ACCESSORIES

☐ T9120 Solid Seat Insert \$160.65
This is shipped in a separate carton.
Solid Seat Insert is manufactured in a vinyl upholstery only.

ACCESSORIES

The below selections are shipped in a separate carton.

☐ T1496 IV/Oxygen Holder Package Assembly
..... \$169.22 K0105 + E2208

☐ T1497 Oxygen Holder Package Assembly
..... \$161.72 E2208

☐ T1322 IV Rod and Bracket, Telescopic \$158.51 K0105
This selection is not available with super hemi seat-to-floor.

MISCELLANEOUS ACCESSORIES

☐ T1358 Rear Anti-Tippers \$92.11 E0971

☐ A1567 Pull to Lock Clamp on Whl Locks (Std on Hemi) ...
..... No Charge

☐ 1721 Hi Push to Lock Clamp on Whl Locks (Std on Adult)
..... No Charge
Only available with HDPG option.

☐ 1305L Crutch and Cane Carrier \$46.05 E2207

☐ A1576 Wheel Lock Extension \$41.77 E0961

☐ 1335XL Non Folding Device - Topaz \$53.55



Yes, you can.®

Invacare® IVC™ 9000 Topaz™ Wheelchair



Seat to Floor Matrix 6" Casters

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
15 1/2"	1 Position Only	20"	Top

Seat to Floor Matrix 8" Casters

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
17 1/2"	Top	24"	Top
19 1/2"	Bottom	24"	Bottom

Weight Limit 700 Pounds

Invacare recommends all wheelchairs to be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular limitations of the product.



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-678-4682

SUBMIT

CLEAR FORM

PRINT FORM