INVACARE® 9000 SL WHEELCHAIR





The **Invacare® 9000 SL Wheelchair** is the perfect high-strength lightweight rental product for those who want quality and durability at a moderate price. This low-maintenance wheelchair also has a full complement of options and accessories. Its contemporary design and high-performance frame make the 9000 SL wheelchair a superior value.

Features

Stylish, contemporary carbon steel frame is long-lasting

Multiple frame color variations for customizing

Dual axle positions allow variations of seat-to-floor heights

Urethane rear tires mounted on "no flex" wheels offer superior performance

Urethane casters for better ride and lighter weight

Low-maintenance, dual-sealed precision bearings

Heavy-duty inner liners keep seat and back from stretching

Space-saver arms reduce overall frame width without limiting seat width

Meets RESNA standards for manual wheelchairs in product performance

HCPCS Code: K0004











- 1. Dual axle positions
- 2. Urethane casters
- 3. Space saver arms
- 4. Seating options available
- 5. Bulider options available



Invacare Corporation

www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668 - 5324

©2012

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 11-226 111241















Model # Overall Width Overall Height 9SL 23", 25"

36" 30"

Overall Depth Folded w/o riggings Width

Seat Width Seat (2" increments) Depth

Seat-to-Floor Height (1" increments) 16", 18" 15.5"-19.5"







15"



Shipping

16"-20"



Carton





Back Height 9SL 15"-17", 17"-19"

Caster Height 6", 8"

Wheel Height 20", 22", 24" Weight Capacity 250 lb.

Weight 43 lb.

Dimensions 32"x13"x 37"

Limited Warranty 5 Years

Yes, you can.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.