



AERO T



THE TILITE AERO T. PURE GEOMETRY.

AERO T INNOVATION.

The soul of a chair lies in how it rolls, how it handles the real world, how it puts the power to the ground, how it manages the next curb, and how it tracks down a hall or a sidewalk. How a chair does all these things is a function of the chair's geometry, and that is the heart of the innovation of the Aero T. The lengths of its tubes and the angles at which they interface are the keys to the exceptional performance of the Aero T.

- **New, Swept-In, Adjustable Dual-Tube Frame.**
 - + Advance frame geometry insures:
 - Extraordinary handling on any surface
 - Superior rigidity
 - Optimal energy transfer
 - + Sweeping curve of the lower frame tube creates a unique look as the lower frame blends in and fades away.
- Curved titanium backrest tube further enhances the clean new look.
- **TiShaft** back release bar makes folding and unfolding the backrest easy from any angle.
- **Adjustable Camber Mount** incorporates **TiLite's Tru-Fit System** with up to 3.5" of Rear Seat Height adjustment.
- **SlipStream** Single-Sided Fork with Integrated Bumper (patent pending).

TECHNICAL SPECIFICATIONS.

Frame Type	Aluminum Dual-Tube Rigid - Full Adjustability
Frame Styles	SuperLite Edition Armed Forces Edition Power Adaptable Frame Reinforced Frame
Total Chair Width*	20"
Weight Limit	265 lbs.
Seat Widths	12" - 20"
Seat Depths	12" - 20"
Front Seat Heights	15" - 21"
Rear Seat Heights	13" - 21"
Front Frame Angles	70°, 75°, 80°, 85° or 90°
Footrest Width Options	Standard Taper
Cambers	0°, 2°, 4°, 6°, 8°, 12°, 15° or 18°
Front Wheels	23 Options
Rear Wheels	67 Options
Handrims	15 Options
Center of Gravity Adjustment	5"
Transport Weight*	12 lbs.
Aesthetics	6 Aluminum Frame Finishes 19 Frame Colors 6 Color Anodize Packages 7 Ultrasuede Colors 8 Spinergy Wheel Colors
HCPCS Code	K0005 (Code Pending)
Warranty	Lifetime on Frame; Limited Lifetime on Components

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lbs.

CUSTOM CONFIGURATION THE TILITE WAY.

Custom configuration is when you are able to select options and components on a chair to meet your needs. While other companies offer some custom configuration, TiLite offers the widest range of options from which to choose. More choices mean more ability to truly create a custom configuration of options that best meets your needs. Choice Matters™ and that is why TiLite offers more options than any other manufacturer! The choice is yours.



TiLite | 800.545.2266 | tilite.com
© 2013 TiSport, LLC.

All TiLite wheelchairs are shown with available options. CF_AeroT_0113RevA